
Chemistry Tubes

Tubos de bioquímica

Tubes chimiques

Tubos Químicos

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Becton, Dickinson and Company
Franklin Lakes, NJ USA Made in USA www.bd.com
BD, BD Logo and all other trademarks are property
of Becton, Dickinson and Company. ©2010 BD 8363931
03/2011

Chemistry Tubes EN

Intended Use:

BD Microtainer® Tubes are used to collect, transport and store skin puncture blood specimens for tests utilizing serum or plasma.

Product Description:

The BD Microtainer® Tube consists of a plastic reservoir and a color-coded closure (Fig. 1). The upper edge of the reservoir serves as a collector. The closure fits snugly over the upper edge of the reservoir to reduce the potential for sample contact when the cap is removed. The tubes contain additives in varying concentrations dependent upon the fill volumes and the required additive to blood ratio for the tube. Additive choice depends on the analytic test method, and is specified by the manufacturer of the test reagents and/or instrument on which the test is performed.

	Color	Additive	Fill Line
Serum tubes:			
Silicone Coated	Red	No Additive	250-500 µL
SST™ Clear/Amber	Gold	Clot Activator & Polymer Gel	400-600 µL
Plasma tubes:			
Lithium Heparin	Green	Lithium Heparin	200-400 µL
PST™ Clear/Amber	Mint Green	Lithium Heparin & Polymer Gel	400-600 µL
Glycolytic Inhibitor	Gray	NaF/Na2 EDTA	400-600 µL

To remove the closure, grasp the ribbed edge and twist while pulling it off. To reinsert, twist the closure while pushing it onto the reservoir.

An optional extender is available. It fits into the bottom of the reservoir and increases the tube length to about 75 millimeters (mm). With this extender, the tube fits in a standard 13 x 75 mm test tube rack. The extender also provides a larger area for specimen labeling.

Serum Tubes:

A silicone coating on the walls of the serum tubes reduces adherence of red cells to tube walls.

SST™ tubes are coated with micronized silica particles to increase surface area. Particles in the white film on the interior surface activate clotting when tubes are mixed 5 times by inversion.

Plasma Tubes:

Plasma tubes contain the optimum amount of additive to anticoagulate the specified volume of skin puncture blood. Glycolytic Inhibitor tubes also contain the optimum amount of additive to prevent glycolysis for up to 24 hours at room temperature.

Amber Tinted Tubes:

The tinting of the amber tube significantly reduces light transmission and is used when protection is required for light sensitive analytes.

Storage:

Store product below 25°C and use prior to expiration date.

Limitations:

- BD Microtainer® Tubes are designed for skin puncture blood collections.
- Current CLSI guidelines recommend that serum or plasma be physically separated from contact with cells as soon as possible, with a maximum time limit of two hours from the time of collection. In addition, concentrations of chemical constituents in serum or plasma from skin puncture blood have been compared to those in serum from venous blood. Statistically and/or clinically important differences in the concentrations of glucose, potassium, total protein, and calcium have been reported. Except for glucose, the concentration of these analytes is lower in skin puncture blood.
- Blood fill quantity must be within the specified range to ensure proper blood to additive ratio for adequate mixing and accurate test results.
- Glycolytic Inhibitor tubes prevent glycolysis for up to 24 hours at room temperature.

Materials Required But Not Provided:

- Gloves for the person performing the collections
- Alcohol swabs
- Dry sterile gauze pads
- Lancet as appropriate for site and volume of blood required
- Sharps and contaminated materials disposal container

Optional Materials For Specimen Collection:

- Warming device if required, dependent on the volume of blood and tests to be performed.
- Adhesive bandage. Avoid use of bandages with patients likely to place fingers or feet in their mouths, as ingestion / aspiration may occur.

Directions For Use:

- Gloves must be worn during skin puncture blood collection.
- Select puncture site, warm as appropriate, and cleanse with alcohol.
- Allow site to dry. Do not dry by wiping, as disinfection occurs during air drying.
- Remove closure from tube and place on a convenient surface. The closure may be nested on the base of the tube.
- Puncture skin with the appropriate lancet, following instructions supplied by the manufacturer. Dispose of lancet in an approved sharps container. Wipe away the first drop of blood with gauze.
- Hold BD Microtainer® Tube at 30°-45° angle from surface of puncture site (Fig. 2). Touch collector end of reservoir to drop of blood. Avoid scraping skin surface to collect blood sample. After collecting 2 or 3 drops, blood will freely flow down reservoir wall to bottom of tube. CAUTION: "Milking" of skin puncture site may cause hemolysis and adversely affect test result accuracy.
- Fill tubes between the fill marks. Over- or underfilling of tube may result in clotting and/or erroneous test results.
- Replace closure by twisting and pushing cap downward.

9. Invert plasma tubes 8 - 10 times. Invert SST™ tubes 5 times.

10. Fill other BD Microtainer® Tubes as required. NOTE: If capillary blood gas testing is to be done, the capillary tubes used to collect samples must be filled before other samples are collected. Otherwise, the EDTA samples for hematology testing should be filled first.

11. Dispose of the used lancet into an approved sharps container if this has not been done. Dispose of any contaminated materials into appropriate container.

Procedure For Serum/Plasma Separation:

- Serum Separation: allow blood to clot for a minimum of 30 minutes. Plasma Separation: invert tube 8-10 times to assure anticoagulation.
- Place sample in centrifuge or centrifuge adapter as necessary, taking care to balance the system.
- Centrifuge as follows:

	Gel Tubes	Non-Gel Tubes
Minimum Time	90 seconds	3 minutes
Rate	6000-15000 g	Minimum 2000 g

4. Speed in revolutions per minute (rpm) may be related to relative centrifugal force (g) by the following formula: $RCF \text{ (in g forces)} = 1.12 \times 10^{-5} \times (\text{rpm})^2 \times r \text{ (in cm)}$; where "r" is the radial distance from the center of the centrifuge head to the bottom of the tube.

5. Remove tube from centrifuge. Serum/Plasma is ready for use and may be pipetted directly into analyzer cup.

Principle For Gel Separation:

BD Microtainer® Gel Tubes contain a polymer barrier material with a specific gravity between that of the serum/plasma and blood cells. During centrifugation of the blood specimen, the polymer barrier material rises to the interface of the serum/plasma and the clot/cells and forms a physical barrier separating the serum/plasma from the clot/cells.

Caution:

- Practice Universal Precautions. Use gloves, gowns, eye protection, other personal protective equipment, and engineering controls to protect from blood splatter, blood leakage, and potential exposure to bloodborne pathogens.
- Handle all biologic samples and blood collection "sharps" (lancets, needles, luer adapters, and blood collection sets) according to the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury) since they may transmit viral hepatitis, HIV (AIDS), or other infectious diseases. Utilize any built-in needle protector if the blood collection device provides one. BD does not recommend resheilding used needles. However, the policies and procedures of your facility may differ and must always be followed.
- Discard all blood collection "sharps" in biohazard containers approved for their disposal.
- Transferring a sample from a syringe to a tube is not recommended. Additional manipulation of sharps increases the potential for needlestick injury. In addition, depressing the syringe plunger during transfer can displace the sample and cause a potential blood exposure. Using a syringe for blood transfer may also cause over or under filling of tubes, resulting in an incorrect blood-to-additive ratio and potentially incorrect analytic results.
- If blood is collected through an intravenous (I.V.) line, ensure that line has been cleared of I.V. solution before beginning to fill blood collection tubes. This is critical to avoid erroneous laboratory data from I.V. fluid contamination.
- Underfilling of tubes will result in an incorrect blood-to-additive ratio and may lead to incorrect analytic results or poor product performance.

Analytic Equivalency:

Evaluations of BD Microtainer® Tubes have been performed for an array of analytes over a variety of test methods and time periods. The BD Technical Services Department is available to answer questions regarding these studies. Please contact them to obtain references and technical reports on these evaluations and any other information regarding the use of BD Microtainer® Tubes with your instrument/reagent system.

It is the laboratory's ultimate responsibility to determine reference intervals for all analytes based upon the tubes used for skin puncture blood specimens by that laboratory. The clinical laboratory should establish/verify its reference ranges if changing specimen collection tube types and sizes, as this could potentially affect analytic results from patient samples.

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