Single Patient Use Hand Controlled & Foot Controlled Electrosurgical Suction Coagulators

Sterility guaranteed unless package has been opened, broken or damaged.

Intended Use:

The Suction Coagulator is used for coagulating tissue and removing fluid and debris from the surgical site. Suction and coagulation can be performed independently or simultaneously.

Maximum voltage is not to exceed 4.5 KV peak.

Contraindications:

These devices should never be used when:

- There is visible evidence of damage to the exterior of the device such as cuts, punctures, nicks, abrasions, unusual lumps, or significant discoloration.
- In the presence of flammable gases, flammable prep solutions or drapes, oxidizing gases such as nitrous oxide, or in oxygen-enriched environments.
- Where conventional monopolar electrosurgery is inappropriate or unsafe, such as for small appendages, in circumcision or finger surgery.

Warning:

- Use lowest possible power setting on the associated electrosurgical unit capable of achieving desired surgical effect. Activation time should be as short as possible.
- Do not directly secure the cable with metal instruments to surgical drapes. Activation when in contact with metal instruments may cause burns at the tissue / instrument interface.
- To avoid burns never allow cable associated with this device to be in contact with skin of patient or touching operator
- Do not permit the cables connected to these devices to be parallel and in close proximity to the leads of other electrical devices.
- Do not immerse or saturate suction coagulator or cord in fluid.
- Do not pull on or stretch the cord.
- Do not test by sparking the active electrode to the patient plate or other objects.
- Always place unused electrosurgical accessories in a safe insulated location such as a holster when not in use.
- Inspect and test each device before use.
- Electrosurgery is potentially hazardous for patients with pacemakers, AICDs, neurostimulators or other active implants. High frequency energy may damage or adversely affect active implants. Ventricular fibrillation may occur.
- After activation the suction coagulator may be hot and care should be taken to prevent accidental burns.
- Do not activate the electrosurgical unit if the tip of the instrument is not in a position to deliver energy to target tissue. Doing so may cause inadvertent burns to the patient.

Precaution:

- Refer to the electrosurgical generator operator manual prior to use.
- After patient use, this product should be considered a potential biohazard and disposed of per hospital policy and applicable law.

Inspection:

These devices should be inspected before each use and should not be used if damage is found. Visually examine the devices for obvious physical damage including:

- Cracked, broken or otherwise distorted plastic parts.
- Broken or significantly bent connector contacts.
- Verify that the electrode is fully and securely seated in the handpiece before use.

Instructions for use:

1. Attach suction tubing to suction connector at end of handle.
2. For hand switching suction coagulator, plug the 3-prong plug into the appropriate receptacle on the electrosurgical unit. For foot switching suction coagulator, use CONMED adapter (Catalog No. 60-0823-001) to connect cable end to the foot controlled active jack of the electrosurgical unit as required.
3. Remove the cap and cleaning rod before operation. Activate the coagulation mode by depressing the activation button on the top of the hand controlled device, or by depressing the COAG footswitch for the foot controlled device. For suction, place finger over the suction port. Suction and coagulation may be activated simultaneously.
4. User may bend suction tube to meet surgical requirements. It is recommended that the cleaning rod be in the tube while forming, to prevent kinking of the tube.
5. To clear eschar clogs in tubing, fully insert cleaning rod through device.
6. To remove eschar buildup on end of tube, wipe with moistened surgical sponge.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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© CONMED Corporation、米国で出版
MDSS GmbH