
Revisions

Rev from	Rev to	JOB #
06	07	52-17

Notes:

- BD Catalog Number: 212525, 212526, 212527, 212528, 212531, 212532, 212539, 212542, 212543, 212544, 212545
- Blank (Sheet) Size: Length: 25.5" Width: 22"
- Number of Pages: 28 Number of Sheets: 1
- Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 4.25" x 5.5"
- Ink Colors: No. of Colors: 2 PMS#: 032 Red; Standard Black
- Printed two sides: Yes  No
- Style (see illustrations below): # 5

- Vendor Printed  Online/In House Printed  Web
- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	<b>REVISED BY</b> By Sonia Thompson at 4:48 pm, Oct 20, 2017	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	<b>PROOFING APPROVED BY</b> By Natalie Morio at 1:29 pm, Oct 26, 2017			
Checked By	<b>THIRD EYE BY</b> By Nichole Graham at 2:03 pm, Nov 29, 2017			
Part Number:	8820191JAA	Category and Description Package Insert, Gram Stain Kits and Reagents	Sheet: 1 of 29 Scale: N/A	<b>A</b>

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Gram Stain Kit		1	Cat. No
Gram Crystal Violet	For the differential staining of bacteria.	1 x 250 mL	<b>212539</b>
Gram Iodine (Stabilized)		1 x 250 mL	
Gram Decolorizer		1 x 250 mL	
Gram Safranin		1 x 250 mL	
Gram Crystal Violet	For staining microorganisms by the differential Gram method.	4 x 250 mL	<b>212525</b>
		1 x 3.8 L	<b>212526</b>
Gram Iodine (Stabilized)	For staining microorganisms by the differential Gram method.	4 x 250 mL	<b>212542</b>
		1 x 3.8 L	<b>212543</b>
Gram Decolorizer	For staining microorganisms by the differential Gram method.	4 x 250 mL	<b>212527</b>
		1 x 3.8 L	<b>212528</b>
Gram Safranin	For staining microorganisms by the differential Gram method.	4 x 250 mL	<b>212531</b>
		1 x 3.8 L	<b>212532</b>
Gram Basic Fuchsin	For staining microorganisms by the differential Gram method.	4 x 250 mL	<b>212544</b>
		1 x 3.8 L	<b>212545</b>

## INTENDED USE

Gram Stain Kits and Reagents are used to stain microorganisms from cultures or specimens by the differential Gram method.

## SUMMARY AND EXPLANATION

The Gram stain was devised in 1884 by Christian Gram in an attempt to differentiate bacterial cells from infected tissue. Although Gram observed what is now called the “Gram reaction,” he did not recognize the taxonomic value of his technique. The Gram stain is now used to differentiate intact, morphologically similar bacteria into two groups based on cell color after staining. In addition, cell form, size and structural details are evident. Such preliminary information provides important clues to the type of organism(s) present and the further techniques required to characterize them.

Because inorganic iodine is rapidly oxidized and loses its effectiveness as a mordant,<sup>1</sup> the Gram Stain Kit (Cat. No. 212539) differs from Gram’s original formulation by offering a more stable organic iodine complex, L-polyvinylpyrrolidone-iodine.

## PRINCIPLES OF THE PROCEDURE

The Gram stain procedure<sup>2</sup> consists of:

Staining a fixed smear with crystal violet.

Applying iodine as a mordant.

Decolorizing the primary stain with alcohol/acetone; and, counterstaining with safranin or basic fuchsin.

A crystal violet-iodine complex forms in the protoplast (not the cell wall) of all organisms stained by this procedure. Organisms able to retain this dye complex after decolorization are classified as gram-positive while those that can be decolorized and counterstained are classified as gram-negative.

Upon disruption or removal of the cell wall, the protoplast of gram-positive (as well as gram-negative) cells can be decolorized and the gram-positive attribute lost. Thus, the mechanism of the Gram stain appears to be related to the presence of an intact cell wall able to act as a barrier to decolorization of the primary stain.

Generally, the cell wall is nonselectively permeable. It is theorized that during the Gram stain procedure, the cell wall of gram-positive cells is dehydrated by the alcohol in the decolorizer and loses permeability, thus retaining the primary stain. However, the cell wall of the gram-negative cells has a higher lipid content and becomes more permeable when treated with alcohol, resulting in loss of the primary stain.

The molecular basis for the Gram stain has not yet been determined.

## REAGENTS

Approximate Formula\* Per Liter

### Gram Crystal Violet

#### PRIMARY STAIN

Crystal Violet.....	3.0 g
Isopropanol.....	50.0 mL
Ethanol/Methanol.....	50.0 mL
Distilled Water.....	900.0 mL

### Gram Iodine

#### MORDANT

(Working solution prepared from Gram Diluent and Gram Iodine 100X)

Iodine Crystals.....	3.3 g
Potassium Iodide.....	6.6 g
Distilled Water.....	1.0 L

### Stabilized Gram Iodine

#### MORDANT

Polyvinylpyrrolidone-Iodine Complex.....	100.0 g
Potassium Iodide.....	19.0 g
Distilled Water.....	1.0 L

### Gram Decolorizer

#### DECOLORIZER

Acetone.....	250.0 mL
Isopropanol.....	750.0 mL

### Gram Safranin

#### COUNTERSTAIN

Safranin O Powder (pure dye).....	4.0 g
Ethanol/Methanol.....	200.0 mL
Distilled Water.....	800.0 mL

### Gram Basic Fuchsin

#### COUNTERSTAIN

Basic Fuchsin.....	0.08 g
Phenol.....	2.6 g
Isopropyl Alcohol.....	4.5 mL
Distilled Water.....	993.0 mL

\*Adjusted and/or supplemented as required to meet performance criteria.

**Warnings and Precautions:** For *in vitro* Diagnostic Use.

Over time, a fine precipitate may develop in Gram Basic Fuchsin. Product performance will not be affected.

### BD BBL Gram Crystal Violet (212525, 212526)

#### Warning

**H332** Harmful if inhaled. **H351** Suspected of causing cancer. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

**P201** Obtain special instructions before use. **P202** Do not handle until all safety precautions have been read and understood.

**P261** Avoid breathing dust/fume/gas/mist/vapors/spray. **P271** Use only outdoors or in a well-ventilated area. **P273** Avoid release to the environment. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P304+P340** IF

INHALED: Remove person to fresh air and keep comfortable for breathing. **P308+P313** IF exposed or concerned: Get

medical advice/attention. **P312** Call a POISON CENTER or doctor/physician if you feel unwell. **P405** Store locked up.

**P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

### BD BBL Gram Decolorizer (212527, 212528)

#### Danger

**H225** Highly flammable liquid and vapor. **H319** Causes serious eye irritation. **H336** May cause drowsiness or dizziness.

**P210** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. **P240** Ground/bond

container and receiving equipment. **P241** Use explosion-proof electrical/ventilating/lighting/...equipment. **P242** Use only non-

sparkling tools. **P243** Take precautionary measures against static discharge. **P261** Avoid breathing dust/fume/gas/mist/vapors/

spray. **P264** Wash thoroughly after handling. **P271** Use only outdoors or in a well-ventilated area. **P280** Wear protective gloves/

protective clothing/eye protection/face protection. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated

clothing. Rinse skin with water/shower. **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing.

**P305+P351** IF IN EYES: Rinse cautiously with water for several minutes. **P312** Call a POISON CENTER or doctor/physician if

you feel unwell. **P337+P313** If eye irritation persists: Get medical advice/attention. **P338** Remove contact lenses, if present and easy to do. Continue rinsing. **P370+P378** In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray. **P403+P233** Store in a well-ventilated place. Keep container tightly closed. **P403+P235** Store in a well-ventilated place. Keep cool. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

#### BD BBL Gram Safranin (212531, 212532)

##### Warning

**H226** Flammable liquid and vapor.

**P210** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. **P241** Use explosion-proof electrical/ventilating/lighting/.../equipment. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P240** Ground/bond container and receiving equipment. **P233** Keep container tightly closed. **P242** Use only non-sparking tools. **P243** Take precautionary measures against static discharge. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P370+P378** In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray. **P403+P235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

#### BD BBL Gram Stain Kit (212539)

##### Danger

**H225** Highly flammable liquid and vapor. **H226** Flammable liquid and vapor. **H315** Causes skin irritation. **H319** Causes serious eye irritation. **H332** Harmful if inhaled. **H336** May cause drowsiness or dizziness. **H351** Suspected of causing cancer. **H401** Toxic to aquatic life. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

**P201** Obtain special instructions before use. **P202** Do not handle until all safety precautions have been read and understood. **P210** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. **P233** Keep container tightly closed. **P240** Ground/bond container and receiving equipment. **P241** Use explosion-proof electrical/ventilating/lighting/.../equipment. **P242** Use only non-sparking tools. **P243** Take precautionary measures against static discharge. **P261** Avoid breathing dust/fume/gas/mist/vapors/spray. **P264** Wash thoroughly after handling. **P271** Use only outdoors or in a well-ventilated area. **P273** Avoid release to the environment. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P302+P352** IF ON SKIN: Wash with plenty of soap and water. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P308+P313** IF exposed or concerned: Get medical advice/attention. **P312** Call a POISON CENTER or doctor/physician if you feel unwell. **P321** Specific treatment (see on this label). **P332+P313** If skin irritation occurs: Get medical advice/attention. **P337+P313** If eye irritation persists: Get medical advice/attention. **P362+P364** Take off contaminated clothing and wash it before reuse. **P370+P378** In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray. **P403+P233** Store in a well-ventilated place. Keep container tightly closed. **P403+P235** Store in a well-ventilated place. Keep cool. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

#### BD BBL Gram Iodine (Stabilized) (212542, 212543)

##### Warning

**H315** Causes skin irritation. **H319** Causes serious eye irritation. **H401** Toxic to aquatic life.

**P264** Wash thoroughly after handling. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P321** Specific treatment (see on this label). **P332+P313** If skin irritation occurs: Get medical advice/attention. **P337+P313** If eye irritation persists: Get medical advice/attention. **P350** Gently wash with plenty of soap and water. **P362+P364** Take off contaminated clothing and wash it before reuse. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

**Storage:** On receipt, store at 15–30 °C. The expiration date is for product in unopened bottles stored as directed. Do not open until ready to use.

Use the traditional Gram Iodine working solution within 3 months of preparation, not exceeding the Expiry of either component.

**Product Deterioration:** Unstabilized Gram Iodine when reconstituted may cause variability in the Gram stain when sufficient iodine is no longer available in solution. Protect the iodine solution from undue exposure to air, light and heat, to ensure that the solution is providing proper mordant activity.

## SPECIMEN COLLECTION AND PREPARATION

Apply the test specimen to a clean glass slide in a manner that will yield a thin, uniform smear. Emulsify colonies from an 18–24 h culture in saline to obtain the proper density.

Allow the smear to air dry.

Fix the smear to the slide using one of the following techniques:

1. Heat fix by passing the slide through a low flame 2–3 times. Cool the slide to room temperature before staining.  
*NOTE:* Do not overheat the slide; excessive heating will cause atypical staining.
2. Methanol fix the slide by flooding with absolute methanol for 1–2 min and rinse with tap water before staining.<sup>3</sup>  
*NOTE:* For proper fixation, store absolute methanol in a brown screw-capped bottle and replenish the working supply every two weeks.

## PROCEDURE

### Reagent Preparation

Prepare the traditional Gram Iodine working solution by adding an entire 2.5 mL ampule of Gram Iodine 100X to 250 mL Gram Diluent or an entire 40 mL vial of Gram Iodine 100X to 3.8 L of Gram Diluent; mix thoroughly.

**Materials Provided:** Gram Crystal Violet, Gram Iodine or Stabilized Gram Iodine, Gram Decolorizer and Gram Safranin or Gram Basic Fuchsin.

**Materials Required But Not Provided:** Microscope slides, bunsen burner or methanol, bacteriological loop, swabs, blotting paper, microscope with oil immersion lens and Gram slide.

### Test Procedure:

1. Flood the fixed smear with primary stain (Gram Crystal Violet) and stain for 1 min.
2. Remove the primary stain by gently washing with cold tap water.
3. Flood the slide with mordant (Gram Iodine or Stabilized Gram Iodine) and retain on the slide for 1 min.
4. Remove the mordant by gently washing with tap water.
5. Decolorize (Gram Decolorizer) until solvent running from the slide is colorless (3–60 sec).
6. Wash the slide gently in cold tap water.
7. Flood the slide with counterstain (either Gram Safranin or Gram Basic Fuchsin) and stain for 30–60 sec.
8. Wash the slide with cold tap water.
9. Blot with blotting paper or paper towel or allow to air dry.
10. Examine the smear under an oil immersion lens.

### User Quality Control

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Run controls using **BBL™** Gram Slide (Cat. No. 231401) or 18–24 h cultures of known gram-positive and gram-negative microorganisms. The following test strains are recommended:

Organism	ATCC®	Expected Results
<i>Staphylococcus aureus</i>	25923	gram-positive cocci
<i>Escherichia coli</i>	25922	gram-negative rods

## LIMITATIONS OF THE PROCEDURE

The Gram stain provides preliminary identification information only and is not a substitute for cultural studies of the specimen. Gram stain results must be confirmed with additional procedures such as direct antigen tests and culture on media.

Prior treatment with antibacterial drugs may cause gram-positive organisms from a specimen to appear gram-negative.

Use of an 18–24 h culture is advisable for best results since fresh cells have a greater affinity than old cells for most dyes. This is particularly true of many spore formers, which are strongly gram-positive when examined in fresh cultures but which later become gram-variable or gram-negative.

The Gram stain reaction is altered by physical disruption of the bacterial cell wall or protoplast. The cell walls of gram-positive bacteria interpose a barrier which prevents leaching of the dye complex from the cytoplasm. Cell walls of gram-negative bacteria contain lipids soluble in organic solvents, which are then free to decolorize the cytoplasm. Therefore, a microorganism that is physically disrupted by excess heating will not react to Gram staining as expected.

"Careful adherence to procedure and interpretive criteria is required for accurate results. Accuracy is highly dependent on the training and skill of the microbiologist."<sup>2</sup>

Gram stain results, including organism morphology, can be affected by the age of the isolate, bacteria containing autolytic enzyme systems, cultures transferred from antibiotic-containing media, as well as specimens collected from patients on antibiotics.<sup>4</sup> "Background material and artifacts can also interfere with interpretation. Precipitated gram-positive stain generally appears as irregular coccoid shapes or as asters resembling fungal hyphae."<sup>4</sup>

#### EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS<sup>1-4</sup>

Reaction	Using Gram Safranin	Using Gram Basic Fuchsin
Gram-positive	Purple-black cells	Bright purple to purple-black cells
Gram-negative	Pink to red cells	Bright pink to fuchsia cells

#### REFERENCES

1. Magee, C.M., G. Rodenheaver, M.T. Edgerton, and R.F. Edlich. 1975. A more reliable Gram staining technique for diagnosis of surgical infections. *Am. J. Surg.* 130:341–346.
2. Kruczak-Filipov, P., and R.G. Shively. 1994. Gram stain procedures, p. 1.5.1–1.5.18. *In* H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 1. American Society for Microbiology, Washington, D.C.
3. Mangels, J.I., M.E. Cox, and L.H. Lindley. 1984. Methanol fixation. An alternative to heat-fixation of smear. *Diag. Microbiol. Infect. Dis.* 2:129–137.
4. Chapin, K. 1995. Clinical microscopy, p. 33-51, *In* P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed., American Society for Microbiology, Washington, D.C.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or [www.bd.com](http://www.bd.com).

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For US: "For Investigational Use Only"

Lower limit of temperature / Донен лимит на температурата / Dolni hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κάτω όρο θερμότητας / Limite inférieure de température / Alumine temperatuuripiir / Limite inférieure de température / Najnižja dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inférieure de température / Температурный нижний порог / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite máximo de temperatura / Limitá mínima de temperaturá / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sıcaklık alt sınırı / Минимальна температура / 温度下限

Control / Контроль / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controlle / Controllo / Контроль / kontroll / Контроль / 对照

Positive control / Положительный контроль / Pozitívni kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positivite kontrol / Contrôle positif / Pozitivna kontrola / Pozitiv kontrol / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozitivná kontrolė / Positivite controle / Kontrola dodatna / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивный контроль / 阳性对照试剂

Negative control / Отрицательный контроль / Negatívni kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negativne kontrol / Contrôle négatif / Negativna kontrola / Negativ kontrol / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiama kontrolė / Negativná kontrolė / Negativite controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативный контроль / 阴性对照试剂

Method of sterilization: ethylene oxide / Метод на стерилизация: этиленов оксид / Způsob sterilizace: etylenoxid / Steriliseringsmetode: etylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Steriliseerimismetood: etüleenoksiid / Méthode de stérilisation: oxyde d'éthylène / Metoda sterilizacije: etilen oksid / Sterilizálás módszere: etilén-oxid / Método de esterilização: óxido de etilene / Стерилизация әдісі – этилен тотығы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizācijas metode: etilēnoksīds / Gesteriliseerid met behulp van etylenoxide / Steriliseringsmetode: etylenoxid / Metoda sterilizacji: etilen oksid / Método de esterilização: óxido de etileno / Método de sterilizace: oxid de etilena / Метод стерилизации: этиленоксид / Metoda sterilizacije: etylenoxid / Metoda sterilizacije: etilen oksid / Steriliseringsmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: этиленоксидом / 灭菌方法: 环氧乙烷

Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringsmetode: bestråling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseerimismetood: kiirgus / Méthode de stérilisation: irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Método de sterilizazione: irradiazione / Стерилизация әдісі – сәулө түсіу / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizācijas metode: apstarošana / Gesteriliseerid met behulp van bestraling / Steriliseringsmetode: bestråling / Metoda sterilizacji: napromienianie / Método de esterilização: irradiação / Metoda de sterilizace: iradiere / Метод стерилизации: облучение / Metoda sterilizácie: ožiarenie / Metoda sterilizacije: ozračavanje / Steriliseringsmetod: stråling / Sterilizasyon yöntemi: irradyasyon / Метод стерилизації: опроміненням / 灭菌方法: 辐射

Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικά κίνδυνοι / Riesgos biológicos / Biologiskelised riskid / Risques biologiques / Biološki rizik / Biológiallag veszélyes / Rischio biologico / Биологичны тауекелдер / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biološki rizici / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscurosi biologice / Биологическая опасность / Biologické riziko / Biološki riziko / Biologisk risk / Biyolojik Riskler / Биологична небезпека / 生物学风险

Caution, consult accompanying documents / Внимание, направте справка в придружаващите документи / Pozor! Prostudujte si příloženou dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, ομψβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaustus! Lugeda kaasvatustsiooni / Attention, consulter les documents joints / Uprozeženje, koristí prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлануу, тиэти құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dmesio, žižrèkite pridedamus dokumentus / Piesardzība, skatit ravaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktigt, se vedlagt dokumentasjon / Nalezý zapozná sač s dolozconými dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Vstraha, pozri sprievodné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se meðfylgjande dokumentation / Dikkat, birlikte verilen belgelerle / Увага: див. супутню документацию / 小心, 请参阅附带文档。

Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Övre temperaturgrænse / Temperaturobergrenze / Ανώτερο όρο θερμότητας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superior de temperatura / Температурный верхний порог / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurlimiet / Övre temperaturgrænse / Górná granica temperatury / Limite máximo de temperatura / Limitá máxima de temperaturá / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Övre temperaturgräns / Sıcaklık üst sınırı / Максимальна температура / 温度上限

Keep dry / Παзете сухо / Skladujte v suchém prostředí / Opbeavars tørr / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivais / Conservar au sec / Držati na suhom / Százaz helyen tartandó / Tenere all'asciutto / Құрғақ күйінде ұста / 건조 상태 유지 / Laikykite sausiai / Uzglabāt sausu / Droog houden / Holds tørr / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / Ne dopuskać popadnięcia wplgni / Uchovávať v suchu / Držite na suvom mestu / Förvaras tørr / Kuru bir şekilde muhafaza edin / Берегти від вологи / 请保持干燥

Collection time / Време на събиране / Čas odběru / Opsamlingstidspunkt / Entnahmeuhrzeit / Ωρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинақы уақыты / 수집 시간 / Paemimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de colheita / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingstid / Toplama zamani / Час забору / 采集时间

Peel / Обелете / Olevfete zede / Abn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Υστίρη κάθαίνι αλπι таста / 벗기기 / Plešti čia / Atfimtē / Schillen / Trekk av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odtrhnite / Oluštiti / Dra isār / Айрма / Відкелити / 撕下

Perforation / Перфорация / Perforce / Perforing / Διήτρηση / Perforación / Perforations / Perforacia / Perforalás / Perforazione / Tecik teicy / 점취선 / Perforacia / Perforacia / Perforate / Perforacia / Perforação / Perforare / Perforacion / Perforacia / Perforasyon / Перфорация / 穿孔

Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nerouživejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inihal beschädigter Packungnicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά. / No use si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Eter niet bezuylman bokan, palidandalaba / 패키지가 손상된 경우 사용 금지 / Jei pakuoje pažeista, nenaudoti / Nelietoti, ja terpakujioms bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Nào usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nerouživajte, ak je obal poškozený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包裝破損，請勿使用

Keep away from heat / Παзете от топлина / Nevystavujte přilísnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valguses / Protéger de la chaleur / Držati dalje od izvora topline / Övja a melegtől / Tenere lontano dal calore / Салқын жерде сақта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no karstuma / Beschermer tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródeł ciepła / Manter ao abrigo do calor / A se feri de căldură / Не награвяте / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Берегти від дії тена / 请远离热源

Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Lögata / Découper / Rezi / Vágja ki / Tagliare / Kecişia / 잘라내기 / Kirpti / Nogriet / Knippen / Kutt / Odciąć / Cortar / Decupajti / Отрезать / Odstřihněte / Iseći / Klipp / Kesme / Pozpaziti / 剪下

Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumisküüraev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинақаның тізбектелген / 수집 날짜 / Paemimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de colheita / Data colectării / Дата сбора / Datum odboru / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期

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Rx Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

**Australian Sponsor:**

Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113  
Australia

ATCC is a trademark of the American Type Culture Collection.

© 2017 BD. BD, the BD Logo and BD BBL are trademarks of Becton, Dickinson and Company.