

BD Nexiva™

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EN • Before needle removal Figure B Before needle removal Figure C Shielded needle

Indications for Use: The Nexiva intravascular catheter is inserted into a patient's vasculature to sample blood, monitor blood pressure, or administer fluids intravascularly. The needle-shielding feature is notched near the luer access port, aid in the prevention of needle-stick injuries. Blood is contained within the device during the catheter insertion process aided in the prevention of blood exposure. This catheter may be used for any patient population with consideration given to adequacy of vascular anatomy and appropriateness of procedure.

The 18-22 gauge Nexiva catheters are suitable for use with power injectors rated for a maximum of 300 psi when the luer access port(s) is removed and a direct connection is made.

Device Description: The BD Nexiva Closed IV Catheter System is designed to minimize blood exposure. It includes a passive needle-shielding mechanism designed to reduce accidental needlestick injury. The closed system is designed to keep blood contained within the device throughout the insertion process, which may prevent potential exposure for clinicians and patients. The system consists of a radiopaque Valium™ tube with an enhanced stabilization platform, a luer adapter, a luer connector, and a luer adapter. The system consists of a radiopaque Valium™ tube with an enhanced stabilization platform, a luer adapter, a luer connector, a clamp, a vent plug and a luer connector. The luer connector includes two needleless connectors, with one attached to the luer connector and the other included in the package. The 20-gauge is capable of withstanding high pressure injection procedures. The stabilization platform and luer adapter are color-coded.

RT ONLY CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

The 18-22 gauge Nexiva catheters are suitable for use with power injectors rated for a maximum of 300 psi when the luer access port(s) is removed and a direct connection is made.

Device Description: The BD Nexiva Closed IV Catheter System is designed to minimize blood exposure. It includes a passive needle-shielding mechanism designed to reduce accidental needlestick injury. The closed system is designed to keep blood contained within the device throughout the insertion process, which may prevent potential exposure for clinicians and patients. The system consists of a radiopaque Valium™ tube with an enhanced stabilization platform, a luer adapter, a luer connector, and a luer adapter. The system consists of a radiopaque Valium™ tube with an enhanced stabilization platform, a luer adapter, a luer connector, a clamp, a vent plug and a luer connector. The luer connector includes two needleless connectors, with one attached to the luer connector and the other included in the package. The 20-gauge is capable of withstanding high pressure injection procedures. The stabilization platform and luer adapter are color-coded.

RT ONLY CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

The 18-22 gauge Nexiva catheters are suitable for use with power injectors rated for a maximum of 300 psi when the luer access port(s) is removed and a direct connection is made.

GENERAL GUIDELINES

- For proper use, clinicians must be familiar with and trained in the use of the BD Nexiva Closed IV Catheter System.
- Observe Standard Precautions on ALL patients.
- Aseptic technique, proper skin preparation, and continued protection of the insertion site are essential.
- Non-pyrogenic. Sterile unless package has been opened or damaged.
- This product is DEHP free. This product does not contain natural rubber latex.
- Volume of 18 and 20 gauge is 0.5 mL. Volume of 22 and 24 gauge is 0.3 mL.

CAUTIONS

- Do not bend the needs when using the catheter system.
- Never insert a needle into the catheter septum (8).
- Do not use scissors at or near the insertion site.
- Do not attempt to override or defeat the needle shielding mechanism.
- Report medication immediately and follow established protocol.
- Immediately dispose of any needles that fail to shield, keeping the needlepoint away from the body and fingers at all times.
- Exposure to blood, either through percutaneous puncture with a contaminated needle or via mucous membranes, may lead to serious illness such as hepatitis, HIV (AIDS), or other infectious diseases.
- Immediately dispose of any needles that fail to shield, keeping the needlepoint away from the body and fingers at all times.
- Exposure to blood, either through percutaneous puncture with a contaminated needle or via mucous membranes, may lead to serious illness such as hepatitis, HIV (AIDS), or other infectious diseases.
- 18-22 gauge catheter systems are suitable for use with power injectors set to a maximum pressure of 300 psi and within maximum flow rate recommendations (see tables), when the access port is removed and a direct connection is made.
- The use of the catheter system must be assumed immediately before power injection.
- Measures should be taken to avoid kinking or obstructing the catheter system during power injection to avoid product failure.
- Re-use may lead to infection or other illness/injury.
- It is recommended that the BD Nexiva device be changed according to facility policy, CDC (USA) guidelines, or if the integrity of the device has been compromised.

WARNINGS

- 24 gauge catheter systems should not be used with power injectors.

CATHETER INSERTION

- Inspect the catheter system. Check that the clamp (5) is not engaged and that the vent plug (1) is secure.
- Holding the colored stabilization platform (2), twist and remove the needle cover.
- Hold the catheter (4) at the finger grips (4) and pull back approximately 1/8" (3mm) on the finger grips (4).
- Push the finger grips (4) back to their original position, making sure that the colored stabilization platform (2) and finger grips (4) are snugly together.
- Access the vessel.
- Observe immediate blood return along the catheter tubing (1). NOTE: Blood return will slow at the catheter/extension tube junction before continuing up the extension tube (1).
- Never insert or advance the entire catheter and needle until slightly to ensure that the catheter tip is within the vessel.
- Advance the catheter off the needle and into the vessel using the push tab (3).
- With the opposite hand, gently press the stabilization platform (2) to the skin to stabilize the catheter.
- Ensure that your fingers are on the finger grips (4), and pull the needle completely out of the catheter until the push tab component (3) releases from stabilization platform (2). DO NOT HOLD UNTIL THE PUSH TAB COMPONENT (3) OF THE DEVICE AS THIS WILL PREVENT THE RELEASE OF THE NEEDLE SHIELD.
- Remove the vent plug (1) and attach administration tubing or access port.
- Secure and apply the sterile dressing.
- Engage the clamp (5).
- Discard the shielded needle in a puncture resistant, leak proof sharps container.
- Release the clamp (5) to flush or infuse.

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BD Medical

9450 South State Street
Sandy, Utah 84070 USA.

For information not listed here, see appropriate Buy Specification.

Title: BD Nexiva Single Port IFU

Artwork Drawing Number: D16383

Single Unit Dimensions: 22" x 17" sheet folded to 5.5" x 4.25"

Print Colors: PMS 2755 Blue

Notes: Finished folded size is 5.5" x 4.25" with the BD Logo panel facing up.

Material: 40# White Offset

Revision History

Rev	Description of Changes	ECO No.
1	Initial release for Nexiva single port consolidation.	ECO100629
2	Add "Indications Use" and "Device Description" sections. Add Rx ONLY statement. Add last line in Cautions section and Warnings section. Updated ISO symbols legend. Updated trademark by-line date. Updated address to RA S: drive. Updated SAP and document numbers. Changed page size to 22" x 17" and fold size to 5.5" x 4.25".	ECO151232

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Brands are trademarks of their respective owners.