

	Wavelength	Radiant Power
RED	660 \pm 2nm	1.8mW
IR	940 \pm 10nm	2.0mW

5. Power Requirements

Two AAA alkaline batteries
Power consumption: Less than 30mA
Battery Life: Two AAA 1.5V, 600mAh alkaline batteries could be continuously operated as long as 30 hours.

It is equipped with a function switch, through which the oximeter can be powered off if no finger is present in the oximeter for longer than eight seconds.

6. Outline Dimensions


Length: 49~51mm
Width: 28~30mm
Height: 28~30mm
Weight: 30~50g (including two AAA batteries)

7. Environment Requirements

Operation Temperature: 5 ~ 40°C
Storage Temperature: -20 ~ +55°C
Ambient Humidity: ≤80% no condensation in operation
≤93% no condensation in storage

8. Equipment Response Time

As shown in the following figure.
Response time of slower average is 12.4s.



9. Classification

According to the type of protection against electric shock: INTERNALLY POWERED EQUIPMENT
According to the degree of protection against electric shock: TYPE BF APPLIED PART
According to the degree of protection against ingress of water: IPX1
According to the mode of operation: CONTINUOUS OPERATION

Declaration

Guidance and Manufacturer's declaration – electromagnetic emissions – For all EQUIPMENT and SYSTEMS

Guidance and Manufacturer's declaration – electromagnetic emission		
The Pediatric Finger Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic Environment – guidance
RF emissions CISPR 11	Group 1	The Pediatric Finger Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Guidance and Manufacturer's declaration – electromagnetic immunity – For all EQUIPMENT and SYSTEMS

Guidance and Manufacturer's declaration – electromagnetic immunity			
The Pediatric Finger Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance Level	Electromagnetic Environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	+/- 6kV contact +/- 8kV air	+/- 6kV contact +/- 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and Manufacturer's declaration – electromagnetic immunity
For all EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and Manufacturer's declaration - electromagnetic immunity

The Pediatric Finger Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance Level	Electromagnetic Environment – guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Pediatric Finger Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{E_r} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz} \quad d = \left[\frac{7}{E_r} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and Manufacturer's declaration – electromagnetic immunity
For all EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING – CONTINUED

- Field strengths from fixed transmitters, such as base station for radio (cellular/cordless) telephones and land mobile radio amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Pediatric Finger Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measurements may be necessary, such as reorienting or relocating the Pediatric Finger Pulse Oximeter.
- Over the frequency range 150 kHz to 80 MHz, fields strengths should be less than 3 V/m

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEMS – For all EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and Pediatric Finger Pulse Oximeter		
The Pediatric Finger Pulse Oximeter is intended for use in electromagnetic environments in which radiated RF disturbances are controlled. The customer or the user of the Pediatric Finger Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Pediatric Finger Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	$d = \left[\frac{3.5}{E_r} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$	$d = \left[\frac{7}{E_r} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
0.01	0.1167	0.2334
0.1	0.3689	0.7378
1	1.1667	2.3334
10	3.6893	7.3786
100	11.6667	23.3334
For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.		
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

Possible Problems and Solutions

Problems	Possible reason	Solution
SpO ₂ or PR can not be shown normally	1. Finger is not inserted correctly. 2. Patient's SpO ₂ value is too low to be measured.	1. Retry by inserting the finger. 2. There is excessive illumination. 3. Try some more times. If you can make sure no problem is existing in the product, please go to a hospital in a timely manner for exact diagnosis.
SpO ₂ or PR is unstable	1. Finger might not be inserted deep enough. 2. Excessive patient movement.	1. Retry by inserting the finger. 2. Be calm.
The oximeter cannot be powered on	1. No battery or low battery power. 2. Batteries might be installed incorrectly. 3. The oximeter might be damaged.	1. Please replace the batteries. 2. Please reinstall the batteries. 3. Please contact your local customer service center.
Indication lamps are suddenly off	1. The product is automatically powered off when no signal is detected for longer than eight seconds. 2. The battery power is too low to work.	1. This is normal. 2. Please replace the batteries.
"Error 3" or "Error 4" is displayed on screen	1. "Error 3" means the red emission LED is damaged. 2. "Error 4" means the infra-red emission LED is damaged.	1. Check the red emission LED. 2. Check the infra-red emission LED.
"Error 6" is displayed on screen	"Error 6" means the screen has failed.	Change the screen.
"Error 7 " is displayed on screen	"Error 7" means the emission LED or reception diode is damaged.	Check the emission LED and reception diode.

Symbol Definitions

Symbol	Definition	Symbol	Definition
	Type BF applied part		No SpO ₂ Alarm
	Consult accompanying documents		Storage temperature and relative humidity
	Protected against dripping water		Serial No.
	Oxygen saturation		Manufacturer's information
	Pulse rate (BPM)		Date of Manufacture
	Low power indication		

Please complete the attached warranty card and return for our record. Refer to the Dynarex website for complete warranty details.

Note: The illustrations used in this manual may differ slightly from the appearance of the actual product.