DASELNE® EVALUATION INSTRUMENTS


BASELINE® BODY LEVEL / SCOLIOSIS METER



Baseline® body level

Use to determine whether body parts are properly aligned. Ideal for Scoliosis screening.

How to use: Hold the Baseline® body level with gauge facing you. Extend the moveable arms and place them on two sites. Read the amount of tilt (in degrees) from the scale.



Baseline® scoliosis meter

Measurements ontained with scoliosis meter can detect scoliosis and abnormal anteropostero curve. Use to measure unstable lumbosacral, cervical and thoracic curves. Scale in cm., in. and degrees.

How to use: Hold the Baseline® scoliosis meter in one hand and loosen the screw on the level with the other hand. Allow the scoliosis meter to turn in your hand, then place the level on the side fastening it with the screw. The lighter markings on the rods are two (2) millimeters apart and the heavier markings are ten (10) millimeters.



Manufacturer and Master Distributor of Physical Therapy and Rehabilitation Products

Fabrication Enterprises, Inc. PO Box 1500, White Plains, New York 10602

tel: 800-431-2830 914-345-9300 fax: 800-634-5370 914-345-9800 info@FabricationEnterprises.com

TABLE OF CONTENTS

Introduction	page 3
Adam's Sign	page 3
How to Measure:	
Shoulder tilt in degrees	page 3
Pelvic tilt in degrees	page 3
Shoulder tilt in cm/mm	page 3
Thoracic deviation in degrees and cm/mm	page 4
Normal antero-postero cervical curve - with spine	page 4
Normal antero-postero thoracic curve - with spine	page 5
Detection of unstable lumbo-sacral angle with spine	page 5
School worksheets	page 6
Physician's worksheets	page 7
More Baseline® evaluation instruments	back cover


INTRODUCTION:

The Baseline® body level and the Baseline® scoliosis meter are used to measure body alignment. The Baseline® body level is excellent for measuring scoliosis (lateral curvature of the spine). The Baseline® scoliosis meter can detect an abnormal anteropostero curve, and measures the unstable lumbosacral, cervical and thoracic curves.

USER NOTE Please consult your appropriate training material for more information on testing protocols, procedures and methods. This pamphlet gives a summary of the types of tests that can be performed with the scolioiss meter and body level.

ADAM'S SIGN is a standardized neurological / orthopedic procedure and testing method. A patient with scoliosis, when bending over, will have no straightening of the curve. This result is "POSITIVE." A straightening of the curve indicates a "NEGATIVE" result.


HOW TO MEASURE:



Shoulder tilt in degrees

PELVIC TILT IN DEGREES


The blocks are released and the rods are placed over the crest of the ilium. The blocks are then pressed firmly toward the midline. Read the angle from the level. If the gauge reads over 2%, the result is listed as positive.



Shoulder tilt in centimeters

SHOULDER TILT IN DEGREES

With the patient in a standing position, place the rods on the acromioclavicular articulation. Read the angle from the level gauge. If the gauge reads over 2%, the reading is listed as positive. NOTE: to read the level gauge and obtain the degree measurement, look for the ball on the gauge.




Pelvic tilt in degrees

SHOULDER TILT IN DISTANCE

Placing the ends of the rods on the acromioclavicular articulation, you can get the degree measurement. If you release the rod on the higher side and move the bar down until "0" is indicated on the gauge, you will be able to calculate the distance by counting the markings above the block on the high side. This indicates the distance deviation.

THORACIC DEVIATION IN DEGREES

With the patient bending forward and with the rods locked, place the instrument over the thoracic rib cage with the ends of the rods touching the apex (highest point) of the curve. This point will vary on each individual. In this position you will be able to acquire the degree readings on the gauge by loosening the rod on the high side and pushing the bar down until the gauge reads "0". You will be able to get the distance reading by counting the lines above the block on the high side.



Thoracic deviation in ribs measured in degrees


NOTE:

As above, with the patient in a 90° flexed position, rotation of the lumbar vertebra is measured by bringing the instrument over the lumbar vertebra to the area of the transverse process. The reading will be determined, as demonstrated with thoracic deviation, in degrees, and distance (cm/mm).

NOTE: In order to acquire the following readings, it is necessary to have the Baseline® scoliosis meter in a vertical position with the gauge facing the examiner and on 0°.

CERVICAL CURVE


The top block is positioned so that the rod touches the spinous of C4 (approximately 3-1/4" below the occipital or base of the brain). It is locked in place. The bottom block is positioned so that it touches the spinous process of C7 (most prominent spinous process in the cervical area palpating from the base of the skull downward). It is locked in place. Obtain the distance reading by counting the hash marks. A measurement over 3 cm. indicates lordosis above the normal range and a measurement below 1 cm indicates kyphosis.



Cervical curve C4-C7 Normal = 1-3 cm.

THORACIC CURVE


The top block is positioned so that the rod touches the spinous of C7 (most prominent spinous process in the cervical area palpating from the base of the skull downward). It is locked in place. The bottom block is positioned to allow the rod to touch the spinous process of the vertebra which is the apex of the thoracic curve. The bottom rod is loosened and the reading distance reading is possible. The normal range is from 1cm to 4 cm. A decrease in the number indicates a hypo-kyphotic area. An increase in the number indicates a more severe kyphosis, over it's normal 40°.



Thoracic curve C7 - apex of curve Normal = 1-4 cm.

ANTEROPOSTERO ANGULATION L5 AND APEX OF SACRUM

The base of the spine is the trigger point of spinal disorders. It is very important that this area undergo the closest scrutiny. Due to all the ramifications in the various severe low back conditions, the top block is positioned so that the rod touches the spinous process of L5 (place your index fingers over the crest of the ilium and then place your thumbs together to locate L5). The bottom block is positioned so that the rod will touch the most prominent tubercle at the apex of the sacrum. The distance reading can then be established. If there is an unstable lumbo-sacral angle, the reading would be greater than 3 cm. If the reading is less than 1 cm, a kyphotic spine would exist. The normal sacral angle is 42°, which would correspond to the 1 to 3 cm normal range.



anteropostero angulation of L5 and apex of sacrum measured in centimeters/millimeters Normal = 1-3 cm.

NOTE: The Baseline® scoliosis meter can measure anterior humping by placing rods on the apex of the anterior rib cage to obtain a distance and degree reading.

BASELINE® SCOLIOSIS METER - SCHOOL WORKSHEETS

STUDENT NAME		ADDRESS		SH Degrees	SHOULDER TILT Degrees Centimeters/Millimeters (in)	PELVIC TILT Degrees	LUMBAR ROTATIONAL DEVIATION Centimeters/Millimeters (in)
	Date of Examination	Examination Results	Physician Follow-up	5th			
5th Grade/Age				6th			
				7th			
Remarks:				8th			
6th Grade/Age				9th			
				10th			
Remarks:				11th			
7th Grade/Age				12th			
				THORA	CIC DEVIATION IN POSTERIC	OR RIB CAGE	ANTERIO-POSTERO DEVIATION
Remarks:				Degre	Degrees Centimeters/Millimeters (in)	s/Millimeters (in)	OF L5-APEX OF SACRUM Centimeters/Millimeters (in)
8th Grade/Age				5th			Normal = 1.3 cm (5.11 in)
				6th			
Remarks:				412			
9th Grade/Age				8th			
				र के			
Remarks:				10 T			
10th Grade/Age				1011			
				TITU			
				12th			
Remarks:							
11th Grade/Age					SPINAL ANALYSIS (Centimeters/N	SPINAL ANALYSIS OF ANTERIO-POSTERO CURVES Centimeters/Millimeters (in)	ERO CURVES
					Cervical Curve Normal = 1.3 cm (5.11 in)	Norn	Thoracic Curve Normal = 1.3 cm (5.11 in)
Remarks:							
12th Grade/Age				6th			
				7th			
Remarks:				8th			
				9th			
				10th			
				11th			
				12th			
							School Nurse

BASELINE® SCOLIOSIS METER - PHYSICIANS WORKSHEETS

LUMBAR ROTATIONAL DEVIATION Centimeters/Millimeters (in)					ANTERIO-POSTERO DEVIATION	OF L5-APEX OF SACRUM Centimeters/Millimeters (in)	Normal = 1.3 cm (5.11 in)					ERO CURVES	Thoracic Curve Normal = 1.3 cm (5.11 in)					Dhusioinn's Signature
PELVIC TILT Degrees					OR RIB CAGE	s/Millimeters (in)						OF ANTERIO-POSTI	Norm					
SHOULDER TILT Degrees Centimeters/Millimeters (in)					SIC DEVIATION IN POSTERIC	Degrees Centimeters/Millimeters (in)						SPINAL ANALYSIS OF ANTERIO-POSTERO CURVES Centimeters/Millimeters (in)	Cervical Curve Normal = 1.3 cm (5.11 in)					
SHOI Degrees					THORAC	Degrees							Nom					
AGE	Follow-up Results																	
ADDRESS	Examination Results																	
PATIENT NAME	Date of Eaxamination																	



www.BaselineProducts.net

The Most Extensive Instrument Range for Measuring Progress in Physical Therapy

- Baseline Measure strength, range-of-motion, and more
- Baseline Lightweight and portable
- Baseline Inexpensive and cost effective
- Baseline Accepted by physical and occupational therapists for over 25 years



wrist dynamometer



back-leg-chest dynamometer



pinch gauges



bubble inclinometer



hand evaluation set



taylor hammers



goniometers



wrist inclinometers

For complete line of Baseline products, see pages 56-87 in the 2008 FEI catalog

or view online at www.FabricationEnterprises.com