HemoCue® WBC

Operating Manual

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GB/US

Table of Content

HemoCue® WBC system	5
Components	6
Start Up	8
Set-up	12
Audio signal	12
Measuring	14
Capillary blood	14
Measuring	22
Control material or venous blood	22
Set-up	24
Printer function	24
Maintenance	26
Troubleshooting Guide	28
Specifications	44

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HemoCue® WBC-systeem

GB/US

Thank you for choosing the HemoCue WBC system. The system is indicated for use for quantitative determination of white blood cell (WBC) count in capillary or venous whole blood. The HemoCue WBC system is for *In Vitro* Diagnostic use only. The HemoCue WBC Analyzer is only to be used with HemoCue WBC Microcuvettes. The HemoCue WBC system is indicated for use in clinical laboratories and for point-of-care settings.

All system components are designed and manufactured to provide maximum safety. Any other use of the system than indicated may impair the safety.

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Components
Bestandteile
Éléments
Onderdelen

GB/US

1. HemoCue WBC Analyzer*
2. AC adapter or
3. Six type AA batteries**
4. HemoCue WBC Microcuvettes**
5. HemoCue WBC Operating Manual, HemoCue WBC Quick Reference Guide and Instruction CD
6. HemoCue Cleaner

Open the carton and lift out the analyzer and accessories and place them on a stable surface (non-vibrating).

* Do not open the cover of the analyzer.

** Not included.

For information about HemoCue WBC Microcuvettes, please contact your HemoCue distributor.

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
GB/US

Only use adapters, as listed under Specifications, in the Adapters section.

1. If AC power is available, plug the supplied power adapter into the power inlet at the back of the Analyzer.
2. If no power is available, insert six type AA batteries, 1.5 V. To open the battery lid on the left side of the Analyzer, press on the ribbed marking engraved on the battery lid and slide it backwards until it opens. Gently take out the battery holder. Place the batteries in the battery holder observing the indication of polarity. Put the battery holder back in the battery compartment and close the lid. Consult local environmental authorities for proper disposal of batteries.
3. Pull the cuvette moving arm out to the loading position.

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


GB/US

4. Press and hold the button until the display is activated (all symbols appear on the display). The Analyzer performs a self test, and after approximately 10 seconds the display will show three flashing dashes and the HemoCue symbol. This indicates that the Analyzer is ready for use.
5. If the cuvette moving arm is in the measuring position, the most recent result will be displayed.
6. To turn the Analyzer off, press and hold the button until the display reads OFF and then goes blank. If the Analyzer is operating on battery power, but not being used, it will automatically turn off after approximately five minutes and if operating on AC power, it will turn off after two hours.

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
Audio signal
Tonsignal
Signal sonore
Geluidssignaal

GB/US

If the audible signal is activated, a signal will be heard when the measurement is completed, or if an error code is displayed. To activate or deactivate the audible signal, follow the steps below.

1. Make sure that the analyzer is turned off. Press the button for approximately 10 seconds.
2. The display now shows a flashing bell symbol.
3. Press the button rapidly to change between signal on and off.

When the setting is completed, press the button for approximately 3 seconds until the bell stops flashing. The analyzer will return to the ready mode.

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
Capillary blood
Kapillarblut
Sang capillaire
Capillair bloed

GB/US

Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal. Always wear protective gloves when handling blood specimens. The microcuvette is for single use only.

1. To perform a test, the cuvette moving arm should be in its loading position. The display will show three flashing dashes and the HemoCue symbol. Take a HemoCue WBC Microcuvette from the vial.
2. Make sure the patient's hand is warm and relaxed. Use only the middle or ring finger for sampling. Avoid fingers with rings on.
3. Clean fingertip with disinfectant and allow to dry completely or wipe off with a dry, lint-free wipe.
4. Using your thumb, lightly press the finger from the top of the knuckle towards the tip.

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



GB/US

5. Sample at the side of the fingertip.
6. While applying light pressure towards the fingertip, puncture the finger using a lancet*.
7. Wipe away the first two or three drops of blood.
8. Re-apply light pressure towards the fingertip until another drop of blood appears.

* Spring loaded lancets with a puncture depth of at least 2 mm are preferred to produce a sufficient blood flow.

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


GB/US

- When the blood drop is large enough, fill the microcuvette in one continuous process. Do NOT refill! NOTE: Make sure that the microcuvette is filled from the tip, placed at about a 45 degree angle towards the blood drop according to the picture on page 18.
- Wipe off excess blood from the outside of the microcuvette with a clean, lint-free wipe. Do not touch the open end of the microcuvette.
- Look for air bubbles in the filled microcuvette. If present, discard the microcuvette and fill a new microcuvette from a new drop of blood. Small bubbles around the edge can be ignored. NOTE: Make sure that the microcuvette is filled according to picture 9 on page 18 since an improper filling angle might cause air bubbles to be introduced.

NOTE: If a second sample is to be taken, it is important that this is done after the measurement of the first sample is complete. Wipe away the remains of the drop of blood and fill the second microcuvette from a new drop of blood as per steps 7–11 above.

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



GB/US

12. Place the filled microcuvette in the cuvette holder within 40 seconds after filling.
13. Gently push the cuvette moving arm towards the measuring position. It will automatically slide to the measuring position and the measurement starts.
14. During the measurement “”, three fixed dashes and the HemoCue symbol will be shown.
15. After approximately 3 minutes, the WBC value is displayed. The result will remain on the display as long as the cuvette moving arm is in the measuring position. Do not remeasure the filled microcuvette.
16. Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal.

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
GB/US

Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal. Always wear protective gloves when handling blood specimens. The microcuvette is for single use only.

1. If the venous blood samples has been stored in a refrigerator, allow it to warm up to room temperature (15–35 °C, 59–95 °F) before mixing. Mix the venous sample tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 10–20 times by hand. For control materials always follow instructions for use provided by the manufacturer.
2. Place a drop of blood or control material onto a hydrophobic plastic or glass slide using a pipette or other suitable transfer device.
3. Fill the microcuvette in one continuous process. Do NOT refill! NOTE: Make sure that the microcuvette is filled from the tip, placed at about a 45 degree angle towards the blood drop according to the picture on page 22.

Perform the analysis as per steps 10–16 on page 19–21.

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
Maintenance
Wartung
Maintenance
Onderhoud

GB/US

Before performing maintenance, make sure that the analyzer is turned off and that the display is blank.

1. The cuvette holder should be cleaned after each day of use. Pull the cuvette moving arm to its loading position.
2. Remove the cuvette holder by lifting it straight up.
3. Clean the cuvette holder with alcohol (20–70 %) or mild detergent. NOTE: The cuvette holder cannot be autoclaved.
4. If the optical parts becomes dirty, an error code will be displayed. To clean the optical parts, push a HemoCue Cleaner* into the opening of the cuvette holder. Move from side to side 5–10 times. If the cleaner is stained, repeat with a new cleaner. Make sure that the cleaner reaches the parts according to 4a.
5. Wait 15 minutes before replacing the cuvette holder. Make sure that the cuvette moving arm is in its loading position before replacing the cuvette holder. The cover may be cleaned with alcohol (20–70 %) or a mild detergent.

* One HemoCue Cleaner is included with the Analyzer. Contact your local distributor to order more.

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Troubleshooting Guide

If you are unable to resolve the problem by following this Troubleshooting Guide, please contact your local HemoCue distributor or HemoCue AB. The Analyzer should be cleaned as described in the Maintenance section prior to service or disposal. Consult local environmental authorities for proper disposal. The Analyzer has no serviceable parts. NOTE! Do not open the cover of the Analyzer.

Symptom	Explanation	Action
Err01	A portion of the image area is unable to be analyzed. 1. Due to measurement error. 2. Due to abnormal sample.	1. Take a new microcuvette and repeat the measurement, as described in the Measuring section. 2. NOTE! The blood sample should be verified with a suitable laboratory method and be questioned as to the pathological condition of the patient.
Err02	Uneven spatial distribution of detected cells.	Take a new microcuvette and repeat the measurement, as described in the Measuring section.
Err03	Image, or part of the image area is detected as out-of-focus.	Take a new microcuvette and repeat the measurement, as described in the Measuring section.
Err30	1. Optical parts dirty. 2. Optical parts wet after cleaning.	1a. Turn off the analyzer and clean the optical parts as described in the Maintenance section. 1b. If the problem continues, the analyzer needs service. Contact your local distributor. 2. Wait 15 minutes before turning on the analyzer after cleaning to make sure that the optical parts are dry.
Err33	Empty microcuvette, not filled with sample.	Take a new microcuvette and repeat the measurement, as described in the Measuring section and make sure that the microcuvette is filled with sample.
Err34	Stray light detected.	a. Turn off the analyzer and make sure the analyzer is not exposed to any bright light sources. b. If the problem continues, the analyzer needs service. Contact your local distributor.
Err35	The battery power is too low.	Turn off the analyzer. a. Replace the batteries, six type AA, as described in the Start Up section. b. Use the power adapter as described in the Start Up section.

Symptom	Explanation	Action
Err60	General hardware error.	<p>a. Turn off the analyzer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement, as described in the Measuring section.</p> <p>b. If the problem continues, the analyzer needs service. Contact your local distributor.</p>
Err61	Self test error during start up of analyzer.	<p>a. Turn off the analyzer and turn it on again after 30 seconds. Take a microcuvette and perform the measurement, as described in the Measuring section.</p> <p>b. If the problem continues, the analyzer needs service. Contact your local distributor.</p>
Err62	Blanking test failed for other reasons than Err30, Err34.	<p>a. Turn off the analyzer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement, as described in the Measuring section.</p> <p>b. If the problem continues, the analyzer needs service. Contact your local distributor.</p>
WbC	Empty cuvette holder.	Take a microcuvette and perform the measurement, as described in the Measuring section.
LLL	<p>1. Measured value is below $0.3 \times 10^9/L$ ($300/mm^3$, $300/\mu L$).</p> <p>2. Empty microcuvette, not filled with sample.</p>	<p>1. NOTE! Results exceeding the LLL limit should be verified with a suitable laboratory method and be questioned as to the pathological condition of the patient.</p> <p>2. Take a new microcuvette and repeat the measurement, as described in the Measuring section and make sure that the microcuvette is filled with sample.</p>
HHH	Measured value exceeds $30.0 \times 10^9/L$ ($30000/mm^3$, $30000/\mu L$).	NOTE! Results exceeding the HHH limit should be verified with a suitable laboratory method and be questioned as to the pathological condition of the patient.

Symptom	Explanation	Action
No characters on the display.	<ol style="list-style-type: none"> 1. Power is not received. 2. If on battery power, the batteries need to be replaced. 3. The display is out of order. 	<ol style="list-style-type: none"> 1a. Check that the power adapter is properly connected to the analyzer and the AC power supply, as described in the Start Up section. 1b. Check that the cable is not damaged. 2. Replace the batteries, six type AA, as described in the Start Up section. 3. The analyzer needs service. Contact your local distributor.
The display gives erroneous characters.	<ol style="list-style-type: none"> 1. The display is out of order. 2. The microprocessor is out of order. 	<ol style="list-style-type: none"> 1, 2. The analyzer needs service. Contact your local distributor.
The display shows "FIR".	General analyzer software error.	<ol style="list-style-type: none"> a. Remove and replace all cables and/or batteries, and restart, as described in the Start Up section. b. The analyzer needs service. Contact your local distributor.
The display shows "▣▣▣".	<ol style="list-style-type: none"> 1. The batteries need to be replaced. 2. If on AC power, the power adapter or the circuit board is out of order. 	<ol style="list-style-type: none"> 1. Turn off the analyzer and replace the batteries, six type AA, as described in the Start Up section. 2a. Check that the correct power adapter is used and properly connected and working, as described in the Start Up section. 2b. The analyzer needs service. Contact your local distributor.
The display does not switch from "Hb" and "WbC" to three flashing dashes and "Hb" (ready for measuring).	The cuvette holder sensor is out of order.	The analyzer needs service. Contact your local distributor.
The cuvette holder is not moving into the correct position.	The magnet in the cuvette holder is missing.	The analyzer needs service. Contact your local distributor.
Measurement on patient samples are higher or lower than anticipated.	<ol style="list-style-type: none"> 1. Improper sampling technique. 2. The microcuvettes are beyond their expiration date, damaged or have been improperly stored. 	<ol style="list-style-type: none"> 1. Take a new microcuvette and repeat the measurement, as described in the Measuring section. 2. Check the expiration date and the storage conditions of the microcuvettes.

Symptom	Explanation	Action
<p>Measurement on control materials are out of range, either too high or too low.</p>	<ol style="list-style-type: none"> 1. Improper sampling technique. 2. The microcuvettes are beyond their expiration date, damaged or have been improperly stored. 3. The control material is beyond its expiration date or has been improperly stored. 4. The control material has not been mixed properly and/or is not at room temperature. 5. The control material is not suitable for use with the HemoCue WBC system. 	<ol style="list-style-type: none"> 1. Take a new microcuvette and repeat the measurement, as described in the Measuring section. 2. Check the expiration date and the storage conditions of the microcuvettes. 3. Check the expiration date and the storage conditions of the control material. If the problem continues, contact the manufacturer of the control material. 4. Make sure that the control material is mixed properly and at room temperature. If the problem continues, contact the manufacturer of the control material. 5. Contact your local distributor for control material information.

Specifications

Intended Purpose/Intended Use

The HemoCue WBC system is indicated for use for quantitative determination of white blood cell (WBC) count in capillary or venous whole blood. The HemoCue WBC system is for *In Vitro* Diagnostic use only. The HemoCue WBC Analyzer is only to be used with HemoCue WBC Microcuvettes. The HemoCue WBC system is indicated for use in clinical laboratories and for point-of-care settings.

IVD Medical Device Directive

The HemoCue WBC system complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark.

Principles of the Method/Procedure

Principle of the method

A hemolyzing agent lyses the red cells in the microcuvette and a staining agent colors the white cells. An image is taken of the stained white cells and the number of cells is counted by image analysis in the analyzer.

Principle of the procedure

The microcuvette serves as a sample container and reaction chamber. The microcuvette is for single-use only. A blood sample of approximately 10 µL is drawn into the cavity by capillary action. The microcuvette is placed in the analyzer. The result is obtained within 3 minutes. The system was designed and developed to establish agreement with the manual light microscopy method for white blood cell count. The system is factory calibrated and needs no further calibration.

Warning and Precaution

The microcuvettes are for *In Vitro* Diagnostic use only. Always handle blood specimens with care as they may be infectious. Consult local environment authorities for proper disposal. Always wear protective gloves when handling blood specimens. The microcuvettes are for single-use only.

Storage and Handling

HemoCue WBC Microcuvettes

The Microcuvettes are to be stored at 15–35 °C (59–95 °F), <90 % non-condensing humidity. Once the seal of the vial is broken, the microcuvettes are stable for 3 months. An unopened vial of microcuvettes can be stored for a shorter period of time (4 weeks) outside the specified storage conditions down to 0 °C (32 °F) and up to 50 °C (122 °F), <90 % non-condensing humidity. Allow the microcuvettes to reach 15–35 °C (59–95 °F) before use. Use the microcuvettes prior to the expiration date that is printed on the package. Keep the vial properly closed. All unused microcuvettes should remain in the original package.

HemoCue WBC Analyzer

The analyzer can be stored in temperatures between 0–50 °C (32–122 °F), <90 % non-condensing humidity. Operating temperature is 15–35 °C (59–95 °F), <90 % non-condensing humidity. Allow the analyzer to reach ambient temperature before use.

Specimen Collection and Preparation

Capillary or venous whole blood may be used. EDTA anticoagulant may be used, preferably in solid form to avoid dilutional effects. Mix all specimen tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 10–20 times by hand. The specimen can be stored at room temperature (15–35 °C, 59–95 °F), or in a refrigerator (2–8 °C, 35–46 °F) for 48 hours. If the specimen has been stored in a refrigerator, it will be viscous and the blood should be allowed to warm up to room temperature before mixing.

Materials Required

- HemoCue WBC Analyzer
- HemoCue WBC Microcuvettes
- Lancet (capillary samples)

- Pipette or other transfer device (venous samples)
- Lint-free tissue (non-fraying)

Quality Control

The HemoCue WBC Analyzer has an internal quality control, the “self test”. Every time the analyzer is turned on, it will automatically verify the measurement performance. When passing the self test, the display will show the HemoCue symbol and three flashing dashes, indicating that the analyzer is ready to perform a measurement. An error code will be displayed if the self test fails. Another part of the built in self test (QC) is performed for each measurement, including both checks of the HemoCue WBC Analyzer but also several condition checks of the HemoCue WBC Microcuvette and the sample itself. The operator’s ability to handle the microcuvette and apply the sample correctly is also included in these self tests. No additional quality controls performed by the operator are required for verification of the system functionality. Note that local, state, or other accreditation agencies may require additional quality control testing. If additional quality control checks are required, R&D Systems R&D HC WBC Control may be used. Please refer to the R&D HC WBC Control, as well as the HemoCue WBC Operating Manual regarding handling and procedures. Please refer to local guidelines for recommended frequency of use.

Expected Values (Dacie and Lewis Practical Haematology)

Adults 4.0–10.0 $\times 10^9/L$
 Children 1 year 6.0–16.0 $\times 10^9/L$
 Children 2–6 years 5.0–15.0 $\times 10^9/L$
 Children 6–12 years 5.0–13.0 $\times 10^9/L$
 Infants 1 month 5.0–19.0 $\times 10^9/L$
 Infants 2 months 5.0–15.0 $\times 10^9/L$
 Infants 3–6 months 6.0–18.0 $\times 10^9/L$

The values above may vary due to a wide range of factors, such as sex, diurnal variations, exercise, physical stress or trauma, pregnancy, indigestion of food, and cigarette smoking.

Measuring Range

Displayed range: 0.3–30.0 $\times 10^9/L$ (300–30000/ mm^3 , 300–30000/ μL). Results above the measuring range will be displayed as HHH. Results below the measuring range will be displayed as LLL.

Limit of Detection

Functional Sensitivity is estimated as the mean concentration for a spiked sample whose CV is 20 %. The Functional Sensitivity has been determined to be 0.3 $\times 10^9/L$. (James O. Westgaard, Basic Method Validation) Limit of Blank is defined as the highest measurement result that indicates that the analyte is not present in the sample. Limit of Blank has been determined to be 0.06 $\times 10^9/L$ (CLSI Document EP17-A).

Limitations of the Method/Procedure

- The measurement should be made as soon as possible, but within 40 seconds after the blood has been drawn into the microcuvette.
- Do not remeasure the filled microcuvette.
- Mixing samples for an extended period may affect the result.
- Results above the measuring range will be displayed as HHH. Results below the measuring range will be displayed as LLL.
- Studies have shown that patient samples with >2 % nucleated red blood cells (NRBCs) may give falsely elevated white blood cell count.

Specific Performance Characteristics

Linearity

The method used in the HemoCue WBC system has according to the CLSI Document EP06-A been demonstrated to be linear between 0.3–30.0 x10⁹/L (300–30000/mm³, 300–30000/μL), within 0.2 x10⁹/L difference in 0.3–3.5 x10⁹/L and within 6 % difference in 3.6–30.0 x10⁹/L.

Within-run and Total Precision

Within-run and Total Precision was determined according to the CLSI Document EP05-A2. The results given below come from four batches of HemoCue WBC Microcuvettes, and five HemoCue WBC Analyzers. Commercially available controls at three different levels were used. The WBC count was measured in duplicate twice a day, morning and afternoon, during 20 consecutive days.

Level	N	x10 ⁹ /L	Within-run Precision CV %	Total Precision CV %
1	400	2.5	4.06	5.4
2	400	7.2	2.92	3.5
3	400	19.0	1.63	1.9

Additionally, a study has been performed to determine the precision using fresh blood samples. The study was performed on venous whole blood samples tested in duplicate.

x10 ⁹ /L	N	CV _{duplicates} %	
		Sysmex XS-1000i	HemoCue WBC
0.3–1.0	12	5.0	10.4
1.1–3.5	52	3.7	4.1
3.6–10.0	209	1.9	3.5
10.1–20.0	88	2.1	2.5
20.1–30.0	35	1.8	2.2
Total	396	2.4	3.2

Validation Studies

Results of the validation studies, based on the intended use of the system, are shown in Figure 1, 2, 3 and 4. Single replicate on the HemoCue WBC system and the mean of duplicate on the comparative method are presented in the plots.


Figure 1
Point-of-Care


Figure 2
Clinical Laboratory

When a capillary skin puncture is performed, several defense systems in the body are activated very quickly. These defense systems cause an increase in the number of WBCs in the blood closest to the wound, leading to greater differences in results from several samples taken from the same finger stick.

A study was performed, on the HemoCue WBC system in which possible influences of capillary sampling technique were included, see Figure 3 below.


Figure 3
Capillary Sampling Technique

An additional study has been performed at a clinical reference laboratory in the U.S. including a large number of blood samples in order to evaluate the HemoCue WBC system for expected demographic variations (age, gender, geographics) in the population, see Figure 4 below.


Figure 4
Demographic Sample Population

Interference Studies

Interference studies were performed according to CLSI Document EP07. The following substances have not been found to interfere with the HemoCue WBC system. The highest concentration or percentages tested is referred to in brackets.

Acetaminophen (20 mg/dL), ascorbic acid (6 mg/dL), ibuprofen (50 mg/dL), caffeine (308 µmol/L), creatinine/urea (30 mg/dL/500 mg/dL), salicylic acid (60 mg/dL), tetracycline (30 mg/dL), unconjugated bilirubin (342 µmol/L), conjugated bilirubin (380 µmol/L), dextran 40 (30 g/L), osmolality (167 mmol/L), HbCO (30 %), methemoglobin (50 %), thrombocytes (2000 x10⁹ g/L), lipemia (intralipid 3100 mg/L) which corresponds to triglycerides (approximately 10.5 mmol/L), hemoglobin (30–250 g/L), reticyocytes (258 x10⁹ L), HbAS, immature leukocytes and platelet clumping have not been found to interfere. pH values between 6.8–8.0 do not interfere with the system.

Technical Specifications

Dimensions: 185 x 133 x 120 mm (7.28 x 5.24 x 4.72 inches)

Weight: 600 g (1.32 pounds) (with 6 AA batteries installed)

Power adapter: CE marked

Only use adapters, as listed under Adapters.

Pollution degree: 2

Overvoltage category: II

Atmospheric pressure: 700 hPa to 1060 hPa.

Equipment not suitable for use in the presence of flammable mixtures.

The HemoCue WBC system has been tested for electrical safety and EMC according to the following standards:

- EN 60601-1:1990, A1:1993, A2:1995, A13:1996 Medical electrical equipment – Part 1: General requirements for Safety.
- EN 60601-1-1:2001 Medical electrical equipment – Part 1-1: General requirements for Safety – Collateral Standard: Safety requirements for medical electrical systems.
- UL 60601:2003 Medical electrical equipment – Part 1-1 General requirements for Safety.
- EN 61010-2-101:2002 Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-101: Particular requirements for *in vitro* diagnostic (IVD) medical equipment.
- UL 61010-1:2004 Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements.
- IEC/EN 60601-1-2+A1 Medical electrical equipment- Part 1: General requirements for Safety – Part 1-2: Collateral Standard: Electromagnetic compability – requirements and tests.

Recommended separation distance between Portable and mobile RF communications equipment and HemoCue WBC Analyzer

The HemoCue systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of HemoCue systems can help prevent electromagnetic interference by maintaining a minimum distance between portable and RF communications equipment (transmitters) and HemoCue systems as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at maximum output power not listed above, the recommended separation distances (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – Electromagnetic immunity

The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or user of the HemoCue systems should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds For explanation of U see NOTE 1	<5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the HemoCue systems image intensifier requires continued operation during power mains interruptions, it is recommended that the HemoCue systems image intensifier be powered from an uninterruptible power supply or a battery.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p> <p>See NOTE 2 and NOTE 3</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the HemoCue systems, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter</p> <p>Recommended separation distance</p> <p>$d=1.2\sqrt{P}$</p> <p>$d=1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> 

NOTE 1 U is the a.c. mains voltage prior to application of the test level.

NOTE 2 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 3 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, MA and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HemoCue systems are used exceeds the applicable RF compliance level above, the HemoCue systems should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the systems.
- b) Over the frequency range 150 KHz to 80 Mhz, field strength should be less than 3 V/m.

Technical specifications (EMC-RF)

Use only cables with the following specification:

USB shielded maximum 2 m

Serial shielded maximum 1.5 m

Guidance and manufacturer's declaration – electromagnetic emissions		
The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or the user of the HemoCue systems should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions	Group 1	The HemoCue systems uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment
RF emissions	Class B	The HemoCue systems are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Adapters

Country: EU/US/GB

Type: FW7333SM/12

Input: 100V~ – 240 V~, 50–60 Hz, 200 mA

Country: EU/US/GB

Type: HCA01

Input: 100 V~ – 240 V~/50–60 Hz/<500 mA

Warning

The device is tested according to standard IEC/EN 60601-1-2 and is found to comply with the standard.

Despite this compliance it is impossible to predict any possible effects by other nearby standing instruments, (stationary, portable or mobile units) or the possible impact of electromagnetic radiance. This is the reason we need to inform users of this Analyzer that disturbance from other equipment might affect the performance of the Analyzer. Should you note that this is the fact, please contact your local HemoCue distributor.

The HemoCue WBC is intended for use in the electromagnetic environment specified in Technical specifications. The customer or user of the HemoCue WBC system should assure that it is used in such an environment. The HemoCue WBC Analyzer uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.

The HemoCue WBC system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards (i.e. IEC 60950-1 for data processing equipment and IEC 60601-1 for medical equipment).

Furthermore all configurations must comply with the system standard, IEC 60601-1-1. Connections of additional equipment to the signal input or signal output connections is deemed to be “configuring a medical system”, and therefore assumes responsibility for the systems compliance within the requirements of IEC 60601-1-1. For further information, please contact your local distributor.

Warranty

The Analyzer carries a 24-month warranty from the day of receipt. After the expiration date of the warranty, service such as maintenance and repairs are offered at a fixed price.

Service and Disposal

The Analyzer should be cleaned as recommended under Maintenance prior to service or disposal. Consult local environmental authorities for proper disposal.

Spare Parts and Accessories

The following accessories and spare parts are available:

- Power adapter
- Cuvette holder
- HemoCue Cleaners
- HemoCue Lancets

Patents

The product is protected by the following patents (or Patent Pending):

SE 0 500 549, US 5 674 457

EP 0 821 784, SE 0 601 576

US 552 500, RCD 00563838

US 11/102 837, EP 06110903.9

US 11/717 675, PCT SE2007/000654

US 60/906 504

Symbols Used

Attention, see instructions for use

CE mark

Class II equipment

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

Serial port

Relative humidity, non-condensing

Temperature

Efficiency Level

References

1. HemoCue WBC Microcuvettes package insert.
2. Dacie and Lewis, Practical Haematology Tenth edition.
3. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline CLSI Document EP05-A2.

4. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline CLSI Document EP06-A.
5. Interference Testing in Clinical Chemistry; Approved Guideline CLSI Document EP07-A2.
6. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline Document CLSI EP09-A2.
7. Protocols for Determination of Limit of Detection and Limit of Quantification; Approved Guideline CLSI Document EP17-A
8. James O. Westgard, Basic Method Validation.

Manufacturer

HemoCue AB
Box 1204
SE-262 23 Ängelholm
Sweden
Phone:
Fax:
E-mail: info@HemoCue.se
www.HemoCue.com

+46 775 70 02 10

+46 775 70 02 12

A Quest Diagnostics Company

HemoCue AB, Box 1204, SE-262 23 Ängelholm, Sweden. Phone +46 775 70 02 10. Fax +46 775 70 02 12.
info@hemocue.se www.hemocue.com