

Non-Latex Biogel®

Skinsense® Indicator® Underglove

Blue synthetic Polychloroprene surgical underglove
with Biogel® coating
Can be used alone or as part of a Puncture Indication System



Sterile, powder-free, non-latex polychloroprene gloves that reduce the possibility of glove-related latex protein sensitization because they are made from a synthetic elastomer.

The underglove is blue in color and can be worn alone or as part of the Puncture Indication System, as the smooth surface is specifically engineered for unparalleled ease of double-donning.

Recommended for use in all general surgeries or surgeries where latex allergies are a concern for patients or clinicians.

Appropriate for use in the presence of cytotoxic agents², unless otherwise specified by the chemical manufacturer's guidelines.

Benefits

- Helps prevent glove punctures from going unnoticed
- Blue underglove worn under any Biogel® non-latex powder-free surgical glove creates the Puncture Indication System
- Blue color spot alerts wearer of punctures in outer latex glove, if a perforation occurs in the presence of fluid
- Unique Biogel coating on the inner surface
 - Makes gloves easy to don even with damp hands
 - Soothes the skin to help prevent moisture loss¹
- Beaded cuff for added security
- Every glove is air inflation tested and visually inspected for quality and user confidence
- Industry leading AQL freedom from holes of 0.65
- Powder-free to eliminate starch powder-related complications
- Low dermatitis potential - this product demonstrated reduced potential for sensitizing users to chemical additives

Features

- Curved finger design
- Smooth surface for easy double donning
- Blue colored

Material Information

- Synthetic polychloroprene copolymer
- Biogel polymer coating
- Non-leachable blue pigment



ACTUAL
COLOR

RE-ORDER # 406

#	SIZE	PAIRS	#	SIZE	PAIRS
40655	5 ½	50/Box	40675	7 ½	50/Box
40660	6	50/Box	40680	8	50/Box
40665	6 ½	50/Box	40685	8 ½	50/Box
40670	7	50/Box	40690	9	40/Box

4 BOXES PER CASE

Essential Information

Responsible company: Mölnlycke Health Care AB, Gamlestadsvägen 3C, SE-402 52 Gothenburg, Sweden
USA: Norcross GA USA 30092 1-800-843-8497
For Inquiries: biogel@molnlycke.com
E-mail address: biogel@molnlycke.com
Website address: www.molnlycke.com
Manufacturer: Made and packed in Malaysia by Regent Hospital Products. Regent Hospital Products is a member of the Mölnlycke Health Care Group.

¹ Source: Testing conducted at PRAC and overseen by a board-certified dermatologist, 2004. Report 04-007
² Breakthrough testing data available upon request. REPR 0618.

Biogel® Skinsense® Indicator® Underglove

PRODUCT SPECIFICATIONS

SIZE	MIN LENGTH (mm)	AVERAGE PALM WIDTH (mm)
5 ½	283 (+20/-10)	71 (±3)
6	285 (+20/-10)	77 (±3)
6 ½	285 (+20/-10)	85 (±3)
7	288 (+20/-10)	91 (±3)
7 ½	298 (+20/-10)	96 (±3)
8	299 (+20/-10)	103 (±3)
8 ½	301 (+20/-10)	109 (±3)
9	301 (+20/-10)	115 (±3)

TYPICAL THICKNESS PROFILE SINGLE WALL

Cuff	6.5 mils	0.165 mm
Palm	7.5 mils	0.190 mm
Finger	8.3 mils	0.210 mm

PHYSICAL GLOVE PROPERTIES

		STANDARD REQUIREMENT	BIOGEL
Force at Break (N):	Initial	≥9	Typically 12
	Aged	≥6	Typically 13
Typical Accelerator Analysis % w/w:			
Dithiocarbamate (DTC)		n/a	Below detection
Diphenyl thiourea (DPTU)		n/a	<0.25
Diphenyl guanidine (DPG)		n/a	<0.07
AQL Freedom from Holes (1000ml water leak test Post packing and irradiation) - Process Average Typically		1.5	0.65 <0.20%

General Information

Contraindications: Do not use this product if you have a known allergy to chemical additives.

Pyrogenicity: Each batch of Biogel gloves is certified to be either non-pyrogenic or to have a low endotoxin level (<0.5 EU/ml).

Product Standards: Biogel gloves are tested and manufactured to the following standards:
Quality / Environmental- ISO 9001, ISO 13485, ISO 14001
Product- ASTM D3577, EN455-1, EN455-2, EN455-3
Sterilisation- Gamma irradiation (25 kGy (2.5 Mrad)), EN552, EN556
Viral Penetration- Bacteriophage test, ASTM F1671
Allergenicity / Pyrogenicity- ISO 10993, (PART 5 and 10) / US Pharmacopoeia USP: (151 & 85)

Distributed by: Mölnlycke Health Care, LLC, 5550 Peachtree Parkway, Suite 500, Norcross, GA 30092, USA.

Storage: Store in a cool, dry place away from sources of heat or direct sunlight.

Packaging: One pair per pack, in a high quality inner wrap, packed into a film pack (constructed of a laminate of polyester and low-density polyethylene). 50 pairs per collation case for sizes 5.5 – 8.5; 40 pairs for size 9.0; 200 pairs per transit case for sizes 5.5 – 8.5; 160 pairs for size 9.0.

Disposal: Gloves & Outer Wrap dispose of as for clinical waste. Paper Inner Wrap, Collation Case & Transit Case can be recycled as paper or disposed of as for clinical waste.

Country of Origin: Malaysia

Date of issue: October 2010

For further information on any of the above, please contact your local dealer or Mölnlycke Health Care representative directly.

