

FOR USE WITH SOFIA ONLY



The Sofia Strep A FIA employs immunofluorescence technology to detect Group A Streptococcal antigens from throat swabs of symptomatic patients. All negative test results should be confirmed by bacterial culture because negative results do not preclude Group A Strep infection and should not be used as the sole basis for treatment. The test is intended for professional and laboratory use as an aid in the diagnosis of Group A Streptococcal infection.

SUMMARY AND EXPLANATION

Group A Streptococcus is one of the most common causes of acute upper respiratory tract infection. Early diagnosis and treatment of Group A Streptococcal pharyngitis has been shown to reduce the severity of symptoms and serious complications such as rheumatic fever and glomerulonephritis.^{1, 2} Conventional procedures for identification of Group A Streptococcus from throat swabs involve the culture, isolation, and subsequent identification of viable pathogen at 24 to 48 hours or longer for results.^{3, 4}

PRINCIPLE OF THE TEST

The Sofia Strep A FIA employs immunofluorescence technology that is used with the Sofia analyzer (Sofia) to detect Group A Streptococcal antigen.

The Sofia Strep A FIA involves the extraction of the antigenic components of the Group A Streptococcus (GAS) bacteria. The patient's Swab specimen is placed in the Reagent Tube containing the Reagent Solution, during which time the bacterial antigens are extracted, making them more accessible to the specific antibodies. An aliquot of the extracted specimen is dispensed into the Cassette sample well. From the sample well, the specimen migrates through a test strip containing various unique chemical environments. If Group A Streptococcal antigens are present, they will be bound by antibodies coupled to fluorescent microparticles that migrate through the test strip. The fluorescent microparticles containing bound antigen will be captured by antibodies at a defined location on the test strip where they are detected by Sofia. If antigens are not present, the fluorescent microparticles will not be trapped by the capture antibodies nor detected by Sofia.

Note: Depending upon the user's choice, the Cassette, now containing the specimen, is either placed directly inside Sofia for automatically timed development (WALK AWAY Mode) or placed on the counter or bench top for a manually timed development and then placed into Sofia (READ NOW Mode).

Sofia scans, measures, and interprets the immunofluorescent signal, using on-board method-specific algorithms. Sofia will then report the test results to the user (Positive, Negative, or Invalid) on its display screen, and it can print out the results via an integrated printer or transmit the results via an LIS connection.

Sofia Strep A FIA Page 1 of 15

REAGENTS AND MATERIALS SUPPLIED

25-Test Kit:

- Individually Packaged Cassettes (25): Polyclonal rabbit anti-Group A Streptococcus antibodies
- Reagent Tubes (25)
- Reagent Solution Bottles (25): 4M Sodium Nitrite and 0.2M Acetic Acid inside glass ampoule
- Fixed Volume Pipettes (25)
- Sterile Rayon Throat Swabs (25)
- Positive Control Swab (1): Swab is coated with heat-inactivated, non-infectious Group A Streptococcus
- Negative Control Swab (1): Swab is coated with heat-inactivated, non-infectious Group C Streptococcus
- Package Insert (1)
- Quick Reference Instructions (1)
- QC Card (located on kit box)
- Printer Paper (1)

MATERIALS NOT SUPPLIED IN KIT

- Timer or watch for use in READ NOW Mode
- Sofia analyzer
- Calibration Cassette (supplied with Sofia Installation Pack)

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use.
- Do not use the kit contents beyond the expiration date printed on the outside of the box.
- Use appropriate precautions in the collection, handling, storage, and disposal of patient samples and used kit contents.
- Use of Nitrile or Latex (or equivalent) gloves is recommended when handling patient samples.⁵
- Dispose of containers and used contents in accordance with Federal, State and Local requirements.
- Do not reuse any used Cassettes, Reagent Tubes, Fixed Volume Pipettes, solutions, or Control Swabs.
- The user should never open the Foil Pouch of the test Cassette exposing it to the ambient environment until the Cassette is ready for immediate use.
- Discard and do not use any damaged Cassette or material.
- The Reagent Solution contains an acidic solution. If the solution contacts the skin or eye, flush with copious amounts of water.
- Testing should be performed in an area with adequate ventilation.
- For more information, consult the Safety Data Sheet available on quidel.com.
- The Reagent Solution Bottle contains glass, break cautiously.
- If the Reagent Solution Bottle is missing the glass ampoule or the solution is green prior to the breaking of the ampoule, discard and use another Reagent Solution Bottle.
- To obtain accurate results, the Package Insert instructions must be followed.
- The Calibration Cassette must be kept in the provided storage pouch between uses.
- Inadequate or inappropriate specimen collection, storage, and transport may yield false test results.
- Specimen collection and handling procedures require specific training and guidance.
- If transport media will be used, use only the Transport Media and configuration recommended in this Package Insert.
- Use rayon-tipped swabs to collect throat specimens. The performance claims in the Performance Characteristics section were obtained with the Swabs provided in the kit. Use of the provided Swabs is recommended. Do not use calcium alginate, cotton-tipped or wooden shaft swabs.

Sofia Strep A FIA Page 2 of 15

- Do not write on the barcode of the Cassette. This is used by Sofia to identify the type of test being run and to identify the individual Cassette so as to prevent a second read of the Cassette by the same Sofia.
- Once a Cassette has been successfully scanned by Sofia, do not attempt to scan the Cassette again in the same Sofia. The barcode on the Cassette contains a unique identifier that will prevent Sofia from performing a second read on a previously scanned Cassette.
- As the detection reagent is a fluorescent compound, no visible results will form on the test strip. Sofia must be used for result interpretation.

KIT STORAGE AND STABILITY

Store the kit at room temperature, 59°F to 86°F (15°C to 30°C), out of direct sunlight. Kit contents are stable until the expiration date printed on the outer box. Do not freeze.

QUALITY CONTROL

There are three types of Quality Control for Sofia and Strep A FIA: Sofia Calibration Check Procedure, built-in procedural control features, and External Controls.

Sofia Calibration Check Procedure


Note: This is a "Calibration Check" procedure.

The Calibration Check Procedure should be performed every 30 days. Sofia is set to remind the user to complete the Calibration Check Procedure.

The Calibration Check is a required function that checks the Sofia optics and calculation systems using a specific Calibration Cassette. This Calibration Cassette is shipped with the Sofia Installation Pack. Refer to the Sofia User Manual for details regarding the Calibration Check Procedure.

Important: Ensure that the Calibration Cassette is stored in the provided storage pouch between uses to protect it from exposure to light.


1. To check the calibration of Sofia, select "Calibration" from the Main Menu.



2. Following the prompts, insert the Calibration Cassette into Sofia and gently close the drawer. Sofia performs the Calibration Check automatically with no user input required.

Sofia indicates when the Calibration Check is completed. Select **OK** to return to the Main Menu.

NOTE: If the Calibration Check does not pass, notify the on-site Supervisor or contact Quidel Technical Support for assistance from 7:00 a.m. to 5:00 p.m. Pacific Time at 800.874.1517 (in the U.S.); 858.552.1100 (outside the U.S.); Fax: 858.455.4960; custserv@quidel.com (Customer Service); technicalsupport@quidel.com (Technical Support) or contact your local distributor.



Sofia Strep A FIA Page 3 of 15

Built-in Procedural Controls

The Sofia Strep A FIA contains two built-in procedural control features. The manufacturer's recommendation for daily control is to document these built-in procedural controls for the first sample tested each day.

A control of the extraction procedure is provided by a color change from clear to green as the Reagent Solution is mixed. The color change is an indication of Reagent Solution integrity and is also an indication that the extraction procedure was performed correctly.

Each time a test is run in Sofia, a procedural control is interpreted by Sofia and the result is displayed on the Sofia screen. This information is automatically logged in Sofia with each test result. A valid result obtained with the procedural control demonstrates that the extracted specimen flowed correctly and the functional integrity of the Cassette was maintained. This procedural control is interpreted by Sofia after the Cassette has developed for 5 minutes. If the specimen has not flowed correctly, Sofia will indicate that the result is invalid. Should this occur, review the procedure and repeat the test with a new patient sample and a new test Cassette.



For example: This result shows that an invalid result had occurred.

External Quality Control

External Controls may also be used to demonstrate that the reagents and assay procedure perform properly. Quidel recommends that Positive and Negative Controls be run:

- once for each untrained operator;
- once for each new shipment of kits-provided that each different lot received in the shipment is tested; and
- as deemed additionally necessary by your internal quality control procedures, and in accordance with Local, State and Federal regulations or accreditation requirements.

The user must first select Run QC on the Main Menu of Sofia and then, when prompted, scan the QC Card (located on the kit box). This card provides information specific to the kit lot, including lot number and expiration date.

Sofia will prompt the user to select the desired mode (WALK AWAY or READ NOW) and then to run the External Control Swabs.

External Positive and Negative Control Swabs are supplied in the kit and should be tested using the Test Procedure provided in this Package Insert or in the Quick Reference Instructions. Additional External Control Swabs may be obtained separately by contacting Quidel's Customer Support Services at 800.874.1517 (in the U.S.) or 858.552.1100 (outside the U.S.).

When the QC test is complete, each result will be displayed as "Passed" or "Failed" for the Positive Control and the Negative Control.

Do not perform patient tests or report patient test results if the QC test does not produce the expected results. Repeat the test or contact Quidel Technical Support before testing patient specimens if a "Failed" result is obtained with the External Controls.

Sofia Strep A FIA Page 4 of 15

SPECIMEN COLLECTION AND HANDLING

SPECIMEN COLLECTION

Use the Rayon-tipped Swabs provided in the kit to collect throat specimens. The performance claims listed in the Performance Characteristics section were obtained with the Swabs provided in the kit. Do not use calcium alginate, cotton-tipped or wooden shaft swabs. Collect throat specimens by standard clinical methods. Depress the tongue with a tongue blade or spoon. Rub the Swab on the back of the throat, both of the tonsils, the uvula and the posterior pharynx. Be careful not to touch the tongue, sides or top of the mouth with the Swab. Consult standard reference procedures such as the collection method described by Facklam.⁶

SPECIMEN TRANSPORT AND STORAGE

It is recommended that Swab specimens be processed as soon as possible after collection. Swabs can be held in any clean, dry plastic tube or sleeve up to 24 hours at room temperature (15°C to 30°C), or refrigerated (2°C to 8°C) up to 48 hours. The following transport media and storage conditions have been tested and are also acceptable (Table 1):

Table 1
Recommended Transport Media

	Recommended Storage Condition	
Transport Media	2°C to 8°C	Ambient Temperature
BD BBL CultureSwab with Liquid Stuarts Media (#220109)*	48 hours	24 hours
Remel BactiSwab with Liquid Amies Media (#R723095)*	48 hours	24 hours

^{*}These transport media systems preserve the sample on the Swab tip via contact with a media-moistened sponge.

If a culture is desired, lightly streak the Swab on a 5% sheep blood agar plate before using the Swab in the Sofia Strep A FIA. Do not perform the Strep A FIA before streaking the Swab, as the Reagent Solution will destroy the bacteria on the Swab, thereby rendering the organism incapable of successful culturing. Alternatively, two throat Swab specimens can be obtained; in this case one can be used separately for culture and the other for the Sofia Strep A FIA.


TEST PROCEDURE

Important:

- All specimens must be at room temperature (15°C to 30°C) before beginning the test.
- Gloves should be worn when handling human samples.
- Do not use the Reagent Solution if it is green prior to breaking the ampoule.
- Do not open the foil pouch of the test Cassette until it is ready for immediate use.

Expiration date: Check expiration on outer box before using. Do not use any test Cassette past the expiration date on the label.

- 1. Verify that Sofia is set to the desired Mode: **WALK AWAY** or **READ NOW.** See the "Using Sofia" section for more information.
- **2.** Squeeze **ONCE** to break the glass ampoule inside the Reagent Solution Bottle prior to running the assay.



Sofia Strep A FIA Page 5 of 15

3. Vigorously shake the Bottle 5 times to mix the Solutions. Solution should turn green after the ampoule is broken.



4. Remove the cap. Holding Bottle vertically, fill the Reagent Tube to the line [approximately 6 drops].



5. Immediately add the patient Swab sample to the Reagent Tube. Vigorously mix the solutions by plunging the Swab **5 times** in an up and down motion in the Tube.

NOTE: Best results are obtained when the specimen is vigorously extracted in the solution.




Vigorously plunge up and down 5x

6. Leave the Swab in the Reagent Tube for 1 minute.




7. **Vigorously** mix the solution again by plunging the Swab 5 times in an up and down motion in the Tube.



Sofia Strep A FIA Page 6 of 15

8. Express as much liquid as possible from the Swab by **squeezing the sides of the Tube as the Swab is withdrawn.**


Discard the Swab in accordance with your biohazard waste disposal protocol.



9. Fill the provided Yellow 100 µL Fixed Volume Pipette with the sample:


To fill the Fixed Volume Pipette with the sample:

- a) FIRMLY squeeze the top bulb.
- b) Still squeezing, place the Pipette tip into the sample.
- c) With the Pipette tip still in the sample, release pressure on the bulb to fill the Pipette.



10. Firmly squeeze the top bulb to empty the contents of the Fixed Volume Pipette into the Cassette sample well. Extra liquid in the overflow bulb is OK.

NOTE: The Fixed Volume Pipette is designed to collect and dispense the correct amount of patient sample. Discard the Pipette in your biohazard waste.



11. Proceed to the next section, "Using Sofia," to complete the test.

USING SOFIA

WALK AWAY/READ NOW Modes

Refer to the Sofia User Manual for operating instructions.

Sofia may be set to two different modes (WALK AWAY and READ NOW). The procedures for each mode are described below.

WALK AWAY Mode

In WALK AWAY Mode, the user immediately inserts the Cassette into Sofia. The user then returns after 5 minutes to get the test result. In this mode, Sofia will automatically time the test development before scanning and displaying the test result.

READ NOW Mode

Allow the test to develop for the full 5 minutes BEFORE placing it into Sofia.

The user must first place the Cassette onto the counter or bench top for 5 minutes (outside of Sofia) and manually time this development step. Then, the user inserts the Cassette into Sofia. In READ NOW Mode, Sofia will scan and display the test result in less than 1 minute. Note: Results will remain stable for an additional 10 minutes after the recommended development time of 5 minutes.

Sofia Strep A FIA Page 7 of 15

Tips for Batch Testing


In order to make batch testing easier, the user can prepare one or more Reagent Solution Bottles in advance of testing samples. The user can break the ampoule inside each Reagent Solution Bottle, shake to mix the solutions, and then store the capped Bottles on the bench top at room temperature for up to 12 hours without loss of activity before using with Swab sample(s).

Critically important, the user should never open the foil pouch thus exposing the test Cassette to the ambient environment, until it is ready for immediate use.

Run Test


1. Input the User ID using the handheld barcode scanner or manually enter the data using the key pad.

NOTE: If you mistakenly scan the wrong barcode, use the Arrow Buttons on the Sofia key pad to re-highlight the field. Then simply rescan using the correct barcode, and the previous one will be overwritten with the correct barcode.





2. Input Patient ID and/or Order # using the handheld barcode scanner or manually enter the data using the key pad.





3. Press Start Test and the Sofia drawer will automatically open.



4. Verify that the correct development mode, WALK AWAY or READ NOW, has been selected. Insert the prepared patient test Cassette into the drawer of Sofia and close the drawer.



Sofia Strep A FIA Page 8 of 15

5. Sofia will start automatically and display the progress as shown in example below. In the WALK AWAY Mode, the test results will be displayed on the screen in approximately 5 minutes. In the READ NOW Mode, the test results will be displayed on the screen in less than 1 minute. See "Interpretation of Results" section.



For example: This display shows that the test in WALK AWAY mode has 4 minutes, 13 seconds remaining. Sofia will read and display the results in about 5 minutes.

INTERPRETATION OF RESULTS

When the test is complete, the results will be displayed on the Sofia screen. The results can be automatically printed on the integrated printer, if this option is selected. Test Lines, which are fluorescent, will never be visible to the naked eye.

The Sofia screen will display results for the procedural control as being "valid or invalid," and will provide a positive or negative result for Strep A. If the procedural control is "invalid," retest with a new patient sample and a new test Cassette.

Positive Result:



For example: This display shows a valid positive result for Strep A.

NOTE: A positive result does not rule out co-infections with other pathogens.

Negative Result:



For example: This display shows a valid negative result for Strep A.

NOTE: A negative result does not rule out possible other infections.

Invalid Result:



For example: This result shows that an invalid result was obtained.

Invalid Result: If the test is invalid, a new test should be performed with a new patient sample and a new test Cassette.

Sofia Strep A FIA Page 9 of 15

LIMITATIONS

- The contents of this kit are to be used for the qualitative detection of Group A Streptococcal antigens from throat Swab specimens.
- The test detects both viable and nonviable Group A Streptococcus bacteria and may yield a positive result in the absence of living organisms.
- Respiratory infections, including pharyngitis, can be caused by Streptococcus from serogroups other than Group A, as well as other pathogens.
- The Sofia Strep A FIA will not differentiate asymptomatic carriers of Group A Streptococcus from those exhibiting Streptococcal infection.⁷
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test or if the sample was collected, transported, or stored improperly.
- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- Patients with symptoms and an antigen negative test should have a follow-up culture.¹
- Test results must be evaluated in conjunction with other clinical data available to the physician.
- Negative test results do not rule out possible other infections.
- Positive test results do not rule out co-infections with other pathogens.

EXPECTED VALUES

Group A Streptococcus bacteria are responsible for about 19% of all upper respiratory tract infections.⁸ Infection is most prevalent in winter and early spring, with most cases arising in patients living in highly populated areas. Consistent with these figures, in the multi-center clinical study conducted by Quidel during 2011 and 2012, 17.4% (128/736) of the patients presenting with pharyngitis were found to be culture positive for Strep A. Nearly half of these subjects, 46%, were male. The subjects' ages ranged from 3-72 years and 88% (647/736) were children (3-17 years of age).

PERFORMANCE CHARACTERISTICS

Sofia Strep A FIA Performance vs. Cell Culture

The performance of the Sofia Strep A FIA was compared to standard bacterial culture and identification in a multi-center clinical field study. This study was conducted by health care personnel during 2011 and 2012 at eight (8) distinct sites in various geographical regions within the United States and two (2) sites in Australia. In this multi-center, point-of-care (POC) field trial, two (2) throat Swabs were collected from 736 patients with symptoms suggestive of bacterial pharyngitis.

One (1) throat Swab was transported on cold ice packs to a central Reference Laboratory, streaked on a sheep blood agar plate (SBA) and cultured for up to 48 hours. Immediately after streaking, this same Swab was tested in the rapid Sofia Strep A FIA. The performance of the Sofia Strep A FIA was determined by comparison of the rapid test result to the corresponding culture result. The results from these analyses are presented in Tables 2, 3a, and 3b.

Table 2 Sofia Strep A FIA Results: Combined

Culture			
	Pos	Neg	Total:
Sofia Pos	116	24	140
Sofia Neg	12	584	596
Total:	128	608	736

Sens. = 90.6% (116/128) (95% CI: 84.3%-94.6%)

Spec. = 96.1% (584/608) (95% CI: 94.2%-97.3%)

PPV = 82.9% (116/140)

NPV = 98.0% (584/596)

Prev. = 17.4% (128/736)

Sofia Strep A FIA Page 10 of 15

Table 3a Sofia Strep A FIA Results: READ NOW Mode

Culture			
	Pos	Neg	Total:
Sofia Pos	100	23	123
Sofia Neg	12	549	561

572

684

Sens. = 89.3% (100/112) (95% CI: 82.2%-93.8%)

Spec. = 96.0% (549/572)

(95% CI: 94.0%-97.3%)

PPV = 81.3% (100/123)

NPV = 97.9% (549/561)

Prev. = 16.4% (112/684)

Table 3b Sofia Strep A FIA Results: WALK AWAY Mode

	Cult	ture	
	Pos	Neg	Total:
Sofia Pos	16	1	17
Sofia Neg	0	35	35
Total:	16	36	52

112

Total:

Sens. = 100% (16/16) (95% CI: 80.6%-100%)

Spec. = 97.2% (35/36)

(95% CI: 85.8%-99.5%)

PPV = 94.1% (16/17)

NPV = 100% (35/35)

Prev. = 30.8% (16/52)

Reproducibility Studies

The reproducibility of the Sofia Strep A FIA was evaluated at three (3) different laboratories. Two (2) different operators at each site tested a series of coded, contrived samples, prepared in negative clinical matrix, ranging from low negative to moderate positive Group A Streptococcus. The inter-laboratory agreement (Table 4) for negative samples was 96.7%-100% and 96.7%-100% for positive samples. The intra-laboratory agreement (Table 5) for all samples ranged from 97.5%-99.2%.

Table 4
Sofia Strep A FIA Reproducibility Study Inter-laboratory Agreement

Site	Low Negative (no bacteria) (0 cfu/test)	High Negative (C₅) (1.5x10³ cfu/test)	Low Positive (C ₉₅) (3.0x10 ³ cfu/test)	Mod. Positive (C _{3x}) (2.8x10 ⁴ cfu/test)
1	30/30	30/30	28/30	30/30
2	30/30	29/30	30/30	30/30
3	30/30	28/30	29/30	30/30
Total	90/90	87/90	87/90	90/90
% Overall Agreement (95% CI)	100% (90/90) (95.9%-100%)	96.7% (87/90) (90.7%-98.9%)	96.7% (87/90) (90.7%-98.9%)	100% (90/90) (95.9%-100%)

Table 5
Sofia Strep A FIA Reproducibility Study Intra-laboratory Agreement

Site	Low Negative (no bacteria) (0 cfu/test)	High Negative (C₅) (1.5x10³ cfu/test)	Low Positive (C ₉₅) (3.0x10³ cfu/test)	Mod. Positive (C3x) (2.8x104 cfu/test)	% Overall Agreement (95% CI)
1	30/30	30/30	28/30	30/30	98.3% (118/120) (94.1%-99.5%)
2	30/30	29/30	30/30	30/30	99.2% (119/120) (95.4%-99.9%)
3	30/30	28/30	29/30	30/30	97.5% (117/120) (92.9%-99.1%)

Sofia Strep A FIA Page 11 of 15

Limit of Detection

The limit of detection (LOD) for the Sofia Strep A FIA was determined using three (3) strains of Group A Streptococcus pyogenes. The LOD ranged from $9x10^3-2x10^4$ colony forming units (cfu)/test (Table 6).

Table 6
Sofia Strep A FIA Limits of Detection for Three Streptococcus pyogenes Strains

Strain	Minimum Detectable Level*
Bruno [CIP 104226]	1.86x10 ⁴ cfu/test
CDC-SS-1402	9.24x10³ cfu/test
CDC-SS-1460	2.34x10 ⁴ cfu/test

cfu/test = colony forming units/test

Analytical Reactivity

Analytical reactivity for the Sofia Strep A FIA was determined using 15 strains of *Streptococcus pyogenes*. Each strain listed below in Table 7 produced positive results in the assay.

Table 7
Analytical Reactivity

Streptococcus pyogenes Strain	Test Quantity*
Strain #1 (ATCC 19615)	6.5x10 ⁴ cfu/test
Strain #2 (ATCC 700942)	7.4x10 ⁴ cfu/test
Strain #3 (ATCC 700952)	8.3x10 ⁴ cfu/test
Strain #4 (Field Clinical Isolate)	3.1x10 ⁴ cfu/test
Strain #5 (Field Clinical Isolate)	7.6x10 ⁴ cfu/test
Strain #6 (Field Clinical Isolate)	7.1x10⁵ cfu/test
Strain #7 (Field Clinical Isolate)	6.3x10 ⁴ cfu/test
Strain #8 (Field Clinical Isolate)	6.3x10 ⁴ cfu/test
Strain #9 (Field Clinical Isolate)	5.3x10 ⁴ cfu/test
Strain #10 (ATCC 700482)	6.5x10 ⁴ cfu/test
Strain #11 (ATCC BAA 1315)	7.2x10 ⁴ cfu/test
Strain #12 (ATCC 700459)	5.4x10 ⁴ cfu/test
Strain #13 (ATCC 12203)	6.9x10 ⁴ cfu/test
Strain #14 (ATCC 700944)	5.3x10 ⁴ cfu/test
Strain #15 (Field Clinical Isolate)	7.0x10⁴ cfu/test

cfu/test = colony forming units/test

Sofia Strep A FIA Page 12 of 15

^{*}The levels of bacteria were determined by limiting dilution, bacterial culture, and colony counting to give cfu/test.

^{*}The levels of bacteria were determined by limiting dilution, bacterial culture, and colony counting to give cfu/test.

Analytical Specificity

Cross Reactivity

The cross reactivity of the Sofia Strep A FIA was evaluated with a total of 61 non-Group A Streptococcus bacterial and fungal microorganisms, and 26 viral isolates. None of the organisms or viruses listed below in Table 8 showed any sign of cross reactivity in the assay. When the same organisms in Table 8 were pre-mixed with Group A Strep and tested in the Sofia Strep A FIA, all results were positive also indicating that the potential cross-reactants did not interfere with the detection of Strep A.

Table 8
Analytical Specificity and Cross Reactivity

Organism/Virus	Test Quantity*	Organism/Virus	Test Quantity*
Arcanobacterium haemolyticum	3x10⁵ cfu/test	Streptococcus sp. Group C strain #4	3x10 ⁶ cfu/test
Bacteroides fragilis	3x10 ⁷ cfu/test	Streptococcus sp. Group C strain #5	3x10 ⁵ cfu/test
Bordetella pertussis	3x10 ⁷ cfu/test	Streptococcus sp. Group D strain #1	3x10 ⁶ cfu/test
Candida albicans	3x10⁴ cfu/test	Streptococcus sp. Group D strain #2	3x10 ⁶ cfu/test
Corynebacterium diphtheria	3x10⁵ cfu/test	Streptococcus sp. Group D strain #3	3x10 ⁶ cfu/test
Corynebacterium pseudodiphtheriticum	3x10 ⁶ cfu/test	Streptococcus sp. Group D strain #4	3x10 ⁶ cfu/test
Enterococcus faecalis	3x10 ⁶ cfu/test	Streptococcus sp. Group D strain #5	3x10 ⁶ cfu/test
Enterococcus faecium	3x10 ⁶ cfu/test	Streptococcus sp. Group F strain #1	3x10⁵ cfu/test
Escherichia coli	1.5x10 ⁷ cfu/test	Streptococcus sp. Group F strain #2	3x10 ⁶ cfu/test
Fusobacterium necrophorum	3x10 ⁶ cfu/test	Streptococcus sp. Group F strain #3	3x10 ⁶ cfu/test
Haemophilus influenzae	3x10 ⁷ cfu/test	Streptococcus sp. Group F strain #4	3x10 ⁵ cfu/test
Haemophilus parahaemolyticus	3x10 ⁶ cfu/test	Streptococcus sp. Group F strain #5	3x10 ⁵ cfu/test
Klebsiella pneumoniae	3x10 ⁷ cfu/test	Streptococcus sp. Group G strain #1	3x10 ⁷ cfu/test
Moraxella catarrhalis	3x10 ⁶ cfu/test	Streptococcus sp. Group G strain #2	3x10 ⁶ cfu/test
Neisseria lactamica	3x10 ⁶ cfu/test	Streptococcus sp. Group G strain #3	3x10 ⁶ cfu/test
Neisseria gonorrhoeae	3x10 ⁶ cfu/test	Streptococcus sp. Group G strain #4	3x10 ⁶ cfu/test
Neisseria meningitidis	3x10 ⁶ cfu/test	Yersinia enterocolitica	3x10 ⁶ cfu/test
Neisseria sicca	3x10 ⁷ cfu/test	Adenovirus Type 1	3x10 ¹¹ TCID ₅₀ /test
Neisseria subflava	3x10 ⁷ cfu/test	Adenovirus Type 3	3x10 ⁵ TCID ₅₀ /test
Proteus vulgaris	3x10 ⁷ cfu/test	Adenovirus Type 4	1.5x10 ² TCID ₅₀ /test
Pseudomonas aeruginosa	3x10 ⁶ cfu/test	Adenovirus Type 5	3x10 ⁵ TCID ₅₀ /test
Serratia marcescens	3x10 ⁷ cfu/test	Adenovirus Type 11	3x10 ⁴ TCID ₅₀ /test
Staphylococcus aureus	3x10 ⁶ cfu/test	Coronavirus 229E	3x10 ⁴ TCID ₅₀ /test
Staphylococcus epidermidis	3x10 ⁶ cfu/test	Coronavirus OC43	3x10 ⁴ TCID ₅₀ /test
Staphylococcus haemolyticus	3x10⁵ cfu/test	Coxsackievirus B5 (Faulkner)	3x106 TCID50/test
Staphylococcus intermedius	3x10⁵ cfu/test	Cytomegalovirus	3x10 ³ TCID ₅₀ /test
Staphylococcus saprophyticus	3x10 ⁶ cfu/test	Echovirus Type 3	1.5x10 ⁴ TCID ₅₀ /test
Streptococcus anginosus	3x10 ⁶ cfu/test	Epstein Barr virus	3x10 ⁷ TCID ₅₀ /test
Streptococcus gordonii	3x10 ⁴ cfu/test	Herpes Simplex virus 1	3x10 ⁴ TCID ₅₀ /test
Streptococcus mitis	3x10 ⁴ cfu/test	Herpes Simplex virus 2	3x10 ⁴ TCID ₅₀ /test
Streptococcus mutans	3x10 ⁶ cfu/test	Influenza A H1N1	3x10 ⁴ TCID ₅₀ /test
Streptococcus oralis	3x10 ⁶ cfu/test	Influenza A H3N2	3x10 ⁴ TCID ₅₀ /test
Streptococcus parasanguis	3x10 ⁶ cfu/test	Influenza B Hong Kong	3x10 ⁴ TCID ₅₀ /test
Streptococcus pneumoniae	3x10 ⁶ cfu/test	Influenza B Panama	1.5x10 ⁴ TCID ₅₀ /test
Streptococcus salivarius	3x10 ⁵ cfu/test	Influenza C Taylor	1.5x10 ⁴ TCID ₅₀ /test
Streptococcus sanguinis	3x10 ⁶ cfu/test	Measles (Edmonston)	3x10 ⁴ TCID ₅₀ /test
Streptococcus sp. Group B strain #1	3x10 ⁶ cfu/test	Mumps (Enders)	3x10³ TCID50/test
Streptococcus sp. Group B strain #2	3x10 ⁶ cfu/test	Parainfluenza virus 1	3x10 ⁴ TCID ₅₀ /test
Streptococcus sp. Group B strain #3	3x10 ⁶ cfu/test	Parainfluenza virus 2	1.2 TCID50/test
Streptococcus sp. Group B strain #4	3x10 ⁶ cfu/test	Parainfluenza virus 3	3x10 ⁶ TCID ₅₀ /test
Streptococcus sp. Group B strain #5	3x10 ⁶ cfu/test	Parainfluenza virus 4A	3x10 ⁴ TCID ₅₀ /test
Streptococcus sp. Group C strain #1	3x10 ⁶ cfu/test	Rhinovirus Type 15	3x10 ⁴ TCID ₅₀ /test
Streptococcus sp. Group C strain #2	3x10 ⁶ cfu/test	Rhinovirus Type 1B	3x103 TCID50/test
Streptococcus sp. Group C strain #3	3x10 ⁶ cfu/test		

cfu/test = colony forming units/test TCID50/test = 50% tissue culture infectious dose

Sofia Strep A FIA Page 13 of 15

^{*}The levels of bacteria were determined by limiting dilution, bacterial culture, and colony counting to give cfu/test. Virus concentrations were determined by standard virology methods, Reed-Muench.

Interfering Substances

Several over-the-counter (OTC) products, whole blood, and blood agar were evaluated and did not interfere with the Sofia Strep A FIA at the levels tested (Table 9).

Table 9
Non-interfering Substances

Substance	Concentration
Crest Pro-Health Night Mint (Cetylpyridnium chloride)	25% v/v
Listerine Antiseptic (Eucalyptol, Menthol, Methyl salicylate, and Thymol)	25% v/v
Listerine Cool Mint (Eucalyptol, Menthol, Methyl salicylate, and Thymol)	25% v/v
Cepacol Dual Relief Spray (Benzocaine and Menthol)	25% v/v
Chloraseptic Max: Sore Throat Relief (Phenol and Glycerin)	25% v/v
Children's Dimetapp DM Cold & Cough Elixir (Brompheniramine maleate, Dextromethorphan HBr, and Phenylephrine HCl)	25% v/v
Children's Wal-Tap Elixir Cold & Allergy (Brompheniramine maleate and Phenylephrine HCl)	25% v/v
Children's Wal-Tap DM Elixir Cold & Cough (Brompheniramine maleate, Dextromethorphan HBr, and Phenylephrine HCl)	25% v/v
Rite Aid Tussin CF (Dextromethorphan HBr, Guaifenesin, and Phenylephrine HCl)	25% v/v
Robitussin Cough & Cold-CF Max (Dextromethorphan HBr, Guaifenesin, and Phenylephrine HCl)	25% v/v
Robitussin Nighttime Cough, Cold, & Flu (Acetaminophen, Diphenhydramine HCl, and Phenylephrine HCl)	25% v/v
Cepacol Sore Throat: Cherry Flavor (Benzocaine and Menthol)	25% w/v
Halls Cherry Mentholyptus (Menthol)	25% w/v
Halls Mentholyptus (Menthol)	25% w/v
Ricola Mountain Herb Throat Drops-Sugar Free (Menthol)	25% w/v
Sucrets Complete-Vapor Cherry (Dyclonine Hydrochloride and Menthol)	25% w/v
Sucrets Complete-Cool Citrus (Dyclonine Hydrochloride and Menthol)	25% w/v
Chlorasceptic Throat Drops-Cherry (Phenol and Glycerin)	25% w/v
BreathSavers 3 Hour Mint-Spearmint (Cetylpyridnium chloride)	25% w/v
Tic Tac Freshmints (Eucalyptol, Menthol, Methylsalicylate, and Thymol)	25% w/v
Whole Blood	5% v/v
Sheep Blood Agar (5% Sheep Blood)	2.16 mg/mL
Horse Blood Agar (5% Horse Blood)	1.67 mg/mL

ASSISTANCE

If you have any questions regarding the use of this product or if you want to report a test system problem, please call Quidel's Technical Support Number 800.874.1517 (in the U.S.) or 858.552.1100, Monday through Friday, from 7:00 a.m. to 5:00 p.m., Pacific Time. If outside the U.S., contact your local distributor or technical support@quidel.com.

REFERENCES

- 1. American Academy of Pediatrics. [Group Streptococcal Infection]. In: Pickering L.K., Baker C.J., Kimberlin D.W., Long S.S., eds. Red Book: 2009 Report of the Committee on Infectious Diseases. 28th ed. Elk Grove Village, IL: America Academy of Pediatrics, 2009, pp. 616-628.
- 2. Youmans G.P., Paterson P.Y., and Sommer H.M., Upper Respiratory Tract Infection: General Considerations, in the Biological and Clinical Basis of Infectious Disease, W.B. Saunders Co., Philadelphia, 1980, pp. 177-183.
- **3.** Versalovic J., Carroll K.C., Jorgensen J.H., Funke G., Landry M.L., Warnock D.W., eds. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, DC. Chapter 20, Streptococcus, 2011, pp. 331-349.
- 4. Youmans G.P., Paterson P.Y., and Sommers H.M. The Biological and Clinical Basis of Infectious Disease, 1975, pp. 172-185.
- **5.** Biosafety in Microbiological and Biomedical Laboratories, 5th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (2007).
- **6.** Facklam R.R. and Washington J.A. Manual of Clinical Microbiology 5th Edition, 1991, p. 241.
- 7. Rammelkamp C.H., Jr. Principles of Internal Medicine, 8th Edition, 1977, pp. 814-820.
- 8. Lauer B.A., Reller L.D., and Mirrett S., Journal of Clinical Microbiology, 17:338-340, 1983.

Sofia Strep A FIA Page 14 of 15



20231 - Sofia Strep A FIA - 25 Tests (multi-language) 20253 - Sofia Strep A FIA - 25 Tests





EC REP

MDSS GmbH Schiffgraben 41 30175 Hannover, Germany



Quidel Corporation 10165 McKellar Court San Diego, CA 92121 USA **quidel.com**

1243801EN01 (03/15)

REF CE mark of conformity Catalogue number EC REP LOT Authorized Representative in Batch code the European Community Manufacturer Use by Intended use Temperature limitation IVD Consult instructions for use For In Vitro diagnostic use CONT Contains sufficient for 25 determinations Contents/Contains CONTROL + CONTROL Positive control Negative control

Sofia Strep A FIA Page 15 of 15