

C. difficile molecular testing with connectivity for easy results monitoring

The Solana C. difficile Assay is an *in vitro* diagnostic test for the direct, qualitative detection of toxigenic *Clostridium difficile* in unformed stool specimens of patients suspected of having *Clostridium difficile*-infection (CDI). Change the way you monitor *C. difficile* results with Virena®. Virena provides an electronic patient de-identified data management system to monitor results on a daily basis. This allows healthcare providers and healthcare facilities the ability to quickly track *C. difficile* infections in their patient population and facilities.

Prompt initiation of patient management and therapy


- Results in just 35 minutes
- High volume throughput capabilities;
 Batch up to 12 tests per runs
- Monitor results in near real-time

Ease and Affordability

- Simple procedure
- Sustainable reagent pricing
- Flexible instrument placement options



Procedure



Clinical Performance

Performance Characteristics of the Solana C. difficile Assay Compared to Broth Enhanced Toxigenic Bacterial Culture (Combined Sites – Raw Specimen)

| Enhanced Toxigenic Culture | | | | | | | | |
|----------------------------|-------|-----|-----|-------|-------------|-------|--|--|
| | | POS | NEG | Total | Sensitivity | 93.0% | | |
| Solana C. difficile | POS | 107 | 6* | 113 | Specificity | 99.2% | | |
| | NEG | 8** | 731 | 739 | | | | |
| | Total | 115 | 737 | 852 | | | | |

^{*}Three (3) of the six (6) specimens were positive for C. difficile toxin gene DNA by an alternative molecular device, three (3) were negative.

Performance Characteristics of the Solana C. difficile Assay Compared to FDA-cleared molecular device (Combined Sites – Raw Specimen)

| FDA-cleared Molecular Device | | | | | | | | |
|------------------------------|-------|-----|-----|-------|-------------|-------|--|--|
| | | POS | NEG | Total | Sensitivity | 97.0% | | |
| Solana C. difficile | POS | 97 | 16* | 113 | Specificity | 97.9% | | |
| | NEG | 3** | 736 | 739 | | | | |
| | Total | 100 | 752 | 852 | | | | |

^{*}Twelve (12) of the sixteen (16) specimens were positive for C. difficile toxin gene DNA by an alternative molecular device, four (4) were negative.





^{**}Six (6) of eight (8) specimens were found positive for C. difficile toxin gene DNA by an alternative molecular device, two (2) were negative. Refer to Package Insert for additional performance claims.

^{**}Two (2) of three (3) specimens were found positive for *C. difficile* toxin gene DNA by an alternative molecular device, one (1) was negative. Refer to Package Insert for additional performance claims.