ri-thermo[®] fastPRObe

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Instruction for use

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Table of Content


1. Introduction
 - 1.1 Important Information read prior to Start-up
 - 1.2 Safety Information and Electromagnetic Compatibility
 - 1.3 Packaging Symbols
 - 1.4 Intended Use
 - 1.4.1 Indication
 - 1.4.2 Contraindication
 - 1.4.3 Intended patient population
 - 1.4.4 Intended operators / users
 - 1.4.5 Required skills / training of operators
 - 1.4.6 Environmental conditions
 - 1.5 Warnings / Attention
 - 1.6 User Responsibility
2. Using the Device for the first time
 - 2.1 Scope of Supply
 - 2.2 Device Function
3. Information about Body Temperature
4. Operation and Function
 - 4.1 Commissioning
 - 4.1.1 Insert the Batteries
 - 4.2 Icon Identification
 - 4.3 Probe covers - Applying & Removing
 - 4.4 Change Probes
 - 4.5 Fast Measuring Mode
 - 4.6 Direct Measuring Mode
 - 4.7 Oral & Axillary Temperature Taking
 - 4.8 Rectal Temperature Taking
 - 4.9 Setup Mode
 - 4.10 Recall last Temperature
 - 4.11 Pulse Timer Mode
5. Cleaning and Disinfection
 - 5.1 General Information
 - 5.2 Cleaning and Disinfection
6. Troubleshooting and Device Error Codes
 - 6.1 Troubleshooting
 - 6.2 Device Error Codes
7. Technical Measurement Check
 - 7.1 Germany
 - 7.2 European Union outside of Germany
8. Technical Data
9. Electromagnetic Compatibility accompanying Documents according to EC 60601-1-2
10. Accessories
11. Disposal
12. Warranty


1. Introduction

1.1. Important Information read prior to Start-up

You have purchased a high quality Riester **ri-thermo**[®] fastPRObe Thermometer, which has been manufactured according to the Directive 93/42 EEC to the Council Directive 93/42/EEC and is subject to the strictest quality controls at all times. Read these instructions for use carefully before putting the unit into operation and keep them in a safe place. If you should have any questions, we are available to answer queries at all times. Our address can be found in these instructions for use. The address of our sales partner will be given upon request. Please note that all instruments described in these instructions for use are only to be used by suitably trained personnel. The perfect and safe functioning of this device is only guaranteed when original parts and accessories from Riester are used.

1.2. Safety Information and Electromagnetic Compatibility

Symbol	Symbol Note
	Follow the instructions in the operation manual. The symbol is printed in blue color on the device. The symbol is printed in black color on the probe cover box.
	Medical Device
	Type B applied part
IPX1	IP- Protection art Protection against vertically falling water (dripping water)
	Warning! The general warning sign indicates a potentially hazardous situation which could result in serious injury.
	Attention! Important note in this instruction use The attention symbol indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices
	Direct current
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician (licensed healthcare practitioner).
	Latex free
	DEHP free

	Non Sterile
	Single use
	Use by YYYY-MM / (Year-Month)
	Manufacturing date YYYY-MM-DD / (Year-Month-Day)
	Manufacturer
	Manufacturer serial number
	Lot number
	Reference number
	Unique Device Identifier
	Temperature for transport and storage condition
	Relative Humidity for transport and storage condition
	Air pressure for transport and storage Air pressure for Ambient Operating
	CE Mark
	Symbol for the marking of electrical and electronics devices according to Directive 2002/96/EC.
	Non-ionizing radiation
	Position
	Recall
	Timer

1.3 Packaging Symbols

Symbol	Symbol Note
	Fragile. Show transport package contents fragile, so handling should be handled with care.
	Beware the package from getting wet.
	Upward. It shows the correct position to transport the package.
	Keep away from sunlight
	„Grüner Punkt“ (country-specific)
	# Indicates the number of pieces in the package

The device satisfies the requirements for electromagnetic compatibility. Please note that under the influence of unfavorable field strengths, e.g. during the operation of wireless telephones or radiological instruments, adverse effects on function cannot be excluded.

The electromagnetic compatibility of this device has been verified by test according to the IEC 60601-1-2 requirements.

1.4 Intended Use

The **ri-thermo**[®] fastPRObe Thermometer is used to measure the temperature of the body in the oral (oral) anus (rectal) and in the armpit (axilla), thus detecting, diagnosing vital body function.

1.4.1 Indication

The **ri-thermo**[®] fastPRObe Thermometer is used for direct body temperature measurement in the mouth (orally) in the armpit (axillary), as well as in the anus (rectal) and thus for detection and diagnosis of vital body function.

1.4.2 Contraindication

The device is not designed, sold, or intended for use except as indicated.

1.4.3 Intended patient population

The device is intended for all patients.

1.4.4 Intended operators / users

The device could be used by a doctor, nurse in hospitals, medical facilities, clinics, doctors' offices. No use in MR environment!

1.4.5 Required skills / training of operators

The operator must have basic knowledge for measurement of vital signs. All connectors and connections are clearly explained in the instructions for use.

The user must comply exactly with the instructions in the instructions for use.

1.4.6 Environmental conditions

The device is determined to be used in a controlled environment .

The device must not be exposed to harsh / rough environmental conditions.

1.5 Warnings / Attention

	<p>Warnings!:</p> <p>The general warning sign indicates a potentially hazardous situation which could result in serious injury.</p>
	<p>No use in MR environment!</p>
	<p>There is a possible danger of inflammation of gases if the device is operated in the presence of inflammatory mixtures or mixtures of pharmaceuticals and air or oxygen or laughing gas!</p> <p>The appliance shall not be operated in premises where inflammatory mixtures or mixtures of pharmaceuticals and air or oxygen or laughing gas are present, e.g. Operating rooms.</p>
	<p>Electric shock!</p> <p>The ri-thermo[®] fastPRObe Thermometer housing may only be opened by authorized persons</p>
	<p>Damage to the device due to falling, or to high ESD influence!</p> <p>If the display information on the display is incomplete, the device must not be reused and must be returned to the manufacturer for repair.</p>
	<p>The device is to determine to be used in a controlled environment.</p> <p>The device must not be exposed to harsh / rough environmental conditions.</p>
	<p>ri-thermo[®] fastPRObe Thermometer probe and probe cover are designed for use only with this ri-thermo[®] fastPRObe Thermometer.</p>
	<p>Do not use this thermometer without first installing a ne ri-thermo[®] fastPRObe Thermometer probe cover.</p>
	<p>Use only ri-thermo[®] fastPRObe Thermometer probe covers with this device.</p>
	<ul style="list-style-type: none">- Only use the blue probe (oral, axillary) with the blue probe insert!- Only use the red probe (rectal) with the red probe insert!- To limit cross contamination, only use blue devices for taking oral and axillary temperatures.- Use red devices only for rectal temperatures.
	<p>Use of any other probe cover will result in erroneous temperature readings.</p>

	The device and probe covers are non-sterile. Do not use on injured tissue.
	To limit cross contamination, only use blue devices for taking oral and axillary temperatures. To aid infection control, use ever a separate probe cover box for blue and a separate for red probe.
	To limit cross contamination, only use red devices only for rectal temperatures. To aid infection control, use ever a separate probe cover box for blue and a separate for red probe.
	Do not open unit. No user-serviceable parts inside. Opening device may affect calibration and voids warranty.
	The batteries must be inserted in the correct position.

	Attention! The caution symbol indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.
	Thoroughly dry all electrical contacts on both probe and thermometer after cleaning or device may fail to function properly.
	For re-calibration, service or integrity checks refer to a qualified biomedical technician or return to manufacturer.
	Disposal of used probe covers must be performed in accordance with current medical practices or local regulations regarding disposal of infectious, biological medical waste.
	Removal of the batteries is recommended if the unit is not going to be used for an extended period of time.
	Disposal of batteries must be consistent with local environmental and institutional policy for NiMH or Alkaline battery disposal.
	Disposal of old battery-powered electronic equipment must be consistent with institutional policy for expired equipment disposal.
	Cleaning frequency and practices must be consistent with institutional policy for cleaning of non-sterile devices.

	<ul style="list-style-type: none"> - We recommend to remove the batteries of ri-thermo[®] fastPRObe Thermometer before cleaning or disinfection. - Take care when cleaning and disinfecting the Thermometer and the probe, so that no liquid penetrates into the interior. - Never place the ri-thermo[®] fastPRObe Thermometer or removable parts of the ri-thermo[®] fastPRObe Thermometer (probe, cables) in liquids! - The thermometer is shipped non-sterile. DO NOT use ethylene oxide gas, heat, autoclave or any other harsh methods to sterilize the unit. - The devices are not meant to undergo machine-processed maintenance and sterilization. This could lead to irreversible damage! - The probe cover is suitable for single use only!
	<p>The patient are not the intended Operator. The device could be used by trained personnel. Trained personnel is a doctor, nurse in hospitals, medical facilities, clinics, doctors' offices.</p>

1.6 User Responsibility

	<p>Attention! It is your responsibility:</p> <ul style="list-style-type: none"> - Before each use, the user must check the integrity and completeness of the thermometer, the probe with spiral cable and the probe cover. All components must be compatible with each other. <p>Incompatible components can result in degraded performance.</p> <ul style="list-style-type: none"> - Check calibration of the device every two years. - Never knowingly use a defective device. - Immediately replace parts that are broken, worn, missing, incomplete, damaged or contaminated. - If the probe is damaged, it must be replaced immediately with a new probe. The damaged probe must not be used any further! - Contact the nearest factory approved service center should repair or replacement become necessary. - Further, the user of the device bears sole responsibility for any malfunction that results from improper use, faulty maintenance, improper repair, damage or alteration by anyone other than Riester or authorized service personnel.
	<p>Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.</p>

2. Using the Device for the first time

2.1 Scope of Supply

<p>Artikel- No. 1830</p>	<p>ri-thermo[®] fastPRObe Thermometer Oral / Axillary Probe cover box (20 pc / Box) User Manual 4 AA Alkaline Batteries</p>
<p>Artikel- No. 1831</p>	<p>ri-thermo[®] fastPRObe Thermometer Rectal Probe cover box (20 pc / Box) User Manual 4 AA Alkaline Batteries</p>

2.2 Device Function

The **ri-thermo**[®] fastPRObe Thermometer are working in Direct Mode with a Fast Measuring Mode and a Direct Measuring Mode.

Measuring range 34° - 43°C (93,2°F to 109°F)

Fast Measuring Mode: Accuracy $\pm 0.3^{\circ}\text{C}$ ($\pm 0.5^{\circ}\text{F}$)

Direct Measuring Mode: Increased measurement accuracy $\pm 0.1^{\circ}\text{C}$ ($\pm 0.2^{\circ}\text{F}$)

Adjustment of body temperature measuring points

Timer: 60 seconds time (15, 30, 45, 60 sec.).

Switching: ° C in ° F

Recall: Calling up the last measured value.

Interchangeable thermometer probe insert with color coding:

- Blue for oral, axillary measuring, - Red for rectal measuring.

Replaceable thermometer probe with color coding:

- Blue for oral, axillary measuring, - Red for rectal measuring.

Power supply by means of four normal "AA" batteries.

Battery status indicator.


Automatic switch-on function when removing the probe from its compartment.

Automatic switch-off function.


Easy-to-read liquid crystal display (LCD) with backlight.

Easily accessible and replaceable probe box.


ri-thermo[®] fastPRObe Thermometer

Thermometer probe blue Oral / Axillary

Thermometer probe red Rectal

1. Recall button
2. Position button
3. Display
4. Thermometer probe
5. Thermometer probe insert
6. Thermometer probe connector
7. Timer button
8. Mount for probe cover box
9. Battery compartment
10. Button for probe cover release
11. Spiral cable
12. Thermometer probe connector
13. Thermometer probe tip
14. Probe cover

3. Information about Body Temperature

It is a common misconception that 37 °C is the 'normal' body temperature.

It is a mistake to assume that the "normal" body temperature is 37 °C. It's the average body temperature. Age, gender, measurement point and time of day determine normal body temperature. It is usually lower in the morning, higher in the afternoon and something again usually lower in the evening. Patient's particular activity, metabolic rate or medications taken are factors that affect body temperature include. With increasing age, the normal body temperature also tends to drop. Due to the different measuring points on the body, the measured temperatures are not directly comparable.





In the following table according to the age of the patient and measuring point the normal temperatures are listed.

Body temperature measuring points	Normal body temperatures according to patient age			
	Patient age	0-2 years	3-10 years	11-65 years
Mount / Oral measuring	- -	35,5 °C - 37,5 °C 95,9 °F - 99,5 °F	36,4 °C - 37,6 °C 97,6 °F - 99,6 °F	35,8 °C - 36,9 °C 96,4 °F - 98,5 °F
Armpit / Axillary measuring	34,7 °C - 37,3 °C 94,5 °F - 99,1 °F	35,9 °C - 36,7 °C 96,6 °F - 98,0 °F	35,2 °C - 36,9 °C 95,3 °F - 98,4 °F	35,6 °C - 36,3 °C 96,0 °F - 97,4 °F
Rectum / Rectal measuring	36,6 °C - 38,0 °C 97,9 °F - 100,4 °F	36,6 °C - 38,0 °C 97,9 °F - 100,4 °F	37,0 °C - 38,1 °C 98,6 °F - 100,6 °F	36,2 °C - 37,3 °C 97,1 °F - 99,2 °F
Ear	36,4 °C - 38,0 °C 97,5 °F - 100,4 °F	36,1 °C - 37,8 °C 97,0 °F - 100,0 °F	35,9 °C - 37,6 °C 96,6 °F - 99,7 °F	35,8 °C - 37,5 °C 96,4 °F - 99,5 °F
Heart	36,4 °C - 37,8 °C 97,5 °F - 100,0 °F	36,4 °C - 37,8 °C 97,5 °F - 100,0 °F	36,8 °C - 37,9 °C 98,2 °F - 100,2 °F	35,9 °C - 37,1 °C 96,9 °F - 98,8 °F


4. Operation and Function

4.1 Commissioning

4.1.1 Insert the Batteries

<p>1.  Attention! Please note sequence! [1] Push on the button! [2] Remove the battery cover!</p>	<p>2. [1] Remove the battery cover!</p>	<p>3.  Attention! [1] Please note the symbols for the correct insert of the batteries!</p>
		
<p>4.  Attention! [1] Please insert correct!</p> 	<p>5.  Attention! Please insert correct!</p> 	<p>6.  Attention! Please insert correct!</p> 
<p>7.  Attention! Please note sequence! [1] Put on the battery cover! [2] closed the battery cover!</p> 		
	<p>Dispose of batteries in a manner consistent with local environmental and institutional policy for NiMH or Alkaline battery disposal.</p>	


4.2 Icon Identification

	
	<p>Oral Measurement position Works only with BLUE probe</p>
	<p>Axillary Measurement position Works only with BLUE probe</p>
	<p>Rectal Measurement position Works only with RED probe</p>
	<p>8bit- LEDs Display probe serial number and other information</p>
	<p>Measurement progress bar Display measurement progress in fast mode</p>
	<p>Temperature unit</p>
	<p>Timer 60 seconds time</p>
	<p>Direct Measuring Mode Speed: slowly Accuracy: High</p>
	<p>Probe Flashing probe when the probe is wrong</p>
	<p>Install / remove probe cover Flashing probe cover reminder to install probe cover prior to use and remove probe cover after use</p>
	<p>High ambient temperature Ambient temperature is higher than the range of equipment. The symbol shows the user that the ambient temperature is outside the ambient temperature range of the thermometer. The device is to determine to be used in a controlled environment. The device must not be exposed to harsh / rough environmental conditions.</p>

	<p>Low ambient temperature Ambient temperature is lower than the range of equipment. The symbol shows the user that the ambient temperature is outside the ambient temperature range of the thermometer. The device is to determine to be used in a controlled environment. The device must not be exposed to harsh / rough environmental conditions.</p>
	<p>Battery Battery status indicator</p>
	<p>4bit- LEDs Display temperature and other information</p>
	<p>Below the range of measurement Measurement of low temperature</p>
	<p>Higher than measurement range Measurement of high temperature</p>

4.3 Probe covers – Applying & Removing

- 4.3.1 Open probe cover box by lifting tab at top corner and pulling to remove top panel.
- 4.3.2 Insert box of probe covers into the top of probe insert. (To aid infection control, use ever a separate probe cover box for blue and a separate for red probe).
- 4.3.3 Remove thermometer probe from the probe insert.
- 4.3.4 To help remind the user to apply or remove a probe cover, the probe cover icon will be displayed.
- 4.3.5 Insert the probe end into a probe cover in the box. Push the handle firmly until you feel the cover “snap” into place.
- 4.3.6 Take appropriate temperature measurement (oral, axillary or rectal).
- 4.3.7 Eject the used cover into bio-waste container by pressing top button.
- 4.3.8 Remove, discard and replace box when empty.

Warning!


- Do not use this thermometer without first installing a new **ri-thermo**® fastPRObe Thermometer probe cover.
- Use only **ri-thermo**® fastPRObe Thermometer probe covers with this device.
- Use of any other probe cover will result in erroneous temperature readings.
- The device and probe covers are non-sterile! Do not use on injured tissue.
- To limit cross contamination, only use blue devices for taking oral and axillary temperatures.
- To limit cross contamination, only use red devices only for rectal temperatures.
- To aid infection control, use ever a separate probe cover box for the blue and a separate for the red probe.


Attention!

Disposal of used probe covers must be performed in accordance with current medical practices or local regulations regarding disposal of infectious, biological medical waste.

4.4 Change Probes

<p>1. [1] Initial situation</p> 	<p>2. [1] Remove the probe</p> 	<p>3.  Attention! Please note sequence! [1] Push the button! [2] Remove the probe insert! [3] Remove the probe connector!</p> 
<p>4.  Attention! Please note sequence! [1] Putting the probe connector in! [2] Putting the probe insert in!</p>	<p>5.  Attention! [1] Putting the probe in!</p>	

Warning!

- Only use the blue probe (oral, axillary) with the blue probe insert!
- Only use the red probe (rectal) with the red probe insert!
- To limit cross contamination, only use blue devices for taking oral and axillary temperatures. To aid infection control, use always a separate probe cover box for the blue and a separate for the red probe.
- Use red devices only for rectal temperatures.

4.5 Fast Measuring Mode

- 4.5.1 The electronic thermometer normally operates in fast measuring mode to provide fast and accurate temperature measurements.
- 4.5.2 When the probe is removed out from the probe insert, the probe preheats to 33,5 °C (92,3°F).
- 4.5.3 To help remind the user to apply a probe cover, the probe cover icon  will be displayed and blinks for a few seconds. To remind the user to remove the probe cover after the measurement, the probe cover icon  remains displayed.
- 4.5.4 Insert the probe end into a probe cover in the box. Push the handle firmly until you feel the cover “snap” into place.
- 4.5.5 Take appropriate temperature measurement (oral, axillary or rectal).
- 4.5.6 The Measurement progress bar  shows the progress.
- 4.5.7 A “long beep” is sounded when measurement is complete and the final temperature is displayed.
- 4.5.8 After returning the probe to the probe insert, the temperature is stored. If the probe is returned to the probe insert before the “long beep” is heard, no temperature will be stored for the recall function
- 4.5.9 Fast measuring mode is always functional regardless of red or blue probe.

Attention!

An ▼ or ▲ will appear on the display whenever the current temperature reading is out of range, either high or low, respectively. However, in instances when no temperature is detected or the temperature does not stabilize, the thermometer will automatically switch to direct measuring mode

4.6 Direct Measuring Mode

- 4.6.1 The **ri-thermo**® fastPRObe Thermometer normally operates in fast measuring mode to provide fast and accurate temperature measurements. However, in instances when no temperature is detected or the temperature does not fast stabilize, the thermometer will automatically switch to direct measuring mode.
- 4.6.2 A snail icon 🐌 will be continuously displayed whenever the thermometer is functioning in direct measuring mode.
- 4.6.3 The Measurement progress bar ██████████ are running and shows
- 4.6.4 A reading is only displayed when the body temperature is stable within a range of ± 0.1 °C for 10 seconds.
- 4.6.5 A "long beep" is sounded when measurement is complete and the final temperature is displayed.
- 4.6.6 The direct measuring mode auto feature is always functional regardless of red or blue probe / probe insert.


Attention!

An ▼ or ▲ will appear on the display whenever the current temperature reading is out of range, either High or Low, respectively.

4.7 Oral & Axillary Temperature Taking

- 4.7.1 The Blue probe must be installed.
- 4.7.2 Remove the probe. The thermometer will turn on automatically.
- 4.7.3 The symbol oral or axillary mode shows on the display. Never the symbol rectal.
- 4.7.4 Press the position button on the front Panel. Now to select either the oral or axillary mode. Apply a probe cover.
- 4.7.5 For oral temperatures measuring, the probe tip must insert deep into the sublingual pocket next to the vertical fold of tissue in middle of tongue, (frenulum linguae), on the one side or the other side, toward the back of the mouth.
- 4.7.6 The Patient must be CLOSING the mouth.
- 4.7.7 Hold the probe Securely on the location until a temperature value will displayed.

4.7.8 For axillary temperatures measuring, patient has to raise the arm. Then place the probe tip gently in the axilla. Press gently to assure a good contact of the patients skin.

4.7.9 The patient must lower the arm and press firmly against his body as long as the measurement lasts. The probe must be held as shown in the picture.


4.7.10 A "long beep" is sounded when measurement is complete and the final temperature is displayed.

However, in instances when no temperature is detected or the temperature does not fast stabilize, the thermometer will automatically switch to direct measuring mode. A snail icon  will be continuously displayed whenever the thermometer is functioning in direct measuring mode. A "long beep" is sounded when measurement is complete and the final temperature is displayed.


4.7.11 Eject the used probe cover into a bio-waste container by pushing the top button on the probe.

4.7.12 After the measurement returning the probe to the probe insert, the temperature is stored for recall. After a new measuring erases last temperature memory.

1.

2.

4.8 Rectal Temperature Taking

4.8.1 Make sure that the red probe / probe insert unit is attached. A separate probe cover box for red probe are installed. Withdraw the probe and the Thermometer turns on automatically.

An icon identifying rectal mode is always displayed provided the red probe insert / probe assembly is attached. Pressing the position button on the front panel to change mode has no effect.

4.8.2 Apply a new probe cover.

4.8.3 Insert the probe into the patient's rectum. To ensure proper tissue contact, angle the probe slightly after insertion.

4.8.4 Depth of insertion is recommended at $\frac{1}{2}$ " to $\frac{3}{4}$ " (12 mm – 19 mm) for adults and $\frac{1}{4}$ " to $\frac{1}{2}$ " (6 mm – 13 mm) for children.

4.8.5 A "long beep" is sounded when measurement is complete and the final temperature is displayed.

However, in instances when no temperature is detected or the temperature does not fast stabilize, the thermometer will automatically switch to direct measuring mode. A snail icon  will be continuously displayed whenever the thermometer is functioning in direct measuring mode. A "long beep" is sounded when measurement is complete and the final temperature is displayed.

4.8.6 Eject the used probe cover into a bio-waste container by pushing the top button on the probe.

4.8.7 After the measurement returning the probe to the probe insert, the temperature is stored for recall.

After a new measuring erases last temperature memory.

4.9 Setup Mode

- 4.9.1 To enter setup mode, first press the „Timer“ button and hold it. Second, press and hold the „Position“ button. Hold both together for 4 seconds. A “beep” is heard.
- 4.9.2 The button „Timer“ is used to move forward through the different configuration option. „Recall“ button is used to change an option configuration.
- 4.9.3 In Setup mode, can query software version, change Celsius and Fahrenheit scales, display battery power, display the serial number of the probe.
- 4.9.4 After 30 seconds automatically exit the Setup mode, or move the probe out and back in again.

4.10 Recall last temperature

- 4.10.1 After each temperature measurement, a “long beep” is heard. The “beep” indicates the temperature measurement has been completed and stored and is available for recall. This temperature can be recalled after the probe is returned to the probe insert.
- 4.10.2 To recall the most recent temperature measurement, press and release button „Recall“ on the front panel. The last measurement and “Recall” taken will appear for several seconds.
- 4.10.3 After a new measuring erases last temperature memory.
If the most recent temperature measurement was incomplete or out of range, dashes will appear on the display during the recall operation.

4.11 Pulse timer mode

- 4.11.1 The **ri-thermo**[®] fastPRObe Thermometer may also be used to help measure a patient’s pulse rate.
- 4.11.2 To access, do not remove probe from the probe insert.
Press and release the button „Timer“ on the front panel. The clock icon and 0.0 will be displayed. Press and release the button „Timer“ a second time to start the timer and elapsed time display.
- 4.11.3 A “beep” is sounded at 15 seconds, 2 beeps at 30 seconds, 1 beep at 45 seconds and 2 beeps at 60 seconds. The count stops at 60 seconds. Pressing the button „Timer“ again will restart the counter.
- 4.11.4 The Timer mode will turn off automatically, 30 seconds after stopping at 60 or whenever the probe is removed from the probe insert.

5. Cleaning and Disinfection

5.1 General Information

Cleaning and disinfection of medical products protects patients, users and third parties, lead to value retention of medical products. Due to product design and materials used there is no possibility to define the maximum limit of re-processing cycles. The lifetime of a medical product is determined by its function and how it is used. Before sending back defective products for repair, the following instructions should be followed.

5.2 Cleaning and Disinfection

The **ri-thermo**[®] fastPRObe Thermometer can be cleaned on the outside (with the exception of the display glass cover) using a damp cloth until optical cleanliness is achieved. Use disinfection products only according to the manufacturer’s instructions. Only disinfectants with proven effectiveness according to national guidelines should be used. After disinfection, please wipe the instruments using a damp cloth in order to eliminate any remnants of the disinfectant.

	<p>Attention!</p> <ul style="list-style-type: none"> - We recommend to remove the batteries of ri-thermo[®] fast-PRObe Thermometer before cleaning or disinfection. - Take care when cleaning and disinfecting the thermometer and the probe, so that no liquid penetrates into the interior. - Never place the ri-thermo[®] fastPRObe Thermometer or removable parts of the ri-thermo[®] fastPRObe Thermometer (probe, cables) in liquids! - The thermometer is shipped non-sterile. DO NOT use ethylene oxide gas, heat, autoclave or any other harsh methods to sterilize the unit. - The devices are not meant to undergo machine-processed maintenance and sterilization. This could lead to irreversible damage! - The probe cover is only single use suitable!
	<p>For all reusable devices, if there are any signs of material degradation, the device should no longer be reused and should be disposed/claimed following the procedure mentioned under Disposal / Warranty.</p>

6. Troubleshooting and device error codes

6.1 Troubleshooting

- 6.1.1 If the **ri-thermo**[®] fastPRObe Thermometer does not function properly after cleaning / disinfecting and reconnecting the probe, clean the battery contacts and the contacts of the thermometer probe connector with a dry cloth.
- 6.1.2 If the **ri-thermo**[®] fastPRObe Thermometer breakdown at any time, it will auto enter to alarming mode, and the corresponding error codes will display.

6.2 Device error codes

- 6.2.1 P01 Probe miss
- 6.2.2 P02 Pt > 55°C
- 6.2.3 P03 HT > 55°C
- 6.2.4 P04 In 10s Pt or Ht can't be heated over 30°C
- 6.2.5 E05 Pt calibration error (use in calibration)
- 6.2.6 E06 HT calibration error (use in calibration)
- 6.2.7 E07 System timer error
- 6.2.8 E09 LCD error (We can only hear the alarm. The screen can't be displayed.)
- 6.2.9 E11 System flash error Icon display of ambient temperature too high (40°C) or too low (10°C)
- 6.2.10 E15 Dead Battery (below 4.6 V)

7. Technical measurement check

7.1 Germany

Pursuant to the ordinance on Operating Medical Devices (Medizinprodukte- Betriebsverordnung, MPBetreibV), a technical check must be conducted every two years. The check can only be conducted by the manufacturer, who is the appropriate authorities or persons who comply with the requirements of MPBetreibV §6.

7.2 European Union outside of Germany

The respective legal provisions apply for all EU countries outside of Germany. Technical measurement check for the **ri-thermo**[®] fastPRObe Thermometer: The probe accuracy can be checked using a calibrated water bath (the probe cover must be attached!). Thermometers which are out of the error threshold must be checked in direct measuring mode. The switch to this mode is done automatically when the probe is removed from the

probe insert and about 30 seconds is placed next to the Thermometer.
The measurement check time for the water bath testing last approx. 60 seconds.

	<p>Attention! We recommend a technical check in two-year intervals.</p>
--	---

7.3 Standard

This thermometer conforms to all of the requirements established in
ASTM Standard E1104
ASTM Standard E1112
DIN EN ISO 80601-2-56

8. Technical Data

Medical Equipment:	Electronic thermometer
Electrical protection:	Electrical protection:
Electrical degree of protection:	Type B
Operating mode:	Continuous manual operation
Dimensions:	205 x 70 x 70mm
Weight:	400g with Batteries and Probe cover box
Battery	4x AA, LR6 Alkaline 1,5 V or 4x AA Mignon NiMH 1,2 V
Temperature Measurement Range:	34 °C - 43 °C (93,2 °F to 109 °F)
Ambient Operating Temperature Range:	10 °C to 40 °C (50 °F to 104 °F) at 10% - 85% non-condensing
Ambient Operating Air pressure:	700 to 1050 hPa
Storage and transport conditions	-20 to 55 °C (-4 °F to 131 °F), at 10 - 95 % non-condensing
Storage and transport Air pressure	700 to 1050 hPa
Laboratory Accuracy Water Bath Accuracy (36°C to 42.0°C):	Direct Measuring Mode $\pm 0.1^{\circ}\text{C}$ ($\pm 0.2^{\circ}\text{F}$) Fast Measuring Mode $\pm 0.3^{\circ}\text{C}$ ($\pm 0.5^{\circ}\text{F}$)
Direct measuring mode Accuracy (All measuring points)	$\pm 0,1^{\circ}\text{C}$ ($\pm 0,2^{\circ}\text{F}$)
Fast measuring mode Accuracy (All measuring points)	$\pm 0,3^{\circ}\text{C}$ ($\pm 0,5^{\circ}\text{F}$)
Typical Measurement Times (After position of the probe)	Direct Measuring Mode Oral Mode: 60 seconds Axillary Mode: 60 seconds Rectal Mode: 60 seconds Fast Measuring Mode Oral Mode: 12 seconds Axillary Mode: 12 seconds Rectal Mode: 12 seconds
Probe Cover Material	HDPE

9. Electromagnetic Compatibility accompanying Documents according to IEC 60601-1-2

The instrument satisfies the requirements for electromagnetic compatibility. Please note that under the influence of unfavorable field strengths, e.g. during the operation of wireless telephones or radiological instruments, adverse effects on function cannot be excluded.

The electromagnetic compatibility of this device has been verified by test according to the IEC60601-1-2 requirements.

9.1 EMC (electromagnetic compatibility)

- 9.1.1 During installation and operation of the device, observe the following instructions:
- 9.1.2 Do not use the device simultaneously with other electronic equipment to avoid electromagnetic interference with the operation of the device.
- 9.1.3 Do not use or stack the device near, on, or under other electronic equipment to avoid electromagnetic interference with the operation of the device.
- 9.1.4 Do not use the device in the same room as other electronic equipment, such as life-support equipment that has major effects on the life of the patient and results of treatment, or any other measurement or treatment equipment that involves small electric current.
- 9.1.5 Do not use cables or accessories that are not specified for the device because that may increase the emission of electromagnetic waves from the device and decrease the immunity of the device to electromagnetic disturbance.

Attention!

Medical electrical equipment is subject to special precautions regarding electromagnetic compatibility (EMC). Portable and mobile radio frequency communication devices can affect medical electrical equipment. The ME device is for operation in an electromagnetic environment and intended for professional facilities such as industrial areas and hospitals. The user of the device should ensure that it is operated within such an environment.

Warning!

The ME device may not be stacked, arranged or used directly next to or with other devices. When operation is required to be close to or stacked with other devices, the ME device and the other ME devices must be observed in order to ensure proper operation within this arrangement. This ME device is intended for use by medical professionals only. This device may cause radio interference or interfere with the operation of nearby devices. It may become necessary to take appropriate corrective measures, such as redirecting or rearranging the ME device or shield. The rated ME device does not exhibit any basic performance features in the sense of IEC 60601-1, which would present an unacceptable risk to patients, operators or third parties should the power supply fail or malfunction.

Warning!

Portable RF communications equipment (radios) including accessories, such as antenna cables and external antennas, should not be used in closer proximity than 30 cm (12 inches) to parts and cables of the **ri-thermo**[®] fastPRObe Thermometer specified by the manufacturer. Failure to comply may result in a reduction in the device's performance features.

Directives and manufacturer's declaration - Electromagnetic emissions

The **ri-thermo**[®] fastPRObe Thermometer are intended for use in the electromagnetic environment specified below. The customer or user of the **ri-thermo**[®] fastPRObe Thermometer should ensure that it is used in such an environment.

Guidance and manufacture's declaration – electromagnetic emission

The **ri-thermo**[®] fastPRObe Thermometer intended for use in the electromagnetic environment specified below. The customer of the user of the **ri-thermo**[®] fastPRObe Thermometer should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The ri-thermo [®] fastPRObe Thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The ri-thermo [®] fastPRObe Thermometer is intended for use in all establishments, including residential areas and those directly connected to a public supply network that also supplies buildings used for residential purposes.
Harmonic emissions IEC 61000-3-2	N / A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N / A	

Guidance and manufacture's declaration – electromagnetic immunity

The **ri-thermo**[®] fastPRObe Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of **ri-thermo**[®] fastPRObe Thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Con: ± 8 kV Ais: ± 15 kV	Con: ± 8 kV Ais: ± 15 kV	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	N / A	N / A	The quality of the supply voltage should be that of a typical business or hospital environment.
Surge IEC 61000-4-5	N / A	N / A	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	N / A	N / A	The quality of the supply voltage should be that of a typical business or hospital environment.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Mains frequency magnetic fields should be at a level characteristic of a typical location in a typical commercial hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacture's declaration – electromagnetic immunity

The **ri-thermo**[®] fastPRObe Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the **ri-thermo**[®] fastPRObe Thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	Pass	Pass	<p>Portable and mobile RF communications equipment should not be used closer to any part of the ri-thermo[®] fastPRObe Thermometer, including the cables, than the recommended distance, which is calculated using the equation applicable to the transmitter frequency.</p> <p>Recommended separation distance $d = 1.2\sqrt{P}$ 150 KHz to 80 MHz $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2,7 GHz</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	
Proximity fields from RF wireless communications equipment	Pass	Pass	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> <p>Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and the recommended distance is given in meters (m). Field strengths from fixed RF transmitters determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of devices marked with the following symbol:</p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a.) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ri-thermo**[®] fastPRObe Thermometer is used exceeds the applicable RF compliance level above, the **ri-thermo**[®] fastPRObe Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the **ri-thermo**[®] fastPRObe Thermometer.

b.) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the **ri-thermo[®] fastPRObe Thermometer.**

The **ri-thermo**[®] fastPRObe Thermometer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **ri-thermo**[®] fastPRObe Thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **ri-thermo**[®] fastPRObe Thermometer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,7 GHz
	$d = 12\sqrt{P}$	$d = 1,2\sqrt{P}$	$d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

10. Accessories

Art.- No. 13334	Probe cover 1 pack of 20 pieces
Art.- No. 13338	Probe cover 25 x 1 pack of 20 pieces
Art.- No. 13335	Thermometer Probe Oral (blue) / Axillary, Thermometer Probe Insert Oral / Axillary (blue)
Art.- No. 13336	Thermometer Probe rectal (red), Thermometer Probe Insert rectal (red)

11. Disposal

	Attention! Disposal of the used medical device must be performed in accordance with current medical practices or local regulations regarding disposal of infectious, biological medical waste.
	Batteries and electrical/electronic devices must be disposed in accordance with locally applicable regulations, not with domestic waste.
	If you have any questions regarding disposal of products, please contact the manufacturer or its representatives.

12. Warranty

This product has been manufactured under the strictest quality standards and has undergone a thorough final quality check before leaving our factory. We are therefore pleased to be able to provide a warranty of 2 years from the date of purchase on all defects, which can verifiably be shown to be due to material or manufacturing faults. A warranty claim does not apply in the case of improper handling. All defective parts of the product will be replaced or repaired free of charge within the warranty period.

This does not apply to wearing parts. For R1 shock-proof, we grant an additional warranty of 5 years for the calibration, which is required by CE-certification. A warranty claim can only be granted if this warranty card has been completed and stamped by the dealer and is enclosed with the product. Please remember that all warranty claims have to be made during the warranty period. We will, of course, be pleased to carry out checks or repairs after expiry of the warranty period at a charge. You are also welcome to request a provisional cost estimate from us free of charge. In case of a warranty claim or repair, please return the Riester product with the completed warranty Card to the following address:

Rudolf Riester GmbH
Dept. Repairs RR
Bruckstr. 31
72417 Jungingen
Germany

Serial number or batch number

Date, Stamp and signature of the specialist dealer

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