EliteVue

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Table of Contents

| | |
|---|----|
| 1. Important information to consider before taking the product into operation | 9 |
| 2. Battery handles | 10 |
| 3. Taking the product into operation (inserting and removing disposable batteries and rechargeable batteries) | 10 |
| 4. Charging the battery handles with rechargeable batteries | 10 |
| 5. Attaching the instrument head | 11 |
| 6. EliteVue | 11 |
| 6.1 Intended purpose | 11 |
| 6.2 Attention / Contraindications | 12 |
| 6.3 Attaching and removing ear specula | 12 |
| 6.4 Focusing wheel | 12 |
| 6.5 Pneumatic test | 12 |
| 6.6 Technical data of the lamp | 12 |
| 6.7 Replacing the lamp | 12 |
| 7. Instruktionen for care | 12 |
| 7.1 General information | 12 |
| 7.2 Sterilization | 12 |
| 7.3. Cleaning and disinfection | 12 |
| 7.4. Spare parts and accessories | 12 |
| 7.5. Maintenance | 12 |
| 7.6. Notices | 12 |
| 7.7. Electromagnetic compatibility | 13 |

1. Important information to consider before taking the product into operation

You have purchased a high-quality **Riester** diagnostic instrument set which was manufactured in compliance with Directive 93/42 EEC and is subject to constant quality checks of the strictest kind. The outstanding quality will ensure that your diagnoses are reliable. In these instructions for use, the use of the **Riester** battery handles as well as the instrument head and its accessories are described. Read the instructions for use through carefully before taking the product into operation and keep them safely. Should you have any questions, please feel free to contact either us or your responsible representative for **Riester** products anytime. You can find our address on the last page of this manual. We will be glad to give you the address of our representative upon request. Please note that all of the instruments described in this manual are intended solely for use by appropriately trained personnel. Please also note that the correct and safe function of our instruments is only ensured if both the instruments used and their accessories are exclusively from **Riester**.

Safety notices:

-  Manufacturer
-  CE marking
-  Temperature limits in °C for storage and transport
-  Temperature limits in °F for storage and transport
-  Humidity limitation for storage and transport
-  Fragile, handle with care
-  Store in a dry place
-  „Grüner Punkt“ (country-specific)
-  Attention! Heed operating instructions!
-  Device of protection class II
-  Type B applied part
-  For single use only
-  Used electrical and electronic devices should not be disposed of with normal household waste, but should be disposed of separately in accordance with national EU guidelines.
-  Batch code
-  Serial number
-  Follow the operating instructions

2. Battery handles

2.1. Intended purpose / Indications

The **Riester** battery handles serve to supply the instrument heads with energy (the lamps are included in the respective instrument heads). They furthermore serve as holders.

2.2. Range of battery handles

The instrument head described in these instructions for use fits the following battery handles and can thus be individually combined.

This instrument head furthermore fits the handles of the **ri-former**® wall model.

PLEASE NOTE:

LED instrument heads are only compatible from a certain serial number of the **ri-former**® diagnostic station onwards. We will be glad to provide you with information on the compatibility of your diagnostic station upon request.

For EliteVue:

2.3. Type C battery handle with rheotronic®, 2.5 V.

To operate this battery handle, you require:

- 2 commercially available alkaline Type C Baby batteries (IEC-standard designation LR14).

2.4. Type C battery handle with rheotronic®, 3.5 V (for ri-charger® L)

To operate this battery handle, you require:

- 1 rechargeable battery from **Riester** with 3.5 V (Item No. 10691, **ri-accu**® L).
- 1 **ri-charger**® L charging unit (Item No. 10705, Item No. 10706)

2.5. Type C battery handle with rheotronic®, 3.5 V

for charging in a 230 V or 120 V wall socket

To operate this wall socket handle, you require:

- 1 rechargeable battery from **Riester** with 3.5 V (Item No. 10692, **ri-accu**® L).

2.6. Type C battery handle with rheotronic®, 3.5 V

(for plug charging unit) To operate this battery handle, you require:

- 1 rechargeable battery from **Riester** with 3.5 V (Item No. 10694, **ri-accu**® L).
- 1 plug charging unit (Item No. 10707).

2.7. Type AA battery handle with rheotronic®, 2.5 V. rheotronic®

To operate this battery handle, you require:

- 2 commercially available alkaline Type AA Mignon batteries (IEC-standard designation LR6)

2.8. Type AA battery handle with rheotronic®, 3.5 V (ri-charger® L)

To operate this battery handle, you require:

- 1 rechargeable battery from **Riester** with 3.5 V (Item No. 10690, **ri-accu**® L).
- 1 **ri-charger**® L charging unit (Item No. 10705, Item No. 10706)

3. Taking the product into operation (inserting and removing disposable batteries and rechargeable batteries)

PLEASE NOTE:

Only use the combinations described in Sections 2.3 to 2.8!

3.1. Inserting the batteries:

Type C and Type AA battery handles (2.3 and 2.7) with **rheotronic**® 2.5 V :

- Twist off the cap of the battery handle on the lower part of the handle by turning the cap counterclockwise.
- Insert the commercially available alkaline batteries required for this battery handle into the battery handle with the plus side facing toward the upper part of the handle. Screw the cap of the battery handle firmly back onto the battery handle.

3.2. Removing the batteries:

Type C and Type AA battery handles (2.3 and 2.7) with **rheotronic**® 2.5 V:

- Twist off the cap of the battery handle on the lower part of the battery handle by turning the cap counterclockwise.
- Remove the batteries from the battery handle by holding the opening of the battery handle slightly downward and shaking a little if necessary. Screw the cap of the battery handle firmly back onto the battery handle.

PLEASE NOTE:

The red safety foil on the plus side must be removed from all rechargeable batteries already inserted at the factory or separately included in delivery before taking the battery into operation for the first time!

PLEASE NOTE:

Only applies to Type C battery handle (2.5) with **rheotronic**®, 3.5 V for charging in a 230 V wall socket 120 V: When using the new **ri-accu**® L, Item No. 10692, care must be taken to ensure that no insulation is affixed to the spring of the battery-handle cap. If the old **ri-accu**® L, Item No. 10692, is used, an insulation must be affixed to the spring (risk of short circuit!).

new **ri-accu**® L

old **ri-accu**® L

3.3. Inserting the rechargeable batteries:

Type C and Type AA battery handles (2.4 and 2.8) with **rheotronic**®, 3.5 V (for **ri-charger**® L).

Type C battery handle (2.6) with **rheotronic**®, 3.5 V (for plug charging unit).

Type C battery handle (2.5) with **rheotronic**®, 3.5 V for charging in a 230 V or 120 V wall socket.

CAUTION! Heed safety notice!

- Twist off the cap of the battery handle on the lower part of the handle by turning the cap counterclockwise.
- Remove the red safety foil on the plus side of the rechargeable battery when taking the battery into operation for the first time.
- Insert the rechargeable battery that is permissible for your battery handle (see 2.2) into the battery handle with the plus side facing toward the upper part of the handle.

Beside the plus sign, you will additionally see an arrow that shows you the direction for insertion into the battery handle. Screw the cap of the battery handle firmly back onto the battery handle.

3.4. Removing the rechargeable batteries:

Type C and Type AA battery handles (2.4 and 2.8) with **rheotronic**[®], 3.5 V (for **ri-charger**[®] L).
Type C battery handle (2.6) with **rheotronic**[®], 3.5 V (for plug charging unit).
Type C battery handle (2.5) with **rheotronic**[®], 3.5 V for charging in a 230 V or 120 V wall socket.

CAUTION! Heed safety notice!

- Twist off the cap of the battery handle on the lower part of the battery handle by turning the cap counterclockwise.
- Remove the rechargeable battery from the battery handle by holding the opening of the battery handle slightly downward and shaking a little if necessary. Screw the cap of the battery handle firmly back onto the battery handle.

4. Charging the battery handles with rechargeable batteries

4.1. Type C and Type AA battery handles (2.4 and 2.8) with **rheotronic**[®], 3.5 V (for **ri-charger**[®] L).

- Can only be charged with charging unit **ri-charger**[®] L (Item No. 10705, Item No. 10706) from **Riester**.
- A separate instruction manual is enclosed with the **ri-charger**[®] L charging unit. This must be heeded.

4.2. Type C battery handle (2.6) with **rheotronic**[®], 3.5 V (for plug charging unit).

- This can only be charged with the plug charging unit (Item No. 10707) from **Riester**. To do this, the small round plug on the lower part of the battery handle is inserted into the rechargeable battery (Item No. 10694 **ri-accu**[®] L) via the opening in the cap of the battery handle. Now connect the mains plug of the plug charging unit to the mains supply. The state of charge of the rechargeable battery is indicated by the LED on the plug charging unit. Red light signifies that charging is in progress, green light signifies that the rechargeable battery is fully charged.

4.3. Type C battery handle (2.5) with **rheotronic**[®], 3.5 V for charging in a 230 V or 120 V wall socket.

- Twist off the lower part of the wall socket handle by turning it counterclockwise. The wall socket contacts appear. Round contacts are for 230 V mains operation, flat contacts are for 120 V mains operation. Now insert the lower part of the handle into the wall socket for charging.

PLEASE NOTE:

Before the wall socket handle is taken into operation for the first time, it should be charged for up to 24 hours in the wall socket.

CAUTION!

The wall socket handle must not be charged for longer than 24 hours.

CAUTION!

The handle must never ever be inserted in the wall socket while the rechargeable battery is being replaced!

Technical data:

Optionally 230 V or 120 V version

PLEASE NOTE:

- Should you not use the device for a longer period of time, or should you take it along when travelling, remove the disposable batteries and rechargeable batteries from the handle.
- New batteries should be inserted when the light intensity of the instrument becomes weaker.
- To get an optimal light yield, we recommend that you always insert new high-quality batteries when replacing batteries (as described in 3.1 and 3.2).
- If it is suspected that fluid or moist condensation may have penetrated into the handle, the handle must under no circumstances be charged. Particularly in the case of wall socket handles, this can lead to a life-threatening electric shock.
- In order to prolong the service life of the rechargeable battery, the battery should only be charged when the light intensity of the instrument becomes weaker.

4.4. Disposal:

Please note that disposable batteries and rechargeable batteries must be disposed of in a specific way. You can obtain information in this regard from your local community or from your responsible environmental advisor.

5. Attaching the instrument head

Fit the instrument head onto the receptacle on the upper part of the handle in such a way that the two recesses of the lower part of the instrument head are seated on the two projecting guide cams of the battery handle. Press the instrument head gently onto the battery handle and turn the handle clockwise up to the stop position. To remove the instrument head, turn it counterclockwise.

5.1. Switching on and off

Type C and Type AA battery handles. Switch the instrument on by tapping the switch ring on the upper part of the handle in clockwise direction. To switch the instrument off, press the ring in counterclockwise direction until the device switches off.

5.2. **rheotronic**[®] for adjusting the light intensity

The light intensity can be adjusted on the Type C and Type AA battery handles via the **rheotronic**[®] system. The light intensity becomes weaker or stronger, depending on how often you tap the switch ring in counterclockwise or clockwise direction.

CAUTION!

The light intensity is always 100% when the battery handle is switched on. Automatic safety switch-off after 180 seconds. **Explanation of the sign on the wall socket handle:**

 Attention! Heed operating instructions!

6. EliteVue

6.1. Intended purpose

The **Riester** EliteVue device was produced to illuminate and examine the auditory canal and the tympanic membrane in combination with the **Riester** ear specula.

6.2. Attention / Contraindications

There may be a risk of ignition of gases when the instrument is used in the presence of combustible mixtures or mixtures of pharmaceuticals. The product and the ear funnels are non-sterile. Do not use on injured tissue. Use new or disinfected ear funnels to limit the risk of cross-contamination. Disposal of used ear funnels must be done in accordance with current medical practices or local regulations regarding the disposal of infectious, biological medical waste. Use only accessories / consumables from Riester or Riester approved accessories / consumables. Cleaning frequency and flow must comply with the regulations for the cleaning of non-sterile products in the respective facility. Cleaning and disinfection instructions of the instructions for use must be observed. The product may only be used by trained personnel.

6.3. Attaching and removing ear specula

To equip the otoscope head, single-use ear specula from **Riester** (in black), or alternatively reusable ear specula from **Riester** (in black) can be selected. The size of the ear speculum is indicated at the back of the speculum. The specula from Welch Allyn can also be used.

L1 and L2 specula

Turn the speculum clockwise until you can feel a resistance. To remove the speculum, twist off the speculum by turning it counterclockwise.

L3 speculum

Fit the selected speculum onto the chrome-plated metal holder of the otoscope until you can feel it snap into place. To remove the speculum, press the blue eject button. The speculum is ejected automatically.

6.4. Focusing wheel

The focusing wheel enables you to adjust the focusing range in the auditory canal (tympanic membrane). Thanks to the high-quality lens system, you have 5.5-fold magnification and a field of view with a diameter of 11 mm at a distance of approx. 20mm.

6.5. Pneumatic test

In order to be able to perform the pneumatic test (= an examination of the tympanic membrane), you need a ball which is not included in the normal scope of delivery, but which can be additionally ordered. The tube of the ball is fitted to the connection. You can then carefully introduce the required amount of air into the auditory canal.

6.6. Technical data of the lamp

Otoscope XL 2,5 V 750 mA mean life span 15h

Otoscope XL 3,5 V 720 mA mean life span 15h

Otoscope LED 2,5 V 300 mA mean life span 10,000h

Otoscope LED 3,5 V 300 mA mean life span 10,000h

6.7. Replacing the lamp

6.8. EliteVue

Twist off the instrument head from the battery handle. The lamp is located in the lower part of the instrument head. Withdraw the lamp from the instrument head with your thumb and index finger or with a suitable tool. Firmly fit the new lamp in place.

7. Instructions for care

7.1. General information

Cleaning and disinfection of the medical devices serves to protect the patient, the user and third parties and to preserve the value of the medical devices. Due to the product design and the materials used, no defined limit can be specified for the maximum number of reprocessing cycles that can be carried out. The life span of the medical devices is determined by their function and by gentle handling of the devices. Defective products must undergo the reprocessing procedure described before being returned for repair.

7.2. Cleaning and disinfection

The instrument head with handle can be cleaned externally with a moist cloth until visually clean. Wipe-disinfection as specified by the disinfectant manufacturer. Only disinfectants with proven efficacy should be used, taking into account the national requirements. After disinfection, wipe the instrument down with a moist cloth to remove possible disinfectant residues.

CAUTION!

- Never immerse the instrument head with handle in any liquid! Take care to ensure that no liquids get inside the casing!
- This item is not approved for automated reprocessing and sterilization. These procedures cause irreparable damage!

7.3. Sterilization

a) Reusable ear specula

The ear specula can be sterilized in the steam sterilizer at 134°C with 10 minutes hold time.

b) Single-use ear specula

For single use only

Warning: Repeated use can cause infection.

7.4. Spare parts and accessories

For a detailed list, refer to our catalog „Instruments for E.N.T.“ which you can download from www.riester.de.

7.5. Maintenance

The instruments and their accessories do not require any special maintenance. If the case should arise that an instrument needs to be checked for whatever reason, please send it to us or to an authorized **Riester** specialist dealer in your area, whose details we will be glad to let you know upon request.

7.6. Notices

Ambient temperature:

0° to +40°

Relative humidity:

30% to 70%, non-condensing

Transport and storage temperature:

-10° to +55°

Relative humidity:

10% to 95%, non-condensing

Barometric pressure:

800 hPa - 1100 hPa

7.7. Electromagnetic compatibility

Medical electrical equipment is subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high-frequency communication equipment can influence medical electrical equipment. This ME device is intended for operation in an electromagnetic environment as specified below. The user of the device should ensure that it is operated in such an environment. This ME device must not be used directly next to or arranged in a stack with other devices. If the device has to be operated near to or in a stacked arrangement with other devices, then the ME device should be monitored in order to verify that it operates as intended in this arrangement. This ME device is intended exclusively for use by professional medical staff. This device can cause radio interference or disrupt the operation of equipment nearby. Suitable remedial measures, such as for instance re-alignment, re-arrangement of the ME device or shielding, may become necessary.

| Guidelines and manufacturer's declaration - electromagnetic emissions | | |
|---|----------------|---|
| The EliteVue device is intended for operation in an electromagnetic environment as specified below. The customer or the user of the EliteVue device should ensure that it is used in such an environment. | | |
| Emission measurements | Compliance | Electromagnetic environment - guidelines |
| HF emissions according to CISPR 11 | Group 1 | The EliteVue device employs HF energy solely for an internal function. Therefore, its HF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| HF emissions according to CISPR 11 | Class B | The EliteVue device is intended for use in all facilities, including living quarters and such as are directly connected to a public power supply that also supplies buildings that are used for residential purposes. |
| Harmonics emissions according to EC61000-3-2 | Not applicable | |
| Voltage-fluctuation/flicker emissions according to IEC61000-3-3 | Not applicable | |

| Guidelines and manufacturer's declaration - electromagnetic immunity | | | |
|---|--|--|--|
| The EliteVue device is intended for operation in an electromagnetic environment as specified below. The customer or the user of the EliteVue device should ensure that it is used in such an environment. | | | |
| Immunity tests | IEC 60601 test level | Compliance level | Electromagnetic environment - guidelines |
| Electrostatic discharge (ESD) according to IEC61000-4-2 | ± 6 kV contact discharge ± 8 kV air discharge | ± 6 kV contact discharge ± 8 kV air discharge | Floors should be made of wood or concrete or be covered with ceramic tiles. If the floor is covered with a synthetic material, then the relative humidity must be at least 30%. |
| Fast transient electrical interference/bursts according to IEC61000-4-4 | ± 2 kV for power lines ± 1 kV for input and output lines | Not applicable | The quality of the supply voltage should correspond to that of a typical business or hospital environment. |
| Surges IEC61000-4-5 | ± 1 kV voltage phase-to-phase ± 2 kV voltage phase-to-earth | Not applicable | The quality of the supply voltage should correspond to that of a typical business or hospital environment. |
| Voltage dips, short-time interruptions and fluctuations in the supply voltage according to IEC61000-4-11 | <-5% U_r (→95% drop in U_r) for 0.5 cycles 40% U_r (60% drop in U_r) for 5 cycles 70% U_r (30% drop in U_r) for 25 cycles <-5% U_r (→95% drop in U_r) for 5 s | Not applicable | The quality of the supply voltage should correspond to that of a typical business or hospital environment |
| Magnetic field at the mains frequency (50Hz) according to IEC61000-4-8 | 3 A/m | 3 A/m | If image disturbances occur, the EliteVue device may have to be placed further away from the sources of mains-frequency magnetic fields, or magnetic shielding may have to be installed; the mains-frequency magnetic field should be measured at the intended set-up site in order to ensure that it is small enough. |
| Note - U_r is the alternating supply voltage prior to application of the test level. | | | |

| Guidelines and manufacturer's declaration - electromagnetic immunity | | | |
|---|---|--|---|
| The EliteVue model is intended for operation in the electromagnetic environment specified below. The customer or the user of the EliteVue device should ensure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidelines |
| <p>Conducted HF interference according to IEC61000-4-6</p> <p>Radiated HF interference according to IEC61000-4-3</p> | <p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5GHz</p> | <p>Not applicable</p> <p>10 V/m</p> <p>3 V/m</p> | <p>Portable and mobile radio equipment should not be used within a distance from the EliteVue device, including cables, that is less than the recommended safety distance as calculated by the equation that is appropriate for the transmission frequency.</p> <p>Recommended safety distance:</p> <p>$d = 1.2\sqrt{P}$</p> <p>$d = 1.2\sqrt{P}$ 80 MHz to 1000 MHz</p> <p>$d = 2.3\sqrt{P}$ 1000 MHz to 2.5 GHz</p> <p>where P is the nominal power of the transmitter in Watts [W] as specified by the manufacturer of the transmitter, and d is the recommended safety distance in meters (m).</p> <p>The field strength of stationary radio transmitters should be less than the compliance level^a at all frequencies as verified by an on-site test^a</p> <p>Interference is possible in the vicinity of equipment marked with the following symbol</p>  |
| <p>Note 1: At 80 MHz and 800 MHz, the higher value applies.</p> <p>Note 2 These guidelines may not apply in all situations. The propagation of electromagnetic waves is influenced by reflection and absorption by buildings, objects and people.</p> | | | |
| <p>a The field strength of stationary transmitters, such as base stations of wireless telephones and mobile field radio services, amateur radio stations, AM and FM radio and television transmitters, cannot be precisely determined theoretically in advance. In order to determine the electromagnetic environment due to stationary HF transmitters, an investigation of the location is advisable. If the field strength determined at the location of the EliteVue device exceeds the compliance level indicated above, then the EliteVue device must be monitored with regard to its normal operation at each place where it is used. If unusual performance characteristics are observed, additional measures such as re-alignment of the EliteVue device or its removal to another place may be necessary.</p> <p>b In the frequency range of 150 kHz to 80 MHz, the field strength should be less than 3 V/m.</p> | | | |

| Recommended safety distances between portable and mobile HF communication devices and the EliteVue device | | | |
|--|---|--------------------|---------------------|
| The EliteVue device is intended for operation in an electromagnetic environment in which the radiated HF interference is monitored. The customer or the user of the EliteVue device can help prevent electromagnetic interference by observing minimum distances between portable and mobile HF communication equipment (transmitters) and the EliteVue device as recommended below, depending on the maximum output power of the communication equipment. | | | |
| Nominal power of the Transmitter W | Safety distance that applies to the transmitter frequency m | | |
| | 150 KHz to 80 MHz | 80 MHz to 1000 MHz | 1000 MHz to 2.5 GHz |
| | Not applicable | $d = 1.2\sqrt{P}$ | $d = 2.3\sqrt{P}$ |
| 0.01 | | 0.12 | 0.23 |
| 0.1 | | 0.38 | 0.73 |
| 1 | | 1.2 | 2.3 |
| 10 | | 3.8 | 7.3 |
| 100 | | 12 | 23 |
| For transmitters whose nominal power is not indicated in the table above, the distance can be determined using the equation belonging to the respective column, where P is the nominal power of the transmitter in Watts [W] as specified by the manufacturer of the transmitter. | | | |
| Note 1: At 80 MHz and 1000 MHz, the distance for the higher frequency range applies. | | | |
| Note 2 These guidelines may not apply in all situations. The propagation of electromagnetic waves is influenced by reflection and absorption by buildings, objects and people. | | | |

WARRANTY

This product has been manufactured under the strictest quality standards and has undergone a thorough final quality check before leaving our factory. We are therefore pleased to be able to provide a warranty **of 2 years from the date of purchase** on all defects, which can verifiably be shown to be due to material or manufacturing faults. A warranty claim does not apply in the case of improper handling. All defective parts of the product will be replaced or repaired free of charge within the warranty period. This does not apply to wearing parts. For R1 shock-proof, we grant an additional warranty of 5 years for the calibration, which is required by CE-certification. A warranty claim can only be granted if the Warranty Card has been completed and stamped by the dealer and is enclosed with the product. Please remember that all warranty claims have to be made during the warranty period. We will, of course, be pleased to carry out checks or repairs after expiry of the warranty period at a charge. You are also welcome to request a provisional cost estimate from us free of charge. In case of a warranty claim or repair, please return the **Riester** product with the completed warranty card to the following address:

Serial number of batch number
Date, Stamp and signature of the specialist dealer

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Rudolf Riester GmbH