

CARDIOVIT CS-200 Excellence

Designed for high-volume workload, equipped with the latest technology to save time and boost productivity



Medical Device Depot Inc.

www.medicaldevicedepot.com
1-877-646-3300

CARDIOVIT CS-200 EXCELLENCE

Experience the **CARDIOVIT CS-200 Excellence** and get a view of the future where precision and performance converge.

Thanks to its new processor – 3rd Generation Core i7 CPU – the **CS-200 Excellence** is now one of the fastest and most efficient devices on the market.

Especially designed for enterprise hospitals, the **CS-200 Excellence** masters each task.



ONE-BUTTON OPERATION

Immediate access to the desired function.

STAND-ALONE OR CONNECTED

The CS-200 Excellence optimizes digital workflow with seamless connectivity to EMR, PACS, and SCHILLER's SEMA3 or SEMA-200 Cardiology Information System. Whether used as a stand-alone system or a complete enterprise solution, it is a scalable solution that fits your needs today and tomorrow.



AUTONOMOUS EMERGENCY ECG

In an emergency every second counts. SCHILLER offers the unique ability to print an emergency ECG even before the system boots up.

1 Intuitive user interface

- Monitor ST changes with auto-comparison of current and reference beats.
- Watch for ST level and slope changes using the 12 or 16*-lead average complexes.

2 All wires are integrated internally

3 SCHILLER 12/16-lead thermal printer

4 Large writing surface for clinical documentation


5 Two big accessory drawers

6 System base cover for external printer. Integrated power isolation transformer providing a stable work platform



EASY DATA ENTRY

AUTOMATIC SpO₂, BP SpO₂ and blood pressure are automatically recorded during exercise testing and sent to the CS-200 Excellence. The BP-200 plus and the SCHILLER K-Sound Analysis (SKA) allow unparalleled accurate blood pressure measurements, with both auscultatory and oscillometric techniques.



EASY DATA ENTRY

Data entry as easy and fast as possible. Simply enter the data via bar code scanner.



BOOST PRODUCTIVITY

Thanks to rich on-screen information content, with customizable display layout, integrated resting ECG during Exercise Testing, Wizard, Dynamic Protocol and much more, our system saves time and boosts productivity.

UNIQUE DYNAMIC PROTOCOL

The system suggests the most appropriate exercise protocol, based on patient age, body weight and questionnaire. No modifications during the test are necessary.¹

WIZARD


The WIZARD is an intuitive tutorial which guides the user one step at a time, highlighting the next available function. It is designed for nurses, technicians and all other clinicians to perform the exercise test without a long preparation period.

ECG REVIEW Full Disclosure data enables review of every beat and arrhythmia for enhanced clinical confidence. Even during the exercise it is possible to compare ECG strips from previous stages or even previous tests.

MULTI-VIEW HOOK-UP SCREEN Quickly and easily attach lead wires:

- Colour-coded electrode placement diagram
- Dynamic "noise level" triangle
- Real-time ECG waveform display all help identify hook-up problems before you start the test

WIZARD PLANNER Plan your wizard steps according to your needs by assigning the actions that should be carried out.



¹ J Cardiopulm Rehabil. 2006 Jan-Feb;26(1):16-23. A nomogram to select the optimal treadmill ramp protocol in subjects with high exercise capacity: validation and comparison with the Bruce protocol. Maeder M, Wolber T, Afefy R, Gadza M, Ammann P, Myers J, Rickli H.

Electrode Test

Electrodes

| Lead | Offset |
|---------|---------|
| R (R) | -8 mV |
| L (L) | -11 mV |
| F (F) | -11 mV |
| C1 (C1) | 2678 mV |
| C2 (C2) | -2 mV |
| C3 (C3) | -2 mV |
| C4 (C4) | -2 mV |
| C5 (C5) | -8 mV |
| C6 (C6) | -8 mV |

Lead configuration: Standard

Legend: ▲ weak signal ▲ electrode off

Use reduced electrode set OK

SDS-200 Excellence [Exercise ECG] - User: default / Systemadministrator - 0001:James Smith

File View Function Settings Suppress Ergo Device Screen Help

HR 90 170 53% BP -- / -- ST -1.2 3.0

Recording: Standard Filter: none

Stage / Step: Warm-up 00:00:25

Exercise: 25W

Expected Load: 176 W

METs: Overall Time: 00:03:36

Average METs(i) 125%

Take Reference Auto STmax 3.40

1.2 mm, 0.0 mV/s

06.05.2014 15:21:05 Duration: recording...

Filter: none

10 mm/mV 150 mm/s Page 7 of 12

Start Pause Hold Next Recover Stop

The wizard will guide you through the planned steps by highlighting the buttons accordingly.

SDS-200 Excellence [Patient Selection] Resting ECG - User Exercise ECG Protocol Start Recovery Full Recovery Stop Recording

SUPERIOR DIAGNOSTIC TOOLS

HyperQ™: analysis of high-frequency QRS signals

Sophisticated filters: industry-leading analysis program for artefact and noise management

16-lead acquisition and analysis: in resting and exercise ECGs




HYPER Q™

Unique, scientifically and clinically proven analysis of high-frequency QRS signals, which significantly improves the accuracy and effectiveness of routine exercise and resting ECG

- ❖ More sensitive diagnosis of ischaemic patients
- ❖ Fewer false alarms in non-ischaemic patients
- ❖ Gender-independent accuracy
- ❖ Determined diagnosis in patients with equivocal exercise and resting ECG
- ❖ Reduces health care costs by saving unnecessary expensive diagnostic tests
- ❖ Prevents unnecessary exposure to radioactive or invasive procedures


Exercise ECG: Patient with abnormal ST and negative Hyper Q

75 years, male, CAD: ST segment analysis shows no abnormalities. However, the HyperQ analysis correctly detected significant ischaemia.




Resting ECG: Patient with positive Hyper Q and ACS (Acute Coronary Syndrome)

62 years, female, ACS: as displayed above, ST segment analysis shows no abnormalities. However, the HyperQ analysis correctly detects significant ischaemia.




SAECG (LATE POTENTIAL ANALYSIS)

This analysis is a non-invasive alternative to invasive testing for ventricular arrhythmia. It allows detection of micropotentials that occur after the QRS complex. Gain time thanks to straightforward measurement data recording; pushing one button is enough.



VECTOR 3D

3D view of the cardiac electrical activity, and a complete diagnosis of frontal and posterior walls.



SOPHISTICATED, ADAPTIVE MULTIPURPOSE FILTER

SCHILLER offers an industry-leading analysis program for artefact and noise management. The RNS filter provides the physician with excellent ECG quality and accurate measurements without overfiltering and without compromising clinical significance.

SCHILLER 16-LEAD ECG ACQUISITION AND ANALYSIS

While in the past resting ECG analysis programs were limited to 12 simultaneous leads, the SCHILLER 16-lead algorithm analyses now up to 16 leads of simultaneously acquired ECG waveforms to provide an interpretation of rhythm and morphology for a wide variety of

patient populations. The algorithm reflects newly updated guidelines and recommendations, such as the 2007 AHA/ACCF/HRS Recommendations Part II¹, and the 2009 AHA/ACCF/HRS Recommendations Part VI² for the Standardization and Interpretation of the ECG.

The SCHILLER 16-lead algorithm goes beyond traditional 12-lead interpretation of the resting ECG. It also provides for incremental diagnostic capabilities not associated with analysis programs of the past. The 16-lead integrated analysis takes advantage of optional right chest and back electrodes to provide extended interpretations for adult chest pain.

In Exercise ECG, the use of right precordial leads (16-lead ECG) along with the standard six left precordial leads greatly improves the sensitivity of exercise testing for the diagnosis of coronary artery disease³.

1 AHA/ACCF/HRS Recommendations for the Standardization and Interpretation of the Electrocardiogram, Part II: Electrocardiography Diagnostic Statement List. J Am Coll Cardiol, 2007;49:1128-135.

2 AHA/ACCF/HRS Recommendations for the Standardization and Interpretation of the Electrocardiogram, Part VI: Acute Ischemia/Infarction. Circulation 2009; 100:e262-e270.

3 N Engl. Med 1999;340:340-5.