
**A Guide to the Use of Ophthalmoscopes
in the Eye Examination**

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The Ophthalmoscope

Transparency of the cornea, lens and vitreous humor permits the physician to directly view arteries, veins, optic nerve and the retina.

Direct observation of the structures of the fundus through an ophthalmoscope may show disease of the eye itself or may reveal abnormalities indicative of disease elsewhere in the body. Among the most important of these are vascular changes due to diabetes or hypertension and swelling of the optic nerve head due to papilledema or optic neuritis. In this sense, the eye serves as a window through which many valuable clinical evaluations may be made.

When a preliminary diagnosis of an imminently dangerous eye condition, such as acute glaucoma or retinal detachment, is made by the examiner, prompt referral to an ophthalmologist may prevent irreversible damage. Or, when distressing but less urgent conditions, such as visual impairment due to cataract or vitreous floater are recognized, the patient can be reassured and referred.

- A. Front surface mirror
- B. Crossed linear polarizing/red-free filter switch
- C. Aperture selection disc
- D. On/Off rheostat control
- E. Rubber brow rest
- F. Lens selection disc
- G. Illuminated lens indicator

Acknowledgment

We wish to express our sincere appreciation to the Scheie Eye Institute and Dr. Steven Koenig, Dr. Ralph Eagle, Dr. Ken Spitzer and John Griffin for their contribution to this booklet.

Standard apertures

There is a wide range of practical apertures to select from: micro spot, small spot, large spot, fixation, slit and cobalt blue filter. A red-free filter is also available for use on the apertures. This selection of apertures covers all the physician's needs in an ophthalmologic examination.

A. Micro spot aperture: Allows quick visual entry in very small, undilated pupils.

B. Small aperture: Provides easier view of fundus through an undilated pupil. Always start examination with this aperture, proceed to large aperture as pupil adapts to light.

C. Large aperture: Standard aperture for dilated pupil and general examination of the eye.

D. Fixation aperture: The pattern of an open center and thin lines permits easy observation of eccentric fixation without masking the macula. The graduated cross hairs can be used to estimate either the amount of eccentric fixation relative to the macula or the size or location of a lesion on the retina or choroid.

NOTE: When viewed outside the patient at a distance less than 40", the fixation aperture will be out of focus. The lens of the eye will insure correct focus on the fundus.

E. Slit or streak: Helpful in determining various levels of lesions, particularly tumors and edematous discs.

F. Cobalt blue filter: When fluorescein dye is injected into the vein of a patient, the physician can observe the movement of this fluid within the vessels. When viewed through the cobalt filter of the ophthalmoscope, the dye will appear as a yellow/green color. If any vessel is abnormal, leaking or hemorrhaging, the cobalt filter will reveal this problem. Fluorescein drops in the eye may also help detect corneal abrasions and other lesions.

Other filters

Welch Allyn ophthalmoscopes #11720 and #11730 are equipped with a unique sliding switch that greatly increases their versatility.

Red-free filter: When the switch is positioned to the left (while facing the instrument front) it will be beneath a green dot and the red-free filter will be in place. This can be used in conjunction with any aperture. The red-free filter excludes red rays from the examination field; this is superior to ordinary light in viewing slight alterations in vessels, minute retinal hemorrhages, ill-defined exudates and obscure changes in the macula. The nerve fibers become visible and the observer may note the disappearance of such fibers, as in optic nerve atrophy. The background appears gray, the disc appears white, the macula appears yellow, the fundus reflex is intense and the vessels appear almost black. This filter is also used to help distinguish veins from arteries...veins stay relatively blue, but oxygenated arterial blood makes arteries appear blacker. This makes differentiation easier for the examiner.


Crossed linear polarizing filter: When the switch is positioned to the right (while facing the instrument front) it will be beneath a white circle with crosshairs inside. The crossed linear polarizing filter will be in place. This filter is used to eliminate corneal glare and reflection and can be used with any aperture. For more information on this filter turn to page 9.

Additional uses for the ophthalmoscope

In addition to examination of the fundus, the ophthalmoscope is a useful diagnostic aid in studying other ocular structures. The light beam can be used to illuminate the cornea and the iris for detecting foreign bodies in the cornea and irregularities of the pupil. By placing the +15.00 lens in the scope and looking at the pupil as in fundus examination [2 inches (5cm) distance from patient], the physician may verify doubtful pupillary action.


The practitioner can also easily detect lens opacities by looking at the pupil through the +6 lens setting at a distance of 6 inches (15cm) from the patient. In the same manner, vitreous opacities can be detected by having the patient look up and down, to the right and to the left. Any vitreous opacities will be seen moving across the pupillary area as the eye changes position or comes back to the primary position.

The Eye

With the ophthalmoscope 2 inches (5cm) in front of the eye, the lenses in the rotating wheel produce clear vision at points indicated in the diagrammatic eye.

The hyperopic or far-sighted eye requires more “plus” sphere for clear focus and the myopic or near-sighted eye requires “minus” sphere for clear focus.

- | | |
|----------------------------|--|
| A. Macula | J. Cornea |
| B. Vitreous humor | K. Ciliary Body |
| C. Sclera | L. Zonule (Suspensory Ligament) |
| D. Choroid | M. Conjunctiva |
| E. Retina | N. Lens |
| F. Ora Serrata | O. Hyaloid canal |
| G. Canal of Schlemm | P. Central retinal vein |
| H. Anterior chamber | Q. Optic nerve |
| I. Iris | R. Central retinal artery |

How to Conduct an Ophthalmologic Examination

Position the ophthalmoscope about 6 inches (15cm) in front and 25° to the right side of the patient. (Step 5)

In order to conduct a successful examination of the fundus, the examining room should be either semi-darkened or completely darkened. It is preferable to dilate the pupil when there is no pathologic contraindication, but much information can be obtained through the undilated pupil.

The following steps will help the physician obtain satisfactory results:

1. For examination of the right eye, sit or stand at the patient's right side.
2. Select "0" on the illuminated lens dial of the ophthalmoscope and start with small aperture.
3. Take the ophthalmoscope and start in the right hand and hold it vertically in front of your own right eye with the light beam directed toward the patient and place your right index finger on the edge of the lens dial so that you will be able to change lenses easily if necessary.
4. Dim room lights. Instruct the patient to look straight ahead at a distant object.
5. Position the ophthalmoscope about 6 inches (15cm) in front and slightly to the right (25°) of the patient and direct the light beam into the pupil. A red "reflex" should appear as you look through the pupil.

6. Rest the left hand on the patient's forehead and hold the upper lid of the eye near the eyelashes with the thumb. While the patient holds his fixation on the specified object, keep the "reflex" in view and slowly move toward the patient. The optic disc should come into view when you are about 1½ to 2 inches (3-5cm) from the patient. If it is not focused clearly, rotate lenses into the aperture with your index finger until the optic disc is as clearly visible as possible. The hyperopic, or far-sighted, eye requires more "plus" (black numbers) sphere for clear focus of the fundus; the myopic, or near-sighted, eye requires "minus" (red numbers) sphere for clear focus.

7. Now examine the disc for clarity of outline, color, elevation and condition of the vessels. Follow each vessel as far to the periphery as you can. To locate the macula, focus on the disc, then move the light approximately 2 disc diameters temporally. You may also have the patient look at the light of the ophthalmoscope, which will automatically place the macula in full view. Examine for abnormalities in the macula area. The

red-free filter facilitates viewing of the center of the macula, or the fovea.

8. To examine the extreme periphery, instruct the patient to:
 - a) look up for examination of the superior retina
 - b) look down for examination of the inferior retina
 - c) look temporally for examination of the temporal retina
 - d) look nasally for examination of the nasal retina.This routine will reveal almost any abnormality that occurs in the fundus.
9. To examine the left eye, repeat the procedure outlined above except that you hold the ophthalmoscope in the left hand, stand at the patient's left side and use your left eye.

CAUTION:

Before actuating the red-free filter/crossed linear polarizing filter slide switch, pull the instrument away from the patient's face to prevent contact with finger or switch.

Examine the disc for clarity of outline, color, elevation and condition of the vessels. (step 7)

Overcoming corneal reflection

One of the most troublesome barriers to a good view of the retina is the light reflected back into the examiner's eye by the patient's cornea — a condition known as corneal reflection.

1. On the ophthalmoscope shown in this book, the crossed linear polarizing filter may be used. This filter reduces the corneal reflection by 99%. To engage, simply move the switch on the front of the instrument to the position underneath the white crosshairs. It is recommended that the crossed linear polarizing filter be used when corneal reflection is present.
2. Use the small spot aperture. However, this reduces the area of the retina illuminated.
3. Direct the light beam toward the edge of the pupil rather than directly through its center. This technique can be perfected with practice.

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
Normal fundus

- Disc:** Outline clear; physiological cup is pale area centrally
Retina: Normal red/orange color, macula is dark avascular area temporally
Vessels: Arterial — venous ratio 2 to 3; the arteries appear a bright red, the veins a slightly purplish color

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
Hypertensive retinopathy (advanced malignant)

- Disc:** Elevated, edematous disc; blurred disc margins
- Retina:** Prominent flame hemorrhages surrounding vessels near disc border
- Vessels:** Attenuated retinal arterioles

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
Proliferative diabetic retinopathy

- Disc:** Net of new vessels growing on disc surface
Retina: Numerous hemorrhages, new vessels at superior disc margin
Vessels: Dilated retinal veins

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Papilledema

- Disc:** Elevated, edematous disc, blurred disc margins; vessels engorged
Retina: Flame retinal hemorrhage close to disc
Vessels: Engorged tortuous veins

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
Benign choroidal nevus

Retina: Slate gray, flat lesion under retina; several drusen overlying nevus

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Retinal detachment

Disc: Normal

Retina: Gray elevation in temporal area with folds in detached section

Vessels: Tortuous and elevated over detached retina

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Welch Allyn, Inc.
U.S.A.