WelchAllyn®

Rectal Light Handle

REF 73211

6V Power Supply

REF 73305, 73322, 73324, 73326

English	1-1
Français	2-1
Espanól	3-1
Português	4-1
Italiano	5-1
Deutsch	6-1
Dansk	7-1
Nederlands	8-1
Svensk	9-1
Norsk	10-1

Copyright 2009, Welch Allyn, Inc. All rights are reserved. No one is permitted to reproduce or duplicate, in any form, this manual or any part thereof without permission of Welch Allyn.

Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal or improper use of the product that may result from failure to use this product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

Welch Allyn® is a registered trademark of Welch Allyn. Welch Allyn is protected under various patents and patents pending.

Printed in USA

Table of Contents

Symbol Descriptions	1-1
Warnings and Cautions	1-1
Instructions for Use	1-3
REF 73305, 73322, 73324, 73326 (Power Supplies)	1-3
REF 73211 Rectal Light Handle	1-4
Maintenance	1-7
Preventative Maintenance	1-7
Lamp Replacement	1-7
Cleaning	1-8
Disposal	1-10
Specifications	1-11
Approvals	1-12
Guidance and Manufacturer's Declaration	1-13
Trouble Shooting Guide	1-17
Troubleshooting Guide - Qualified Personnel Only	1-18
Warranty and Service	1-19
Limited Warranty	1-19
Technical Assistance	1-19
Welch Allyn Service Centers	1-21

Thank you for purchasing the Welch Allyn REF 73211 Rectal Light Handle and/or the REF 73305, 73322, 73324, 73326 power supply. Please review this manual prior to using this system.

Only properly trained physicians for endoscopic type examinations should use the 73211 Rectal Light Handle and its appropriate Power Supply when coupled with the Welch Allyn devices listed in the table on page 1-6.

The 73211 Rectal Light Handle is designed for connection only with the power source and appropriate country voltage and outlet style listed within this manual. No other use of these products is intended or advised. Please follow the guidelines set forth in this manual for proper connection and use of the Rectal Light Handle.

Classification

Class II, Type BF Applied Part, Ordinary Equipment

Intermittent Operation

Symbol Descriptions

Attention, Warning, Caution:
See instructions for use.

Type BF Applied Part

Caution, high temperatures

Direct Current

Caution, high-intensity light

Alternating Current

Double Insulation

Not protected from the
ingress of water.

Warnings and Cautions

WARNING: To avoid eye injury, do not look directly into the light beam.

WARNING: Lamp is harmful to skin if touched when hot. Allow lamp to cool a minimum of 5 minutes before replacing.

WARNING: Rectal Light Handle and/or Power Supply is not suitable for use in the presence of a flammable anesthetic mixture with air or oxygen or nitrous oxide. An explosion may result.

WARNING: Lamp generates heat and should remain off until just prior to and immediately after use. Do not place Rectal Light Handle onto a flammable surface if lamp is illuminated or hot.

WARNING: Disconnect power supply from outlet prior to any maintenance.

WARNING: Only use the 6V power supplies described within this manual with the 73211 Rectal Light Handle. Failure to observe this warning may cause damage to the power supply and/or non-approved device

WARNING: The use of accessories, transducers, and/or cables other than those specified, with the exception of those sold by the manufacturer as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

CAUTION: Rectal Light Handle Illumination is intense.

CAUTION: This endoscopic accessory is for use with the Welch Allyn devices listed in the table on page 1-6 only. Do not use the illuminator for any other use.

CAUTION: Use only the cleaning methods described in this manual.

CAUTION: There are no user-serviceable parts inside the device except the lamp. Contact the authorized Service Center listed on page 1-21 for service related issues.

CAUTION: The 73211 Rectal Light Handle is designed for connection only with the power source and appropriate country voltage and outlet style listed within this manual. No other use of these products is intended or advised.

CAUTION: Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information starting on page 1-13.

CAUTION: Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.


Instructions for Use

The Rectal Light Handle's power supply (REF 73305, 73322, 73324, 73326) provides power from an appropriate wall outlet to the Rectal Light Handle (REF 73211) which houses the lamp (REF 07800). The Rectal Light Handle and its appropriate power supply provide an illumination system that is intended for use with Welch Allyn devices and/or accessories (see table on page 1-6).

REF 73305, 73322, 73324, 73326 (Power Supplies)

Verify that the power supply will connect to the appropriate outlet and that the power supplied to the outlet is within the range specified in the table provided.				
REF	Nom. Input Voltage	Input Current at Nom. Input Voltage (mA)	Input Voltage Range (V.)	Plug Style
73305	120V 60 Hz	71	110 - 130	A
73322	230V 50 Hz	30	210 - 230	B
73324	240V 50 Hz	27	230 - 250	C
73326	234V 50 Hz	27	230 - 250	D

Outlet Type


Insert the Power Supply into the receptacle of an appropriate voltage and plug configuration.

Once properly connected, the green LED located on the power supply housing will illuminate indicating power is present to the power supply cord.


REF 73211 Rectal Light Handle

Attach the 73211 Rectal Light Handle to the appropriate device selected from the table on page 1-6.

To assemble/disassemble the Rectal Light Handle System:

1. Select the appropriate device from the table on page 1-6.
2. Insert the 73211 Rectal Light Handle over the device's rectangular "input port" while twisting the 73211's knob in the "release" direction as shown. Release the knob and verify that the 73211 Rectal Light Handle is securely attached to the device.
3. Plug the REF 73305, 73322, 73324, or 73326 Power Supply into an appropriate outlet.
4. Align the plug at the end of the Rectal Light Handle cord to the connector at the end of the Power Supply and connect. The lamp in the 73211 will turn on immediately.
5. To turn the lamp off, disconnect the Rectal Light Handle from the Power Supply.
6. To remove the device from the Rectal Light Handle, twist the 73211's knob in the "release" direction and pull the 73211 off of the device's input port.

The Rectal Light Handle (REF 73211) and power supply (REF 73305, 73322, 73324, 73326) are designed to work with the following Welch Allyn devices for endoscopic procedures:

Welch Allyn 73211 Rectal Handle			
32 Series		38/39 Series Fitted with 38700	
32010	11mm x 10cm Pediatric FO Sig	38108	8 mm x 89mm Std Ano Spec
32020	11mm x 25cm Stricture FO Sig	38114	14mm x 89mm Std Ano Spec
32410	15mm x 15cm Pediatric FO Sig	38119	19mm x 89mm Std Ano Spec
32420	15mm x 25cm FO Sig	38122	22mm x 89mm Std Ano Spec
32810	19mm x 15cm Std FO Sig	38900	20mm x 68mm Rotating Ano Spec
32820	19mm x 25cm Std FO Sig	38614	14mm x 89mm Operating Ano Spec
33830	19mm x 30cm Std FO Sig	38619	19mm x 89mm Operating Ano Spec
32830	10mm x 35cm Std FO Sig	38622	22mm x 89mm Operating Ano Spec
33220	23mm x 25cm Std FO Sig	39614	14mm x 127mm Long Ano Spec
33620	27mm x 25cm Std FO Sig	39619	19mm x 127mm Long ano Spec
37 Series		39622	22mm x 127mm Long Ano Spec
37019	19mm x 7cm FO Anoscope	Disposables Fitted with 36019	
37023	23mm x 7cm FO Anoscope	53130	19mm x 25cm KleenSpec Disp Sig
37027	27mm x 7cm FO Anoscope	53110	Disposable Anoscope


Maintenance

Preventative Maintenance

Examine cords, connections, and transformer case for wear or damage monthly.

After every use, clean the Rectal Light Handle according to the Cleaning Section of this manual.

Lamp Replacement


NOTE: The lamp is the only user-serviceable part on the Rectal Light Handle (REF 73211). Refer to the above figure as a visual aid to replace the lamp.


1. Disconnect the system power supply from the wall outlet and then disconnect the Rectal Light Handle from the power supply cord. Allow lamp to cool for a minimum of 5 minutes.
2. Unscrew the Rectal Light Handle outer sleeve counter clockwise (CCW) to expose the lamp as shown.
3. Unscrew the lamp CCW from the Rectal Light Handle body.
4. Inspect the lamp socket of the Rectal Light Handle for damage or irregularities. Contact the authorized Service Center listed on page 1-21 if any irregularities are found.
5. Using a new Welch Allyn lamp, verify the lamp number is REF 07800 and that the o-ring is present and in good condition.
6. Screw the lamp clockwise (CW) into the Rectal Light Handle body and tighten sufficiently to compress the lamp o-ring.
7. Screw the Rectal Light Handle sleeve CW onto the Rectal Light Handle body and tighten sufficiently to secure to the handle.

8. Connect the power supply to the appropriate wall outlet and reconnect the Rectal Light Handle cord to the power supply cord and verify that the lamp is illuminated.
9. If the lamp fails to illuminate, refer to the Troubleshooting section on page 1-17 or contact your local Service Center (page 1-21) for further assistance.

Cleaning

Part Numbers 73210 and 73211
(including 07800-U lamp)

Part Number 73500
(including 03100-U lamp)


WARNING Allow the lamp to cool for at least 5 minutes before handling.

NOTE Leave the lamp installed while reprocessing 73210, 73211 and 73500.

Assembly and Disassembly Instructions for 73210 and 73211

1. Disconnect the light handle from the power supply cord and the device.
2. Unscrew the outer sleeve of the handle by turning it counterclockwise.

3. Follow one of the reprocessing methods described below.
4. After performing a reprocessing procedure, wipe the lamp with a soft, dry, lint free cloth.
5. Re-attach the clean sleeve to the handle by screwing in a clockwise direction.
6. Re-connect the power supply cord to the light handle and verify that the unit works by plugging it in. If the lamp fails to illuminate follow the troubleshooting instructions located in light handle manual (REF 700791) supplied with your handle.

Cleaning (73210, 73211 and 73500)

1. Prepare an enzymatic detergent formulated for endoscopic instruments according to the manufacturer's instructions. These instructions were validated using ENZOL/CIDEZYME.
2. Submerge the item in the cleaning solution.
3. Using a soft bristle brush, scrub the item submerged in the cleaning solution for a minimum of 5 minutes then rinse with sterile water for 5 minutes.

Intermediate Level Disinfection (73210, 73211 and 73500)

1. Follow the Cleaning instructions above.
2. Using CaviWipes™, wipe all item surfaces and keep wet with CaviWipe solution for a minimum of 5 minutes.
3. Allow the item surface to dry for a minimum of 10 minutes.

Sterilization (73211 ONLY)

Caution Autoclaving the light handle may reduce the product's functional life.

1. Follow the Cleaning instructions above.
2. Perform one autoclave cycle that includes 6 minutes at 270° F (132° C) and a cool down period of at least 15 minutes. These instructions were validated using a gravity autoclave. Validate any other sterilization equipment and load configurations.

Disposal

There are no harmful materials associated with the Rectal Light Handle System.

There are no risks associated with the disposal of this device after its useful life if it is cleaned per the instructions in this manual.

Specifications

Classification:	Class II, Type BF Applied Part
Input:	100 \sim to 240V \sim , 50/60 Hz, 200 mA
Output:	6.0V \equiv , 1150 mA
Category:	Not AP/APG Equipment
Physical Specifications:	
Rectal Light Handle:	96mm (3.775") x 25mm (1.00"), 104 grams (3.67 oz)
Power Supply:	43mm (1.71") L x 52mm (2.06") W x 55mm (2.18") H
Power Supply Cord:	183 cm (71.4")
Illumination:	
Lamp Life:	75 hours
Lamp Voltage:	6 V
Lamp:	Welch Allyn REF 07800
IPX Rating:	IPXØ
Environment	
Operating:	+10° C to +35° C (50° F to 95° F), 95% RH 500 hPa to 1060 hPa
Transport/Storage:	-20° C to +49° C (-2° F to 120° F), 95% RH, 500 hPa to 1060 hPa
 Intermittent Operation Rectal Light Handle	10 minutes on max; 5 minutes off min

Approvals

CERTIFIED TO: CAN/CSA STD C22.2 NO. 601.1

CONFORMS TO: EN/IEC 60601-1, UL 60601-1,
EN/IEC 60601-1-2, EN/IEC 60601-2-18

The CE mark on this product indicates that it has been tested to and conforms with the provisions noted within the 93/42/EEC Medical Device Directive.

European Regulatory Manager

European Regulatory Manager
Welch Allyn Ltd.
Navan Business Park
Dublin Road
Navan, County Meath, Republic of Ireland
Tel.: +353 46 90 67700 • Fax: +353 46 90 67756

Australia AS/NZS 3200.1

Guidance and Manufacturer's Declaration

Emissions for all EQUIPMENT and SYSTEMS

The REF 73211 and its powering accessories is intended for use in the electromagnetic environment specified below. The customer or user of the REF 73211 and its powering accessories should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The REF 73211 and its powering accessories uses RF energy for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The REF 73211 and its powering accessories is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 31000-3-3	Complies	

Electromagnetic Immunity for all EQUIPMENT and SYSTEMS

The REF 73211 is intended for use in the electromagnetic environment specified below. The customer or user of the REF 73211 should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
ESD IEC 61000-4-2	± 6kV Contact ± 8kV Air	Not Applicable, See Note 1	Floors should be wood, concrete, or ceramic tile. If floors are synthetic material, the relative humidity should be at least 30%.
EFT IEC 61000-4-4	±2kV Mains ±1kV I/Os	Not Applicable, See Note 1	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, Short Interruptions, and Voltage Variations on power supply input lines. IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	During each dip in the mains voltage the light will flicker. During the 5-Second dropout, the light shuts off, but returns when mains power is restored. See guidance at the right.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the REF 73211 and its powering accessories requires that the light remain on during the power outage, the system must be powered from an uninterruptible power supply or a battery.
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	3 A/m	Not Applicable, See Note 1	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<p>Note 1: The REF 73211 and its power accessories are a simple power supply and a lamp. The only EMC test events that could affect the performance or environment of the system are Harmonics, Flicker, Radiated and Conducted Emission, Surge Immunity, and Voltage Dips, ESD, Radiated RFI, EFT, Conducted RFI, and Power Frequency Magnetic Immunity should not affect the performance of this product; at the levels required in EN60601-1-2:2001.</p>			

Guidance and Manufacturer's Declaration - Emissions

The REF 73211 and its powering accessories is intended for use in the electromagnetic environment specified below. The customer or user of the REF 73211 and its powering accessories should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80MHz (outside ISM bands)</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>Not Applicable, See Note 1</p> <p>Not Applicable, See Note 1</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the REF 73211 and its powering accessories, including cables, than the recommended separation distances calculated below:</p> $d = \left(\frac{3.5}{V_1}\right) \sqrt{P}$ $d = \left(\frac{3.5}{E_1}\right) \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left(\frac{7}{E_1}\right) \sqrt{P}$ <p>800 Mhz to 2.5 GHz</p> <p>Where P is the max power in watts and d is the recommended separation distance in meters.</p> <p>Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V_1 and E_1). Interference may occur in the vicinity of equipment containing a transmitter.</p>
<p>Note 1: The REF 73211 and its power accessories are a simple power supply and a lamp. The only EMC test events that could affect the performance or environment of the system are Harmonics, Flicker, Radiated and Conducted Emission, Surge Immunity, and Voltage Dips, ESD, Radiated RFI, EFT, Conducted RFI, and Power Frequency Magnetic Immunity should not affect the performance of this product; at the levels required in EN60601-1-2:2001.</p>			

Recommended Separations Distances between portable and mobile RF Communications equipment and the REF 73211 and its powering accessories

The REF 73211 and its powering accessories is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the REF 73211 and its powering accessories can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment (transmitters) and the REF 73211 and its powering accessories as recommended below, according to the maximum output power of the communications equipment.

Max Output Power (Watts)	Separation (m) 150kHz to 80MHz $d = \left(\frac{3.5}{V_1}\right) \sqrt{P}$	Separation (m) 80 to 800MHz $d = \left(\frac{3.5}{E_1}\right) \sqrt{P}$	Separation (m) 800MHz to 2.5GHz $d = \left(\frac{7}{E_1}\right) \sqrt{P}$
0.01	Not Applicable	Not Applicable	Not Applicable
0.1	Not Applicable	Not Applicable	Not Applicable
1	Not Applicable	Not Applicable	Not Applicable
10	Not Applicable	Not Applicable	Not Applicable
100	Not Applicable	Not Applicable	Not Applicable

Trouble Shooting Guide

Problem	Possible Cause	Solution
Lamp does not light.	Rectal Light Handle is not plugged into the power supply.	Verify that the Rectal Light Handle is properly plugged into power supply cord.
	Power Supply is not plugged in.	Verify the power supply is properly plugged into the outlet.
	Faulty Lamp	Replace with a new lamp.
	Loose Lamp	Verify that the lamp is secured into the Rectal Light Handle.
	Wrong Lamp	Verify that the lamp is a REF 07800 Welch Allyn lamp.
	Damaged Lamp	Inspect the lamp and the o-ring for damage. Replace with REF 07800 lamp if damage is found.
Rectal Light Handle does not stay plugged into the transformer cord.	Rectal Light Handle/ transformer interconnect plug is damaged.	Inspect interconnect plug for damage. Replace either or both REF 73211 or REF 73305, 73322, 73324, 73326.
Light output seems low.	Faulty Lamp	Replace with a new lamp.
	Dirty lamp lens	Inspect lamp lens for debris and clean.
This guide is intended to provide assistance towards determining the fault condition of this product. Please contact the authorized Welch Allyn Service Center listed within the product owners manual for additional assistance.		

Troubleshooting Guide - Qualified Personnel Only

Problem	Possible Cause	Solution
Lamp does not light.	No power to the power supply outlet.	Verify that power is supplied to power supply outlet.
	Faulty Power Supply.	Test the voltage supplied through the power supply at the cord end. Verify the voltage reads approximately 6.0 volts DC (REF 73305, 73322, 73324, 73326). Only properly trained personnel shall perform this test.
	Faulty Rectal Light Handle.	Verify that the Rectal Light Handle is faulty by replacing with a new Rectal Light Handle or with power supplied through the transformer (above), remove lamp and measure the voltage between the lamp socket center contact and socket threads. Verify voltage reads approximately 6.0 volts (DC with REF 73305, 73322, 73324, 73326). Only properly trained personnel shall perform this test.
Light output seems low.	Improper voltage supplied to the transformer.	Verify the voltage supplied to the transformer outlet meets the requirements stated in this manual.
This guide is intended to provide assistance towards determining the fault condition of this product. Please contact the authorized Welch Allyn Service Center listed within the product owners manual for additional assistance.		

Warranty and Service

Limited Warranty

Instrumentation purchased new from Welch Allyn, Inc. (Welch Allyn) is warranted to be free from original defects in material and workmanship under normal use and service for a period of one year for Welch Allyn Rectal Light Handle (REF 73211) and the Power Supply (REF 73305, 73322, 73324, 73326) from the date of first shipment from Welch Allyn. This warranty shall be fulfilled by Welch Allyn or its authorized representative repairing or replacing at Welch Allyn's discretion, any such defect, free of charge for parts and labor.

Welch Allyn must be notified via telephone of any defective product and the item must be immediately returned, with a Service Notification number provided by Welch Allyn, securely packaged and postage prepaid to Welch Allyn. Loss or damage in shipment shall be at the purchaser's risk.

Welch Allyn will not be responsible for loss associated with the use of any Welch Allyn product that (1) has had the date of manufacture number defaced, (2) has been repaired by anyone other than an authorized Welch Allyn Service Representative, (3) has been altered, (4) has been used in a manner other than in accordance with the instructions, or (5) has been abused or exposed to extreme environmental conditions.

THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ANY IMPLIED WARRANTY OR MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE, OR OTHER WARRANTY OF QUALITY, WHETHER EXPRESSED OR IMPLIED. WELCH ALLYN WILL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM THE SALE, USE, OR IMPROPER FUNCTIONING OF THE INSTRUMENTATION REGARDLESS OF THE CAUSE. THE DAMAGES FOR WHICH WELCH ALLYN WILL BE RESPONSIBLE INCLUDE, BUT ARE NOT LIMITED TO, LOSS OF REVENUE OR PROFIT, DOWNTIME COSTS, AND LOSS OF USE OF THE INSTRUMENTATION.

Technical Assistance

If you have an equipment problem that you cannot resolve, you may call the Welch Allyn Service Center nearest you for assistance. Technical service support is available by telephone on normal business days at the phone numbers listed under Welch Allyn Service Centers.

Prior to contacting Welch Allyn or its service center, please refer to the following areas of the illuminator and/or transformer to identify the date of manufacture of the component(s).

- Rectal Light Handle: Four digits, month then year, stamped next to the REF 73211 label on the Rectal Light Handle sleeve.
- Power Supply: Four digits, week then year stamped or labeled on the prong side of the transformer plug.

If you are advised to return a product to Welch Allyn for repair or routine maintenance, schedule the repair with the service center nearest you.

Before returning a product for repair you must obtain authorization from Welch Allyn. Service Notification number will be given to you by our service personnel. Be sure to note this number on the outside of your shipping box. Returns without a Service Notification number will not be accepted for service. Prior to shipment, please clean and disinfect the component as instructed in this manual.

Welch Allyn Service Centers

For information about any Welch Allyn product, call Welch Allyn Technical Support:

USA	+1 800 535 6663 + 1 315 685 4560	Australia	+61 2 9638 3000 +800 074 793
Canada	+1 800 561 8797	China	+86 21 6327 9631
European Call Center	+353 46 90 67790	France	+33 1 55 69 58 49
Germany	+49 7477 9271 70	Japan	+81 3 3219 0071
Latin America	+1 305 669 9003	Netherlands	+31 157 505 000
Singapore	+65 6419 8100	South Africa	+27 11 777 7555
United Kingdom	+44 207 365 6780	Sweden	+46 85 853 6551

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Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80MHz (outside ISM bands)</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>Not Applicable, See Note 1</p> <p>Not Applicable, See Note 1</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the REF 73211 and its powering accessories, including cables, than the recommended separation distances calculated below:</p> $d = \left(\frac{3.5}{V_1} \right) \sqrt{P}$ $d = \left(\frac{3.5}{E_1} \right) \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left(\frac{7}{E_1} \right) \sqrt{P}$ <p>800 Mhz to 2.5 GHz</p> <p>Where P is the max power in watts and d is the recommended separation distance in meters.</p> <p>Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V_1 and E_1). Interference may occur in the vicinity of equipment containing a transmitter.</p>
<p>Note 1: The REF 73211 and its power accessories are a simple power supply and a lamp. The only EMC test events that could affect the performance or environment of the system are Harmonics, Flicker, Radiated and Conducted Emission, Surge Immunity, and Voltage Dips, ESD, Radiated RFI, EFT, Conducted RFI, and Power Frequency Magnetic Immunity should not affect the performance of this product; at the levels required in EN60601-1-2:2001.</p>			

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