



ADJUSTABLE EXAM TABLE

with Drawers and
Paper Dispenser

ADI996-01-DR3-BLU ADI996-01-DR3-BLK
ADI996-01-BLU ADI996-01-BLK

ADI996-01-DR3-BLK



**Adjustable Exam Table with Drawers
and Paper Dispenser, Black**

Upgrade your medical examination experience with our unique Viva Comfort Adjustable Exam Tables. These exceptional medical exam tables redefine convenience and functionality, boasting a unique combination of features that set them apart from all others. To create a safe and hygienic environment, each table is equipped with a built-in paper roll dispenser and paper cutter, allowing for quick liner changing between patients. Additionally, it is outfitted with a 2" thick foam padding to provide superior patient comfort and antimicrobial upholstery, helping guarantee a sterile environment for your patients.

Our table's versatility shines through with its adjustable backrest that can be adjusted to multiple positions, allowing healthcare professionals to customize the examination setup to their precise requirements. We've also included an easy-to-install three-drawer set that affix to the bottom of the table. This will provide ample storage for everyday medical supplies, helping streamline examinations. The table's foldable legs are adorned with rubber feet to protect your flooring from being scratched or marked up. Furthermore, they help provide a space-saving storage option when the table is not in use, making it an ideal addition to clinics and medical practices. These one-of-a-kind Adjustable Exam Tables promise to streamline your examination process while enhancing patient comfort, all in a sleek and sophisticated design that complements any medical setting.

FEATURES

Accommodating: Adjustable backrest that can be positioned at different heights.

Build: The stainless-steel frame maintains the longevity of the exam bed and is hassle-free to clean.

Compact: The unique foldable leg design helps provide a space-saving storage option when the table is not in use.

Built-In: The included paper roll dispenser and paper cutter allows for quick changing of the liner.

Storage: An included spacious, three-drawer set with each table gives ample space to hold medical supplies.

Sanitation: The antimicrobial mattress is designed to provide a safe and hygienic environment for patients by reducing the risk of cross-contamination.

SPECIFICATIONS

Adjustable Exam Table with Paper Dispenser	With Drawers		Without Drawers	
	ADI996-01-DR3-BLU	ADI996-01-DR3-BLK	ADI996-01-BLU	ADI996-01-BLK
Upholstery Color	Midnight Express	Black	Midnight Express	Black
Total Weight	84.43 lb.		60.18 lb.	
Overall Size	75.59" L x 27.56" W x 31.5" H		75.59" L x 27.56" W x 31.5" H	
Weight Capacity	440 lb.		440 lb.	
Material	Steel		Steel	
Product Certification	ISO9001, ISO13485, CE, FDA, CARB Compliant		ISO9001, ISO13485, CE, FDA, CARB Compliant	
Compatible Paper Roll Size	22.04" L x 5" H		22.04" L x 5" H	
Backrest Angles	0°, 32.5°, 45°, 54.5°, 62°, 67°		0°, 32.5°, 45°, 54.5°, 62°, 67°	
Height of Bed From Ground To Base	27.95"		27.95"	
Mattress	2-inch Padded Foam Antimicrobial PU Leather Upholstery		2-inch Padded Foam Antimicrobial PU Leather Upholstery	
Mattress Dimensions	Overall Mattress: 74.01" L x 27.55" W		Overall Mattress: 74.01" L x 27.55" W	
	Backrest Mattress: 24.01" L x 27.55" W		Backrest Mattress: 24.01" L x 27.55" W	
	Kneeboard Mattress: 50" L x 27.55" W		Kneeboard Mattress: 50" L x 27.55" W	
Drawer Dimensions	Overall Drawer: 19.02"H x 15.16" W x 17.72" D		Not Included	
	Fully Open Drawer: 14.56" D			
	Top Drawer: 4.72" H x 16.14" W			
	Middle Drawer: 5.71" H x 16.14" W			
	Bottom Drawer: 6.69" H x 16.14" W			


FEATURES



Adjustable Backrest



Three-Drawer Set



Paper Roll Dispenser



Paper Cutter

ADI996-01-BLU

Adjustable Exam Table with Paper Dispenser, Midnight Express



Foldable Legs

UPHOLSTERY COLOR OPTIONS



Black



Midnight Express