



SMALL UP AND DOWN STAIRCASE

ADIME906-01



Get ready to take your physical therapy sessions to new heights with our innovative Small Up & Down Staircase. We've created these stairs exclusively for physical therapy use as they are a game-changer in the world of rehabilitation. Don't let their compact size fool you; they pack a punch with a plethora of features that therapists will adore, and they can fit in small or large-spaced practices. Each staircase is crafted with a sturdy combination of steel fasteners and metal braces to endure the test of time and intensive therapy sessions. Located on each main pressure point are steel upright supports that will ensure the handrails are stable for your patients. We've designed the handrail support system to be one-piece to guarantee a comfortable grip for patients. The solid hardwood handrails, reinforced with a rigid bolt-through mounting and broad bottom support, offer safety and durability.

To ensure the staircase surface is non-slip, we've added two full-width pieces of safety tread to each step and four full-width pieces to the platform. These safety treads will help your patients stay stable on their feet while they ascend and descend. With its straightforward design, you can prevent turning around and have patients walk a straight path instead.

Our staircase is specifically engineered to have a 500 lb. bariatric load capacity, while the powder-coated steel uprights and brace-reinforced underside guarantee stability and support. Whether you have limited space or require seamless, unidirectional movement, invest in our exceptional Small Up & Down Staircase today.

FEATURES

Versatile: They can fit into small and large-spaced physical therapy practices.

Durable: Crafted with a sturdy combination of steel fasteners and metal braces, the staircase is built to last.

Comfortable: The one-piece handrail support system ensures a secure grip for patients during therapy sessions.

Safety: Solid hardwood handrails, reinforced with a rigid bolt-through mounting, offer protection, and long-lasting durability.

Non-Slip: We've added two full-width pieces of safety tread to each step and four full-width pieces to the platform for ultimate safety.

Layout: The design allows patients to ascend and descend without turning around on the staircase.

Enduring: Its hardwood plywood base with a topcoat finish ensures a long-lasting, pristine appearance and can be easily cleaned.

Capacity: With a 500 lb. bariatric load capacity, the staircase accommodates a wide range of patients with varying needs.

Construction: The powder-coated steel uprights and brace-reinforced underside guarantee stability and support during therapy sessions.

SPECIFICATIONS


Small Up And Down Staircase	
Material	Iron, Solid Wood, and Plywood
Overall Dimensions	47.6" H x 33" W x 54.7" L
Weight Capacity	500 lb.
Weight	100 lb.
Clearance	A minimum of an 8' ceiling is required.
Platform Dimensions	11.81" H x 29.92" W x 23.23" D
1 Step Dimensions	5.91" H x 29.92" W x 9.84" D
2 Step 1st Step Dimensions	3.94" H x 29.92" W x 9.84" D
2 Step 2nd Step Dimensions	5.91" H x 29.92" W x 9.84" D

FEATURES

Step Dimensions



SOLID HARDWOOD HANDRAILS, STEEL FASTENERS, & METAL BRACES OFFER OPTIMAL PROTECTION AND DURABILITY.



NON-SLIP SURFACE

TWO PIECES OF SAFETY TREAD ON EACH STEP AND FOUR PIECES ON THE PLATFORM FOR ULTIMATE SAFETY.

