

First Responder Kit

The **LSP First Responder Kit** was developed for the prepared EMT to keep in his/her personal vehicle should they be the first to arrive at an emergency scene.

Each kit comes complete with an aluminum cylinder, oxygen regulator, nasal cannula, and a disposable manual resuscitator ... all conveniently packed in a durable carry bag with a storage compartment and adjustable shoulder strap.

Although the kit is designed to provide a low-cost application for multiple placements, all the components are manufactured to high quality standards to meet the requirements of the EMS environment.

If you're the first on the scene, you can respond with confidence with the LSP First Responder Kit.



SPECIFICATIONS

Supply Pressure to Regulator: 500-2200 psi

Flow Capacity of Regulator: Constant flow outlet: 1, 2, 4, 6, 8, 10, 12, 15, 18, 20, 25 lpm Manual Resuscitator: Disposable; single-patient use; adult, with O₂ reservoir and cuffed mask;

Latex free

Weight: Approximately 8 lbs. complete

ORDERING INFORMATION

L903 First Responder Kit: 21030 oxygen regulator (CGA 870, 0-25 lpm,

hose barb outlet, gauge), L670-040 adult disposable manual resuscitator, 33239 adult cannula with 7' tubing, aluminum C cylinder (empty), carry bag



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA 314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S. Toll Free: 800/444-3940 FAX: 800/477-7701 Canada & Mexico Toll Free: 800/446-0552 FAX: 800/246-6201 International Telephone: 314/268-1683 FAX: 314/771-5183