



AM-100 Treatment Table

- | | | |
|--|---|---|
| <ul style="list-style-type: none"> • 1-Section Treatment Table • 27"W x 76"L | <ul style="list-style-type: none"> • Min. Height: 17" • Max Height: 36" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • 400 lb. Capacity • Non-Skid Adjustable Glides (No Casters) |
|--|---|---|



AM-150 Treatment Table

- | | | |
|--|---|--|
| <ul style="list-style-type: none"> • 1-Section Treatment Table • 27"W x 76"L | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" | <ul style="list-style-type: none"> • 400 lb. Capacity • Power-Assisted Caster System |
|--|---|--|

*Also available:
 AM-134 - 34"W with
 500 lb. capacity
 AM-140 - 40"W with
 500 lb. Capacity*