



### AM-SP300 Treatment Table

- 3-Section Treatment Table
- 27"W x 76"L
- Min. Height: 24"
- Max Height: 35.5"
- 1 1/2" Firm Density Foam
- Head: 15" x 27" with Contoured Face Opening
- Center: 22" x 27"
- Foot: 39" x 27"
- 400 lb. Capacity
- One Step Caster System



### AM-SP400 Traction Table

- 4-Section Traction Table
- 27"W x 76"L
- 99" Overall Length
- Min. Height: 24"
- Max Height: 35.5"
- 1 1/2" Firm Density Foam
- Head: 15" x 13" with Contoured Face Opening
- Thoracic: 17.5" x 27"
- Center: 10" x 27"
- Foot: 33" x 27"
- 400 lb. Capacity
- One Step Caster System

### Also Available:

AM-SP100  
AM-SP200

AM-SP350  
AM-SP450

AM-SP500  
AM-SP1044 4 Caster Base Option

See corresponding AM Series model for Top Configuration Information