

AM-SX SERIES

HI-LO TREATMENT TABLES



AM-SX5000 Treatment Table with Motorized Center

- | | | | |
|--|---|--|--|
| <ul style="list-style-type: none">• 3-Section Treatment Table• 3 Piece Head• 27"W x 72"L | <ul style="list-style-type: none">• Min. Height: 18"• Max Height: 37"• 1 1/2" Firm Density Foam | <ul style="list-style-type: none">• Head: 15" x 13" with Contoured Face Opening• Arms: 15" x 7"• Center: 19" x 27"• Foot: 38" x 27" | <ul style="list-style-type: none">• 400 lb. Capacity• 4 Total Lock Swivel Casters |
|--|---|--|--|

The AM-SX Series Treatment Tables offer a unique design which will work well for smaller areas in the clinic. This design has exceptional stability, full height adjustment and section mobility. The flexibility of these tables allows the patient to be placed in proper position to achieve maximum results with minimum effort.