

MATERIAL SAFETY DATA SHEET

PRODUCT NAME: E-Z SCRUB™ Scrub 7471749 4% CHG
Catalog Nos. 371073, 377477 and 377479
Date: 04/19/10 **Revision #:** 4



PRODUCT AND COMPANY INFORMATION – SECTION 1

Company Identification

Name: Becton Dickinson Medical Systems
Address: 9450 South State Street
Sandy, UT 84070

FOR CHEMICAL EMERGENCY, SPILL, LEAK, FIRE, EXPOSURE, OR ACCIDENT:

In the continental U.S., call: 801.565.2300 (24 hour)
For additional information, call: 888.237.2762
(M-F, 8:00 AM-4:30 PM EST)

COMPOSITION INFORMATION – SECTION 2

Product Name: BD E-Z Scrub with Chlorhexidine Gluconate as Surgical Hand Scrub
Chemical Formula: Proprietary Mixture

<u>Chemical Name</u>	<u>Percent Wt.</u>	<u>CAS Number</u>	<u>RTECS Number</u>
Chlorhexidine Gluconate	4.0	18472-51-0	DU1950000

This is a water-based solution of surfactants and Chlorhexidine Gluconate.

OSHA REGULATORY STATUS

	This material is classified as not hazardous under OSHA regulations.
X	While this material is not classified as hazardous under OSHA regulations, this MSDS contains valuable information critical to the safe handling and proper use of the product. This MSDS should be retained and available for employees and other users of the product.
	This material is classified as hazardous under OSHA regulations.

HAZARDS IDENTIFICATION – SECTION 3

EMERGENCY OVERVIEW

- Clear, pink/red liquid with perfumed odor.
- Presents little or no hazard if spilled.
- Presents little or no fire hazard.
- For external use only: keep out of eyes, ears, and mouth.

POTENTIAL HEALTH EFFECTS

Eyes: May cause permanent damage.
Ears: May cause permanent hearing loss.

MATERIAL SAFETY DATA SHEET

PRODUCT NAME: E-Z SCRUB™ Scrub 7471749 4% CHG
Catalog Nos. 371073, 377477 and 377479
Date: 04/19/10 **Revision #:** 4



POTENTIAL HEALTH EFFECTS

- Skin:** May cause irritation in hypersensitive individuals. Avoid contact with meninges.
- Inhalation:** Overexposure may cause headache, fatigue, dizziness, loss of coordination, and unconsciousness. Vapor has anesthetic properties.
- Ingestion:** Harmful if swallowed.

FIRST AID MEASURES – SECTION 4

Seek medical attention for all cases of overexposure.

FIRST AID MEASURES

- Eyes:** Flush thoroughly with water for at least 15 minutes. Seek medical attention.
- Ears:** Flush with water, repeating several times. Seek medical attention.
- Skin:** Flush with water. Seek medical attention.
- Ingestion:** Wash out mouth with water. Treat with gastric lavage using milk, egg white, or mild soap. Seek medical attention.
- Inhalation:** Remove patient from exposure. Keep warm and at rest. Seek medical attention if ill effects occur.

FIRE FIGHTING MEASURES – SECTION 5

	<u>NFPA Rating</u>	
Health	Fire	Reactivity
1	0	0

Flammability Properties

- Flash Point / Method:** N/A
- Flammability Limits:** LEL: N/A UEL: N/A
(in air % by volume)
- Autoignition Temperature:** U/K
- Hazardous Combustion Products:** Carbon dioxide, carbon monoxide, hydrogen chloride gas, nitrogen oxides, ammonia.
- Extinguishing Media:** Water, foam, dry chemical, carbon dioxide.
- Fire Fighting Instructions:** As in any fire, wear full protective clothing and NIOSH-approved self-contained breathing apparatus with full facepiece operated in the pressure demand or other positive pressure mode.
- Unusual Fire and Explosion Hazards:** Emits toxic fumes under fire conditions.

MATERIAL SAFETY DATA SHEET

PRODUCT NAME: E-Z SCRUB™ Scrub 7471749 4% CHG
Catalog Nos. 371073, 377477 and 377479
Date: 04/19/10 **Revision #:** 4



ACCIDENTAL RELEASE MEASURES - SECTION 6

Absorb spilled material with suitable absorbent. Transfer to container for disposal in accordance with regulatory requirements. Mop area clean.

HANDLING & STORAGE – SECTION 7

Handling: External use only. Avoid contact with eyes, ears, and mouth.

Storage: Store at room temperature. Protect from light. Avoid excessive heat (104°F).

EXPOSURE CONTROLS / PERSONAL PROTECTION – SECTION 8

Engineering Controls: General ventilation for conditions of prescribed use.

Personal Protective Equipment

Eye/Face Protection: Not required for normal use. Wear safety goggles or a full face shield for spill cleanup or if splashing is possible.

Skin Protection: Not required for normal use. Wear impervious gloves and apron for spill cleanup.

Respiratory Protection: Not required for normal use. Use NIOSH approved respirator for organic vapors, if necessary.

Exposure Guidelines

OSHA PEL: N/A

ACGIH TLV: N/A

PHYSICAL & CHEMICAL PROPERTIES – SECTION 9

Appearance, Physical State: Clear, pink / red liquid

Odor: Perfumed

Molecular Weight: N/A

pH: 6-7.5

Boiling Point: Similar to water

Melting Point: Similar to water

Vapor Pressure (mm Hg): Similar to water

Vapor Density (Air=1): Similar to water

Solubility in Water: Soluble

Specific Gravity: 1.06

Viscosity @ 25°C: 50-500 CPS

MATERIAL SAFETY DATA SHEET

PRODUCT NAME: E-Z SCRUB™ Scrub 7471749 4% CHG
Catalog Nos. 371073, 377477 and 377479
Date: 04/19/10 **Revision #:** 4



STABILITY & REACTIVITY – SECTION 10

Stability: Stable under normal conditions
Conditions to Avoid: High heat
Materials to Avoid: Strong oxidizing agents, anionic compounds
Hazardous Decomposition Products: Carbon monoxide, carbon dioxide, nitrogen oxides, ammonia, hydrogen chloride gas
Hazardous Polymerization: Will not occur

TOXICOLOGICAL INFORMATION – SECTION 11

Chlorhexidine gluconate: oral LD50 = 1260 mg/kg (mouse) - 2 g/kg (rat)

No significant human data is available for this material.

ECOLOGICAL INFORMATION – SECTION 12

Harmful to aquatic organisms; may cause long-term adverse effects in the aquatic environment. Low toxicity to sewage microorganisms.

DISPOSAL CONSIDERATION – SECTION 13

Dispose of in accordance with Federal, state and local regulations. Discarded product is not a hazardous waste under RCRA, 40 CFR 261. Empty container retains product residue. Observe all hazard precautions.

TRANSPORT INFORMATION - SECTION 14

DOT CLASSIFICATION: Not regulated in non-bulk packaging

IATA CLASSIFICATION: Not Regulated

IMDG CLASSIFICATION: Not Regulated

REGULATORY INFORMATION – SECTION 15

US REGULATIONS

OSHA Not Regulated
TSCA All ingredients are on the TSCA Chemical Substances Inventory (8(b)).
CERCLA Not Regulated
SARA Title III Not Regulated

MATERIAL SAFETY DATA SHEET

PRODUCT NAME: E-Z SCRUB™ Scrub 7471749 4% CHG
Catalog Nos. 371073, 377477 and 377479
Date: 04/19/10 **Revision #:** 4



STATE REGULATIONS

MA Substance List	Not Regulated
NJ Right-to-Know Hazardous Substance List	Not Regulated
PA Hazardous Substance List	Not Regulated
Canadian WHMIS	Not Regulated

OTHER INFORMATION - SECTION 16

Definitions:

UIK = Unknown **NIA** = Not Applicable

To the best of our knowledge, the information contained herein is accurate. However, neither Becton Dickinson nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist. The possibility exists that the EU will not recognize this MSDS due to the fact that several components of the MSDS are reflective of ANSI Z-400.1-1998. Although ILO (International Labor Organization) has adopted ANSI Z-400.1-1998, ultimate disposition lies with the competent authority.

MATERIAL SAFETY DATA SHEET

PRODUCT NAME: E-Z SCRUB™ Scrub 7471749 4% CHG
Catalog Nos. 371073, 377477 and 377479
Date: 04/19/10 **Revision #:** 4



VERSION HISTORY

Version	Description	ECO Number	Version Prepared By
4	Review and Update: From the Header: Remove Obsolete Catalog Number 377496 add Catalog Number 371073. Update the Date to 4/19/10 and the Revision to 4. Also, corrected QDMS Attributes from WWB104 to WWB108.	ECO123316	Rachel LeBlanc Janet Boggess
3	Review and update document. Show review date in header (10/29/08) and footer (10-08). After review, update section 14 to remove "Combustible liquid" So this section becomes: "Not regulated in non-bulk packaging". CCB# 1599 (10/30/2008).	ECO92782	Janet Boggess Rachel LeBlanc
2	Date of MSDS 14536 and D14554 updated to current date 9/13/05 to comply with regulation in Canada. And change Version number.	ECO35390	George Nolan