

MEDICAL LIGHTS FROM A NAME YOU TRUST



SPECIALTY LIGHTING

Aaron specialty lighting consists of a line of specialty ophthalmic and penlights including the 6802A Fluoro-Dot™ and the Ultraviolet Woods Light (UV59), which is excellent for detecting certain skin disorders and for use in rape detection procedures.



FLEXIBLE & LIGHTS

Aaron manufactures two groups of flexible lighting products.

The first is available in three different flexible shaft lengths 5", 10", & 15" (12.7 cm, 25.4 cm, & 38.1 cm) and is supplied sterile for any surgical procedure requiring additional lighting.

The second group of patented flexible lights are lighted stylets. When using the lighted stylet, the user is assured of the correct placement of the endotracheal tube when a bright glow in the midline of the laryngeal prominence is seen.



7100 30th Avenue North • St. Petersburg, FL 33710-2902
 U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640
 Int'l Phone +727-384-2323 • Fax +727-347-9144
www.aaronmed.com • sales@aaronmed.com

MEDICAL LIGHTS EXPLANATIONS & STYLE NUMBERS

LIGHTED INTUBATION STYLETS & SURGICAL LIGHTS

U.S. Patent #4,517,499 & #4,347,553
Canada Patent #1192171
Far East Patent #36224

These special flexible lights were originally designed for use in general surgery and orthopaedics. They are protected by two United States, one Canada, and one Far East Patent. The non-sterile versions of these lights are for routine non-office examinations or where sterility is not required. They also have the outstanding feature of replaceable battery and bulb.

Surch-Lite™ 5" (non-sterile) No. NS05
Five-inch flexible shaft light with replaceable batteries and bulb. Individually packaged.

Surch-Lite™ 10" (non-sterile) No. NS10
Ten-inch flexible shaft light with replaceable batteries and bulb. Individually packaged.

Surch-Lite™ 15" (non-sterile) No. NS15
Fifteen-inch flexible shaft light with replaceable batteries and bulb. Individually packaged.

Surch-Lite™ 5" (sterile) No. ST05
Five-inch flexible shaft individually packaged sterile in boxes of ten.

Surch-Lite™ 10" (sterile) No. ST10
Ten-inch flexible shaft individually packaged sterile in boxes of three.

Surch-Lite™ 15" (sterile) No. ST15
Fifteen-inch flexible shaft individually packaged sterile in boxes of three.

Aaron's Surch-Lite™ Orotracheal Lighted Intubation Stylets assist in the visualization required for the proper placement of endo-tracheal tubes. By using the lighted stylet, the user is assured of the correct placement of the endotracheal tube when a bright glow in the midline of the laryngeal prominence is seen.

Flexible shafts, a protective coating over the bulb housing, and a one-hand activating on/off switch provide easy operation, safety and convenience.

The Aaron lighted intubation stylet is designed for emergency and trauma situations both in the field and E.R., as well as intubation procedures in the operating room. It also serves as an effective alternative to the standard laryngoscope.

Surch-Lite™ Orotracheal No. SLOT
Lighted Intubation Stylet (shaft length 40cm).
No. NLOT (non-sterile version of SLOT).

Surch-Lite™ Pediatric Orotracheal No. PDOT
Lighted Intubation Stylet (shaft length 25cm).
No. PDON (non-sterile version of PDOT).

SPECIALTY LIGHTS

Tech-Lite™ No. 6800 (6800A = Replaceable Batteries)
Specially designed lens-end bulb gives an intense concentric beam which allows more effective illumination. Larger than the standard penlight, its AA-sized Alkaline batteries give a shelf life of two years and a continuous burn time of four hours.

White Slit™ No. 6801 (6801A = Replaceable Batteries)
Designed with a high-intensity white slit beam for cataract screening to determine the depth of the anterior chamber and to check for corneal edema. The White Slit™ has a two-year shelf life and a four-hour continuous burn time.

Fluoro Dot™ No. 6802 (6802A = Replaceable Batteries)
This AA-sized light combines the outstanding features of a two-year shelf life and four-hour burn time with a special Wratten No. 47A cobalt blue filter. This light is designed for use with fluorescein dye and is excellent for diagnosing abrasions and fitting contact lenses.

Fluoro-Slit™ No. 6803 (6803A = Replaceable Batteries)
Designed with a special bulb and Wratten No. 47A cobalt-blue filter to project a concentrated cobalt-slit image for diagnostic use with fluorescein dye, the Fluoro-Slit™ has a two-year shelf life and a four-hour continuous burn time.

1mm Light No. 6804 (6804A = Replaceable Batteries)
This AA-sized light combines the outstanding features of a two-year shelf life and four-hour burn time. This light is designed for use as an aid for muscle tests with a 1mm beam of light.

Ultraviolet Woods Light™ No. UV59
This hand-held ultraviolet florescent light is excellent for detecting certain skin disorders and its compact size makes it extremely portable. The light includes four AA-sized alkaline batteries, plus the added benefit of a standard krypton flashlight on one end.

PENLIGHTS

Disposable Penlight 6/pack No. 6666
Penlight w/Pupil Gauge - 6/pack No. 66RN
Replaceable Battery Penlight - 6/pack No. 666R
Cobalt Penlight - 3/pack No. 6303

Deluxe Penlight No. A238
Designed with black anodized body and silver accents - 1/box
Disposable Penlight - 1/pack No. 1125
Reusable Penlight - 1/pack No. 1135
Cobalt Filter Only - 1/pack No. 7156

ST05Flexible Sterile Light 5" (NS05 = Non-Sterile)
ST10Flexible Sterile Light 10" (NS10 = Non-Sterile)
ST15Flexible Sterile Light 15" (NS15 = Non-Sterile)
SLOTOral Lighted Stylet Sterile (NLOT = Non-Sterile)
PDOTOral Pediatric Lighted Stylet Sterile (PDON = Non-Sterile)
6800Tech-Lite™ (6800A = Replaceable Batteries)
6801White Slit™ (6801A = Replaceable Batteries)
6802Fluoro Dot™ (6802A = Replaceable Batteries)
6803Fluoro-Slit™ (6803A = Replaceable Batteries)
68041mm Light (6804A = Replaceable Batteries)
UV59Ultraviolet Woods Light

6666Disposable Penlight 6/pack
1125Disposable Penlight 1/pack
66RNDisposable Penlight w/Pupil Gauge 6/pack
666RReplaceable Battery Penlight 6/pack
1135Replaceable Battery Penlight 1/pack
6303Cobalt Blue Penlight 3/pack
7156Cobalt Blue Filter Only 1/pack
A238Deluxe Penlight 1/pack

Distributed By: