

## A single, convenient source for so many solutions.

Brewer understands that you need more than quality power tables to equip your exam and procedure rooms for success. That's why we offer a broad selection of ruggedly constructed clinical support solutions to complete your clinical care space. Each product is manufactured to Brewer's rigorous quality standards for design, materials and craftsmanship to provide the reliable, trouble-free service you need. Plus, all our clinical support products are backed by a standard 5-Year Warranty for added peace of mind.



## Hampers & Carts

### Durable designs you can depend on.

#### 18" Round Linen Hamper (45.7 cm)

Durable steel construction with white powder-coat finish. Dimensions: 35.5" H x 18.25" outside diameter (90.2 cm x 46.3 cm).

**Model 33330**

#### Accessories

- 31331 - 18" Rope style  
10 oz. cotton bag
- 34000 - Bag holder



#### 18" Tilt-Top Round Linen Hamper (45.7 cm)

Durable steel construction with white powder-coat finish. Step-on pedal for hands-free use. Dimensions: 35.75" H x 18" outside diameter (90.6 cm x 45.7 cm).

**Model 33334**

#### Accessories

- 31331 - 18" Rope style  
10 oz. cotton bag
- 34000 - Bag holder



#### Multi Purpose Cart

Strong steel frame with durable white powder coat finish. Total weight capacity: 50 lbs (22.7 kg). Dimensions with handle:

37.5" H x 19.5" W x 29" D (95.3 cm x 49.5 cm x 73.6 cm).

Baskets sold separately.

#### Models

- 63400 - Cart frame only
- 63021 - 12" White powder coat wire basket
- 63020-1 - 6" White powder coat wire basket (1 per carton)
- 63020-2 - 6" white powder coat wire basket (2 per carton)



#### Adjustable Folding Linen Hamper

Holds 18", 22" or 25" diameter bags (45.7 cm, 55.9 cm or 63.5 cm).

Features durable white powder-coat finish. Folds for compact storage.

Dimensions: 37" H x 20" x 20" D (93.9 cm x 50.8 cm x 50.8 cm).

**Model 33395**

#### Accessories

- 31331 - 18" Rope style  
10 oz. cotton bag
- 31341 - 25" Rope style  
10 oz. cotton bag



#### Plastic Utility Cart

Includes three ABS plastic shelves, 16" W x 21.75" L (40.6 cm x 55.3 cm). Total weight capacity: 200 lbs (90.7 kg). Dimensions with handle: 28" L x 17.5" W x 28" H (71.1 cm x 44.5 cm x 71.1 cm).

**Model 63530**



#### Stainless Steel All Purpose Cart

Features non-corrosive stainless steel construction & three 15.5" W x 24" L (39.4 cm x 61 cm) shelves. Total weight capacity: 300 lbs (136 kg). Dimensions with handle: 28" L x 16.5" W x 32" H (71.1 cm x 41.9 cm x 81.3 cm).

**Model 63500**



#### Tilt-Top Square Hamper

Open design for easy unloading. Step-on pedal for hands-free use. White powder-coat finish. Built-in bag holder supports heavy-duty bags. Dimensions: 38" H x 18.5" W x 19.5" D (96.5 cm x 47 cm x 49.5 cm).

**Model 11410**

#### Accessories

- 31331 - 18" Rope style  
10 oz. cotton bag



#### Mobile Equipment Cart

Powder coated steel cabinet with integral handles, locking casters, & rubber bumpers to prevent wall damage. Has three small drawers & one large drawer.

Dimensions: 34.4" H x 24.2" W x 21.6" D (87.4 cm x 61.5 cm x 54.9 cm).

**Model 63100**

#### Accessories

- 100407 - Drawer dividers

