

A single, convenient source for so many solutions.

Brewer understands that you need more than quality power tables to equip your exam and procedure rooms for success. That's why we offer a broad selection of ruggedly constructed clinical support solutions to complete your clinical care space. Each product is manufactured to Brewer's rigorous quality standards for design, materials and craftsmanship to provide the reliable, trouble-free service you need. Plus, all our clinical support products are backed by a standard 5-Year Warranty for added peace of mind.



Durable designs you can depend on.

18" Round Linen Hamper (45.7 cm)

Durable steel construction with white powder-coat finish. Dimensions: 35.5" H x 18.25" outside diameter (90.2 cm x 46.3 cm).

Model 33330

Accessories

- 31331 - 18" Rope style
- 10 oz. cotton bag
- 34000 - Bag holder



Multi Purpose Cart

Strong steel frame with durable white powder coat finish. Total weight capacity: 50 lbs (22.7 kg). Dimensions with handle:

37.5" H x 19.5" W x 29" D (95.3 cm x 49.5 cm x 73.6 cm). Baskets sold separately.

Models

- 63400 - Cart frame only
- 63021 - 12" White powder coat wire basket
- 63020-1 - 6" White powder coat wire basket (1 per carton)
- 63020-2 - 6" white powder coat wire basket (2 per carton)



Plastic Utility Cart

Includes three ABS plastic shelves, 16" W x 21.75" L (40.6 cm x 55.3 cm). Total weight capacity: 200 lbs (90.7 kg). Dimensions with handle: 28" L x 17.5" W x 28" H (71.1 cm x 44.5 cm x 71.1 cm).

Model 63530



18" Tilt-Top Round Linen Hamper (45.7 cm)

Durable steel construction with white powder-coat finish. Step-on pedal for hands-free use. Dimensions: 35.75" H x 18" outside diameter (90.6 cm x 45.7 cm).

Model 33334

Accessories

- 31331 - 18" Rope style
- 10 oz. cotton bag
- 34000 - Bag holder



Adjustable Folding Linen Hamper

Holds 18", 22" or 25" diameter bags (45.7 cm, 55.9 cm or 63.5 cm). Features durable white powder-coat finish. Folds for compact storage. Dimensions: 37" H x 20" x 20" D (93.9 cm x 50.8 cm x 50.8 cm).

Model 33395

Accessories

- 31331 - 18" Rope style
- 10 oz. cotton bag
- 31341 - 25" Rope style
- 10 oz. cotton bag



Stainless Steel All Purpose Cart

Features non-corrosive stainless steel construction & three 15.5" W x 24" L (39.4 cm x 61 cm) shelves. Total weight capacity: 300 lbs (136 kg). Dimensions with handle: 28" L x 16.5" W x 32" H (71.1 cm x 41.9 cm x 81.3 cm).

Model 63500



Tilt-Top Square Hamper

Open design for easy unloading. Step-on pedal for hands-free use. White powder-coat finish. Built-in bag holder supports heavy-duty bags. Dimensions: 38" H x 18.5" W x 19.5" D (96.5 cm x 47 cm x 49.5 cm).

Model 11410

Accessories

- 31331 - 18" Rope style
- 10 oz. cotton bag



Mobile Equipment Cart

Powder coated steel cabinet with integral handles, locking casters, & rubber bumpers to prevent wall damage. Has three small drawers & one large drawer. Dimensions: 34.4" H x 24.2" W x 21.6" D (87.4 cm x 61.5 cm x 54.9 cm).

Model 63100

Accessories

- 100407 - Drawer dividers

