

# Brewer

## *FLEX* ACCESS EXAM TABLE



Provide a better patient experience and drive your bottom-line with a sleek, flexible exam table.



- Unmatched ergonomic patient access
- Industry leading 700lb weight capacity
- Easy, safe and more independent patient movement
- Multiple options to suit your practice



Chair



Table



Patient Transfer Grab Bars



SafeGlide™

### Chair-To-Table Design

- Efficiently and smoothly moves patient from consultation to examination
- 18" low height for easy and safe patient transfers
- Return to chair feature easily takes table back to chair position with one press of a button (5800/5801 models)

### Advanced Safety

- Exclusive SafeGlide™ leg extension increases safety and satisfaction
- Optional patient transfer grab bars provide patients with a stable and secure grab point to use when getting on and off table
- Transfer grab bars conveniently swing out for simplified and more independent lateral wheelchair transfers
- Optional ergonomic stirrups



### Comfortable, Ergonomic Upholstery

- Seamless, medicinal stain resistant
- Available in 12 standard colors
- CAL-133 fire resistant upholstery available in 4 colors
- An array of additional colors available (call for pricing)



### Reliable, Durable Design

- Industry-leading standard 3-year warranty
- Easy to operate and service
- UL 60601-1 Approved
- Meet Brewer's ISO13485:2003 quality process

Brewer

Furnishing solutions for safe practices.

# Brewer

## *FLEX* ACCESS EXAM TABLE

### Models

- #5700 - FLEX Access Exam Table with pneumatic back
- #5701 - FLEX Access Exam Table with pneumatic back, pelvic tilt, front drawer heater and outlet
- #5800 - FLEX Access Exam Table with power back
- #5801 - FLEX Access Exam Table with power back, pelvic tilt, front drawer heater, and outlet

### Accessories and Kits

#### Kits:

- #2102894: Ergonomic Stirrup for 5700 and 5800 (factory installed)
- #2102898: Ergonomic Stirrup for 5700 and 5800 (field installed)
- #2102893: Ergonomic Stirrup for 5701 and 5801 (factory installed)
- #2102897: Ergonomic Stirrup for 5701 and 5801 (field installed)
- #2102891: Safety Grab Bars (factory installed)
- #2102895: Safety Grab Bars (field installed)
- #2102892: Pass-through laptop work surface (factory installed)
- #2102896: Pass-through laptop work surface (field installed)

#### Accessories:

- #98464 stainless steel treatment pan
- #99502 one universal socket clamp
- #99504\* deluxe articulating knee crutches
- #100253 urology drain pan system
- #100407 front drawer dividers (can order up to 7)
- #100540 two support adapters for knee crutches (plant installed)
- #100550 two support adapters for knee crutches (field installed)
- #102031 Brewer Welch Allyn® IV bracket (plant installed)
- #102032 Brewer Welch Allyn® IV bracket (field installed)

\* Needs two #99502 and one of either #100540 or #100550

### Certifications

- UL 60601-1



### Advanced Safety

- Exclusive SafeGlide™ leg extension increases patient and caregiver safety and satisfaction
- Return to chair feature - (5800 and 5801 models)
- Anti-tamper on/off switch prevents accidental movement
- Sensors minimize injuries during table movement
- ADA compliant at 18" seat height
- Easy on / off for patients

### Patient Comfort

- Easy lateral wheelchair transfers
- Multi-density support upholstery

### Large Front Drawer Storage

- Front drawer: 19.8" x 19.2" x 3.8" (50.3 cm x 48.8 cm x 9.7 cm)
- Front drawer storage can accommodate up to 6 - 18" x 3" paper rolls

### Table Specifications

- 700 lb. (317 kg) weight capacity
- Height: 18" - 37" (45.7 cm - 94 cm)
- SafeGlide™ Leg extension length: 21" (53.3 cm)
- Standard foot control
- Standard 3-year warranty
- Standard easy to clean, seamless upholstery
- 12 standard upholstery color options
- CAL-133 fire resistant upholstery available in 4 colors
- Custom colors also available
- Length: 57" extends to 78" (144.8 cm - 198.1 cm)
- Upholstery width: 28" (71.1 cm)
- 115VAC, 60 Hz, 6.0 amps max
- Safe, low voltage actuators
- Safe, single power cord design, detachable

The Power to Advance

The Brewer Company is an ISO 13485:2003 certified company.

Call your local Brewer sales representative to learn about our complete product line.

Brewer

The Brewer Company LLC  
N88W13901 Main Street, Menomonee Falls, WI 53051  
1.888.Brewer.1 | www.brewercompany.com