

Model DS1150

# Mechanical Height Rod

Precision engineered and constructed of heavy gauge aluminum, the Doran Model DS1150 Mechanical Height Rod will deliver accurate height readings for years and years. Flip the measurement arm up and slide the height measurement arm into position. The patient's height is highlighted with an easy to read, bright red line.

Marked in millimeters and inches, the Model DS1150 mechanical height rod can be used to measure adult or pediatric patients.



## MEASURING RANGE

24 - 82" (60 - 210 cm)

## GRADUATIONS

1/16" (1 mm)

## CONSTRUCTION

Aluminum and ABS Plastic

## HEIGHT ROD DIMENSIONS

59" L x 3" W x 1" H

## MEASUREMENT ARM DIMENSIONS

11" L x 8" W x 3.5" H

## WARRANTY

**2** Two Years



The red line on the height rod ensures easy and accurate height measurements.



**Doran Scales, Inc.** 1-800-264-4106  
1315 Paramount Parkway Fax: 630-879-0073  
Batavia, IL 60510 [www.doranmedical.com](http://www.doranmedical.com)