



Black Latex Examination Gloves

Maximum Protection Against Exposure Hazards

Dynarex Black Latex Examination Gloves are ideal for law enforcement, correctional facilities, tattoo shops, and high-end beauty salons. They are also a good choice for light-duty maintenance when protection is required.

- **Maximum protection, puncture resistant, durable, tactile sensitivity, and comfortable fit**
- **Powder free to minimize skin irritation (≤ 2 mg residual powder per glove)**
- **Textured for secure grip in wet or dry use**
- **Polymer coated and ambidextrous**
- **3-5mil. range min. thickness at finger & palm**

Technical Specifications: Dynarex Black Latex Examination Gloves are medical grade gloves that meet or exceed ASTM D3578 specifications. Protein content is 50µg/dm² or less.

CAUTION: This product contains natural rubber latex which may cause allergic reactions. Do not expose this product to any person known or suspected to be latex sensitive. Do not use in latex-free areas or on patients with spina bifida.



Reorder No.	Black Latex – Powder Free	Quantity
2321	Black Latex Exam Gloves - S	10/100/cs
2322	Black Latex Exam Gloves - M	10/100/cs
2323	Black Latex Exam Gloves - L	10/100/cs
2324	Black Latex Exam Gloves - XL	10/100/cs

Visit our website to see our other types of gloves.



For more information contact your Dynarex sales representative 888-396-2739