



WE CARE LIKE FAMILY™



## Nasopharyngeal Airways

### Used To Deliver Oxygen In Higher Concentrations

Recommended for use in the semi-conscious or unconscious patient with an intact gag reflex. Not easily dislodged during transport and better tolerated than an Oropharyngeal Airway. Designed to be inserted into the nasal passageway to secure an open airway. When a patient becomes unconscious, the muscles in the jaw commonly relax and can allow the tongue to slide back and obstruct the airway.

- Sterile
- Flared end prevents the device from becoming lost inside the patient's nose
- Soft – no damage to nasal passage
- Flexible – resists collapse or kinking
- Not made with natural rubber latex



Reorder No.	FR Size #	Spec ID Millimeter	Quantity
4591	20Fr.	5.0 ± 0.15	10/box
4592	22Fr.	5.5 ± 0.15	10/box
4593	24Fr.	6 ± 0.15	10/box
4594	26Fr.	6.5 ± 0.15	10/box
4595	28Fr.	7 ± 0.15	10/box
4596	30Fr.	7.5 ± 0.15	10/box
4597	32Fr.	8 ± 0.15	10/box
4598	34Fr.	8.5 ± 0.15	10/box
4599	36Fr.	9 ± 0.15	10/box

See our website for our complete line of Laryngeal Mask Airways (Reorder No. 5032-5044)



Find us on:



**dynarex**® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)