

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

V-Series, Six Drawer Crash Cart with Breakaway Lock **#V30-6B**

This high capacity crash cart is designed for healthcare facilities under budget pressure requiring a durable, fully equipped emergency cart offering Harloff's quality construction at a value price.



Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand. White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Specifications

Construction:

Steel construction

Features:

- · Break-away lock bar
- 100 break-away seals included
- 30" of vertical drawer space
- Aluminum side rail / push handle
- · Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel two with
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- · Durable, hard-baked powder coat finish

Internal Drawer Dimensions:

• 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D)
- Cart weight: 118 lbs (53.5 kg)
- Approx shipping weight: 141 pounds (63.9 kg)

Drawer Configuration:

- (4) 3" (7.6 cm) (1) 6" (15.2 cm)
- (1) 12" (30.5 cm)

Integrated Pull Handle on all drawers