



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Value Line Medication Cart with Key Lock # VLT216BOX

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Color Shown is Light Gray

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system with 4 keys
- Internal locking narcotics box in bottom 9" drawer
- Capacity of up to 216 unit dose boxes or 60 sections for strip packaging
- Pull out Side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel – one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

Drawer Configuration:

- 1 - ea 3" multipurpose storage drawer
- 6 - ea 3" drawers with internal unit dose dividers
- 1 - ea 9" drawer with 1 divider and locking narcotics box

Internal Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75" H x 29.5" W x 22.5" D
(121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart weight: 170 pounds (77.3 kg)
- Approx. shipping weight: 210 pounds (95.5 kg)

Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>