



A. Model 1427-57-121



Common Features

MAT PLATFORMS — SERIES 1427, 1429:



- Extremely stable lift system elevates from 20" wheelchair height to 30" high.
- "Hands-Free" Foot control.
- Seamless nylon-reinforced vinyl is wear and tear resistant.
- 2" High-Density urethane foam vinyl top.
- Steel frame and 4 corner legs have durable black powder coat finish.
- Legs have levelers for uneven floors.
- 4-Leg "Open" design is patient sling-lift accessible on all four sides.

COLORS: (SEE PAGE 42)

Upholstery – 113 Blue, 119 Gray, 121 Green, 122 Black

**You must specify upholstery color with order.*

A. Series 1427 Bariatric Mat Platforms



SPECIAL FEATURES:

- Electrically activated hydraulic cylinders are built into each (4) corner steel leg.
- Motorized lift system can securely lift 750 lbs. from 20" wheelchair height to 30" high..
- "Hands-Free" foot control operation.

MODEL		WIDTH	LENGTH	HEIGHT
1427-46*	<i>New Product!</i>	4'	6'	20"-30"
1427-47*		4'	7'	20"-30"
1427-57*		5'	7'	20"-30"
1427-68*		6'	8'	20"-30"

B. Series 1429 Bariatric 1100 lb. Capacity Heavy-Duty Mat Platform



SPECIAL FEATURES:

- Electrically activated hydraulic cylinders are built into each (4) corner steel leg.
- Motorized lift system can securely lift 1100 lbs. from 20" wheelchair height to 30" high.
- Static load tested to 4400 lbs.
- Seamless nylon reinforced vinyl is tear resistant.
- Plywood reinforced subtop for added strength.

MODEL	WIDTH	LENGTH	HEIGHT
1429-57	5'	7'	20"-30"
1429-68	6'	8'	20"-30"

B. Model 1429-68-113

