



A. Model 1566



B. Model 1567

## Common Features

EXERCISE STAIRCASES —  
MODELS 1560-65-036, 1566, 1567:



- Handrails and uprights are satin finished premium hardwood and shipped R.T.A.
- Handrails are ergonomically designed for comfortable grip and ease of use.
- Stairs are finished in Natural Oak laminate with "anti-slip" safety treads on all steps and platform.
- BriteStep™ has Yellow Hi-Visibility, Hi-Impact 3mm edging for enhanced safety & durability.

### A. Model 1566 "Mini" Staircase



#### SPECIAL FEATURES:

- Space Saver can fit in 8 ft. high ceiling.
- (2) 6" high steps and (3) 4" high steps. Steps are 30" wide.
- Platform is 30"W x 24"D x 12"H.

MODEL	LENGTH	WIDTH	HEIGHT
1566	52"	30"	12"

### C. Model 1560-65-036 Convertible Staircase



#### SPECIAL FEATURES:

- Sets up and converts to straight or corner.
- (4) 6" H x 36" W steps on one end. (6) 4" H x 36" W steps on the other end.
- 36" x 36" x 24" high platform.
- Requires 9 ft. minimum ceiling.

MODEL	LENGTH	WIDTH	HEIGHT
1560-65-036	109"	36"	24"
Corner	64" x 82"	36"	24"

### B. Model 1567 Straight Staircase



#### SPECIAL FEATURES:

- Handrails enclose platform.
- (4) 6" high x 30" wide steps.
- Platform is 30" wide x 24" deep x 24" high.
- Requires 9 ft. minimum ceiling.

MODEL	LENGTH	WIDTH	HEIGHT
1567	52"	30"	24"



C. Model 1560-65-036  
Straight setup



Corner setup



A. Model 1569



B. Model 1582

*New Product!*

## Common Features

EXERCISE STAIRCASES —  
MODELS 1569, 1580, 1582:



- Handrails and uprights are satin finished premium hardwood and shipped R.T.A.
- Handrails are ergonomically designed for comfortable grip and ease of use and have 12" A.D.A. overhang.
- More compact footprint takes up less space and fits in conventional 8' high ceiling.
- Stairs are finished in Natural Oak laminate with "anti-slip" safety treads on all steps and platform.
- BriteStep™ has Yellow Hi-Visibility, Hi-Impact 3mm edging for enhanced safety & durability.

### A. Model 1569 ADA "Mini" Staircase



**SPECIAL FEATURES:**

- (2) 6" high steps and (3) 4" high steps. Steps are 30" wide.
- Platform is 30"W x 24"D x 12"H.
- A.D.A. 11"D steps x 30"W have anti-slip safety treads.

MODEL	LENGTH	WIDTH	HEIGHT
1569	55"	30"	12"

### B. Model 1582 ADA Straight Staircase



**SPECIAL FEATURES:**

- A.D.A. 11" deep steps x 36" wide have "anti-slip" safety treads.
- (3) 6" high steps on one end.
- Platform is 36" x 36" x 18" high.

MODEL	LENGTH	WIDTH	HEIGHT
1582	58"	36"	18"

### C. Model 1580 ADA Convertible "Space-Saver" Staircase

**SPECIAL FEATURES:**

- Staircase can be installed in straight line or 90° corner arrangement.
- A.D.A. 11" deep Steps x 36" wide have "anti-slip" safety treads.
- (3) 6" high steps on one end, and (4) 4½" high steps on other end.
- Platform is 36" x 36" x 18" high.

MODEL	LENGTH	WIDTH	HEIGHT
1580	91"	36"	18"
Corner	58" x 69"	36"	18"

C. Model 1580



Model 1580  
straight setup