

## Common Features

### BLOOD DRAWING CHAIRS — MODELS 2191, 2196:

- Designed to save time and increase efficiency.
- Gray powder-coated steel frame.
- One-piece contoured Black polypropylene seat 16"W x 16"D x 18"H.
- Armrests are interchangeable and can be set up for left or right hand use.
- One Armrest pivots 360°. The other Armrest is non-pivoting.
- Both Armrests have 10" height range.
- Easy-to-clean Gray laminate finish.
- Shipped ready-to-assemble.

### A. Model 2191 Blood Chair – Standard Height



#### SPECIAL FEATURES:

- Seat is comfortable 19" height.
- Both Armrests adjust from 26" to 36" high with control knobs.

MODEL	LENGTH	DEPTH	HEIGHT
2191	29"	18"	19"

### B. Model 2196 Blood Chair with Drawer



#### SPECIAL FEATURES:

- Same as Model 2191 but equipped with (1) polystyrene supply drawer. 8"W x 19"D x 3½"H.

MODEL	LENGTH	DEPTH	HEIGHT
2196	37½"	18"	19"

A. Model 2191



B. Model 2196



C. Model 2198



### C. Model 2198 Bariatric Blood Chair



#### SPECIAL FEATURES:

- Extra-wide 24½" x 18½" deep x 22" high Black vinyl padded seat and backrest.
- "Flip-up" L-shaped, padded right hand Armrest folds upright 100° for convenient patient access. Folds down flat in use.
- Left Armrest is non-folding. Armrests are not interchangeable. Only available as shown.
- Both armrests are adjustable in height from 27" to 36" high with easy-to-use control knobs.
- Gray powder coated steel legs and frame
- Black vinyl is stocking color.

MODEL	WIDTH	DEPTH	HEIGHT
2198	34"	29"	37½"

### D. Models 2000, 2002

#### Foot Stools



- Durable and attractive Silver Vein powder-coated steel finish.
- Anti-slip treaded 14" x 11" step.
- Metal tubular legs with safety tips.

MODEL	LENGTH	WIDTH	HEIGHT
2000	14"	11"	9"
2002	Same but packed (2) per carton		

D. Model 2000



### E. Models 2030, 2032

#### Foot Stools with Safety Handrail



- Same as Model 2000 but with handrail.
- 33"H steel handrail snaps in without tools.

MODEL	LENGTH	WIDTH	HEIGHT
2030*	14"	11"	33"
2032	Same but packed (2) per carton		

\*Platform is 9" high

E. Model 2030



### F. Models 2216-346, 2216-927 Laminate Foot Stools



- Heavy duty, laminate wood construction in 346 Oak or 927 Folkstone Gray.
- Large 20" x 16" non-slip tread step.

MODEL	LENGTH	WIDTH	HEIGHT
2216-346	20"	16"	8"
2216-927	20"	16"	8"

F. Model  
2216-927



F. Model  
2216-346



SEE BARIATRIC FOOT STOOLS ON PAGE 34.