



MODEL 4062

MODEL 4062
CONVERTIBLE TAPING BENCH/ TREATMENT TABLE



- Multipurpose table functions both as a taping bench with backrest and as a treatment table.
- To use as a taping bench, raise the gas-spring padded backrest up to 75° with recessed control handles on each side of the table.
- Just lower the backrest and the table quickly converts into a full-length treatment table.
- Convenient front full extension pullout for taping supplies.
- Doors have concealed, adjustable, 110° self-closing Euro-hinges.

- (2) 21"W x 24"D open storage sections with adjustable shelves allow access from either side of the table.
- Comfortable 2" thick high-density foam tops in (3) sections – 32"L seating section, 30"L elevating middle section, and 15"L head section.
- 3mm high impact PVC edges on all door and drawer fronts.
- Laminate: Natural Oak.
- Upholstery: Choice of (10) Pro-Form™ vinyl colors.

MODEL	LENGTH	WIDTH	HEIGHT
4062	78"	27"	34"