



A. Model 4343



B. Model 4346



Leaf detail

### Common Features

THERAPY TABLES – MODELS 4343, 4345, 4346, 4347

- Easy-to-operate Crank Hydraulic Lift system adjusts table from 27" to 39" high for sit-down, stand-up or wheelchair use.
- Crank handle folds out of the way when not in use.
- Extruded aluminum telescoping legs with tubular metal base with levelers.



### A. Model 4343 Hand Therapy Table

- Sail White 1 1/8" high-pressure laminate top with T-molding edges has 4" x 24" cutouts in front and back for easy therapist-patient use.

| MODEL | LENGTH | WIDTH | HEIGHT    |
|-------|--------|-------|-----------|
| 4343  | 48"    | 32"   | 27" - 39" |

### B. Model 4345, 4346, 4347 Hi-Lo ADA Expandable Table

- Sail White 1 1/8" thick high pressure laminate tops with T-molding edges.
- Includes (2) 12" Folding Leaves on each end to accommodate (2) additional patients.

| MODEL            | LENGTH    | WIDTH | HEIGHT    |
|------------------|-----------|-------|-----------|
| 4345 (not shown) | 48" - 72" | 32"   | 27" - 39" |
| 4346             | 60" - 84" | 32"   | 27" - 39" |
| 4347 (not shown) | 72" - 96" | 32"   | 27" - 39" |

C. Model 6282



D. Model 6283



E. Model 6284



### Common Features

"SPRING POWER" THERAPY TABLES — MODELS 6282, 6283, 6284:



- Height of "spring-powered" table can be easily adjusted with dual hand control levers.
- 14" height range allows easy use sitting, standing or in wheelchair.
- Folkstone Gray laminate 3/4" tops.
- Black epoxy heavy-duty dual steel uprights and "U" base provide excellent stability.

### C. Model 6282 Hand Therapy Table

- Compact top facilitates use.
- Equipped with (4) leveling legs.

| MODEL | LENGTH | DEPTH | HEIGHT            |
|-------|--------|-------|-------------------|
| 6282  | 32"    | 18"   | 26 1/2" - 40 1/2" |

### D. Model 6283 Hand Therapy Table with Drawer

- Same as 6282 but with Drawer.

| MODEL | LENGTH | DEPTH | HEIGHT            |
|-------|--------|-------|-------------------|
| 6283  | 32"    | 18"   | 26 1/2" - 40 1/2" |

### E. Model 6284 Mobile Therapy Table

- 32" x 24" top has 19"W x 6"D cutout.
- (4) 2" twin-wheel nylon casters.

| MODEL | LENGTH | DEPTH | HEIGHT            |
|-------|--------|-------|-------------------|
| 6284  | 32"    | 24"   | 28 1/2" - 42 1/2" |

### OPTION FOR MODELS 6282, 6283

075 – Set of (4) 2" nylon locking casters (increases height by 2")