

A. MODEL 4061

A. MODEL 4061
TWO-PERSON TAPING TABLE /
CONVERTIBLE TREATMENT TABLE



- Solid (1) piece construction (not modular).
- (2) Padded taping sections 26" wide x 34" deep x 36" high with (1) removable middle top pad.
- Just add middle pad to convert taping table into a treatment table.
- (3) Storage sections with oak laminate interiors and shelf.
- (1) Center drawer on nylon rollers.
- Doors have concealed, adjustable, self-closing Euro-hinges.
- Comfortable 2" thick high-density foam tops.
- 3mm high impact PVC edges on all door and drawer fronts.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak

Upholstery - 10 Pro-Form™ vinyl colors

MODEL	LENGTH	WIDTH*	HEIGHT
4061	74"	34"	36"

*width of cabinet is 30"



Shown as Treatment Table with middle pad in place

OPTIONS

MODELS 4061, 4062

*A945 - Drawer Lock

*A946 - (1) Lock/Bolt for Double Doors

*FACTORY INSTALLED

B. MODEL 4062
CONVERTIBLE TAPING BENCH/
TREATMENT TABLE



- Multipurpose table functions both as a taping bench with backrest and as a treatment table.
- To use as a taping bench, raise the gas-spring padded backrest up to 75° with recessed control handles on each side of the table.
- Just lower the backrest and the table quickly converts into a full-length treatment table.
- Convenient front full extension pullout for taping supplies.
- Doors have concealed, adjustable, 110° self-closing Euro-hinges.
- (2) 21"W x 24"D open storage sections with adjustable shelves allow access from either side of the table.
- 2" thick high-density foam tops in (3) sections - 32"L seating section, 30"L elevating middle section, and 15"L head section.
- 3mm high impact PVC edges on all door and drawer fronts.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak

Upholstery - 10 Pro-Form™ vinyl colors

MODEL	LENGTH	WIDTH	HEIGHT
4062	78"	27"	34"

B. MODEL 4062

