

Advantages

Adjustable Unit Measurement Bar

The mechanical rotating poise bars allow for easy reading of weight in lb or kg

Exclusive 10 Year Limited Warranty

Health o meter® Professional's industry leading warranty protects you in case of damage or repair

Built-in Height Rod

The easy-to-use height rod provides height measurements from 23" to 84"

▼ Powder Coated Fabrication

Powder coating inside and out protects it from rust, extending the life of the scale

Engraved Numerals

Engraved numerals on the poise bars can withstand frequent use, providing long-lasting readability

▼ 500lb Capacity Beam Scale

Markets highest capacity beam scale, providing versatile use





450KL

1

Optional Wheels Item # 55000



Specifications

• Capacity: 500 lb / 200 kg

· Counterweights included

• Graduation: 1/4 lb / 100 g

• Platform Size: 10-1/2" (w) x 14" (d), platform mat: non-slip black

43 lb

892076002169

- Easy to read rotating poise bars for dual reading in lb and kg
- Professional Dual Reading Height Rod: 23" 84" / 60 cm 213 cm
- Optional wheels available (item # 55000)

20.25" x 58.5" x 6.5" / 4.46

- · Factory calibrated to ensure accuracy
- Powder coated, non-rust inside and out for extended life
- 10 Year Limited Warranty

Health o meter®

Professional

We will do everything possible to make it easy for our customers to do business with us.