


Characteristics

On-line chlorination creates excellent donning and low protein content. Pair packaged, inner wrapper in traditional peel-open pouch.



**Latex
Exam Glove
Sterile
Pairs**



DermAssist™

**Powder-Free Latex
Series 104**

PRODUCT DETAILS

SIZE	ITEM NO.	PACKAGING	DESCRIPTION
S	104100	50 pairs/box, 4 boxes/case	Gloves, Exam, Latex, Sterile, Powder-Free, Pairs
M	104200	50 pairs/box, 4 boxes/case	
L	104300	50 pairs/box, 4 boxes/case	
XL	104350	50 pairs/box, 4 boxes/case	

Features

- Online Chlorinated
- Textured
- Low Protein

Benefits

- Softer Feel with Better Grip
- Improved Wet/Dry Grip
- Reduced Chances of Latex Allergic Reactions

Product Solutions You Trust



DermAssist™ Sterile Pairs Latex Examination Gloves

ONE STERILE P AIR
Latex Examination Glove s

Powder-Free ▾ **Ambidextrous**

Item #	Size
104100 (SPF2100)	Small (6.5-7.0")
104200 (SPF2200)	Medium (7.5-8.0")
104300 (SPF2300)	Large (8.5-9.0")
104350 (SPF2350)	X-Large (9.5-10.0")

STERILE EO

Single Use Only

These gloves contain no more than 50 micrograms of protein per gram of glo

CAUTION
This product contains natural rubber latex which may cause allergic reactions. Safe use of these gloves by or on latex sensitized individuals has not been established.
ADVERTENCIA: La seguridad del uso de este guante por o en individuos sensibles al latex no ha sido establecida. Este producto contiene latex de caucho natural y puede causar reacciones alergicas.
MISE EN GARDE: La Sécurité de l'utilisation de ce gant par ou sur des individus sensibles au latex n'a pas été établie. Ce produit contient du latex de caoutchouc naturel pouvant causer des réactions allergiques.

If you suspect that you or your patients are allergic to natural rubber latex, discontinue use immediately and seek the advice of a physician.

LOT NUMBER:
MANUFACTURED:

Sterile Unless Package Is Damaged Or Open

Manufactured in India for Innovative Healthcare Corporation
8222 Chester Page Blvd., Sand Springs, OK 74482

DermAssist is manufactured in compliance with the following international standards

Designation	Standard
ASTM D3578-05	Standard Specification for Rubber Examination Gloves
ASTM D5151-06	Standard Test Method for Detection of Holes in Medical Gloves
ASTM F1671-03	Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens
ASTM D5712-05e1	Standard Test Method for Analysis of Aqueous Extractable Protein in Natural Rubber

Size	Item Number	Total Length +/- 5mm	Thickness Palm +/- 0.03mm	Thickness Finger +/- 0.05mm
Small	104100	240	0.08	0.08
Medium	104200	240	0.08	0.08
Large	104300	240	0.08	0.08
X-Large	104350	240	0.08	0.08

Tensile Strength & Elongation	Before Aging	After Accelerated Aging
Tensile Strength Min. (Mpa)	18	14
ASTM Requirement Min. (Mpa)	18	14
Elongation Minimum Min. (%)	650	500
ASTM Requirement Min. (%)	650	500



Innovative Healthcare Corporation is certified to both ISO 9001:2000 quality management system and ISO 13485:2003 quality management system for medical devices.