



KIMBERLY-CLARK* KC100 Procedure Masks


ASTM F2100-11 Level 1

NEW



All procedure masks are made with So Soft inner layer*

Three Layers of Comfort and Protection



The KIMBERLY-CLARK ADVANTAGE*

KNOWLEDGE NETWORK* Accredited Education
Ongoing Customer Support
Expert Sales Force
Tools & Best Practices
Clinical Research
Commitment to Excellence

Infection prevention website:

www.HAIwatch.com



For more information, please call your sales representative, or visit our website at www.kchealthcare.com.

*Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. All rights reserved. The COLOR LAVENDER is a Trademark of KCWW. ©2011 KCWW. H02460 H1273-11-01

KIMBERLY-CLARK* KC100 Procedure Masks ASTM F2100-11 Level 1 Compliant.

Kimberly-Clark, a world leader in face masks, is the first to offer masks that meet and exceed the ASTM F2100-11 Level I standard. This new standard reinforces the importance of labeling by specific protection levels based on fluid exposure, filtration, and breathability. This greatly simplifies the process of mask selection. The ASTM F2100-11 standard is endorsed by APIC and AORN.

So Soft* Mask Comfort

KIMBERLY-CLARK* KC100 Procedure Masks are made with an inner layer featuring So Soft* bicomponent fabric technology. It's wonderfully breathable and designed to keep you comfortable.

A Rainbow of Protection

Choose from three pastel colors—Yellow, Lavender, or Blue. All offer the same ASTM F2100-11 Level I protection for low fluid exposure.

Test Method	Test Description	ASTM 2100-11 Requirement Level 1	KC100 New Mask*
ASTM F1862	Fluid Resistance with synthetic blood	80 mm Hg	80 mm Hg
MIL-M-36954C	Delta P	< 4.0 mm H ₂ O	3.0 mm H ₂ O
ASTM F2101	Bacterial Filtration Efficiency (BFE)	≥ 95%	99.7%
ASTM F2299	Submicron Particle Filtration (PFE)	≥ 95% @ 0.1 micron	98.4%
16 CFR Part 1610	Flammability	Class 1	Class 1

*Based on the calculated lower 99% confidence interval about the average

KIMBERLY-CLARK* Procedure Mask ASTM F2100-11

NEW Item #	Product Description	ASTM F2100-11 LEVEL 1	Packaging
25867	KC100 Procedure Mask	Yellow	500 eaches/case 10 cartons of 50 each
25868	KC100 Procedure Mask	Lavender	500 eaches/case 10 cartons of 50 each
25869	KC100 Procedure Mask	Blue	500 eaches/case 10 cartons of 50 each



Kimberly-Clark

Trusted Clinical Solutions*