



## Under Buttocks Drape with Fluid Collection Pouch



1. Orient Under Buttocks Drape with patient and insert right and left hands as marked.




2. Slide drape under patient's buttocks.

(89414 depicted within OB IV Pack)



3. Draping is complete.



### Features

- 40" x 44"
- Fluid Collection Pouch
- Absorbent Pad

### Stock #

### Description

### Items Per Case

89414

Under Buttocks Drape with Fluid Collection Pouch,  
Sterile, 1 per package

30

79414 ♦

Under Buttocks Drape with Fluid Collection Pouch,  
Non-Sterile, Bulk Pack

80

♦ Except as noted, all non-sterile products are for custom trays and available only through your custom tray supplier.

† As with any drape or pack recommendation, it is the hospital's decision as to the selection of the appropriate draping technique.

Copyright 2002 Kimberly-Clark Corporation. All Rights Reserved. Your visit to this site and use of the information hereon is subject to the terms of our Legal Statement. \*Registered Trademark or Trademark of Kimberly-Clark Corporation or its affiliates.

<See the following page for Pack Options>

Pack Stock # 88841		Description <i>OB Pack IV</i>	Packs/Case 5
Component Stock #			
<b>Sterile</b>	Non-Sterile ♦	<b>Contents:</b>	
<b>89618</b>	79498	1 Back Table Cover, Zone-Reinforced, 44" x 90"	
		1 Set-Up Cover, 76" x 44"	
<b>95411</b>	70411	1 ULTRA Impervious Surgical Gown, Large, with Towel	
<b>89701</b>	79720	2 Absorbent Towels, 15" x 26"	
		1 Baby Blanket	
<b>89414</b>	79414	1 Under Buttocks Drape, 40" x 44", with Fluid Collection Pouch <ul style="list-style-type: none"> <li>• Absorbent Pad</li> </ul>	
<b>89408</b>	79408	2 Leggings, with 6" Cuffs	
	79810	1 Abdominal Drape, 60" x 38" <ul style="list-style-type: none"> <li>• CONTROL* PLUS Fabric Reinforcement</li> </ul>	
		1 Peri Pad	
		1 HOLLISTER* Cord Clamp	
		1 Vaginal Packing Sponge	
		1 DAVOL* Bulb Syringe – 2 oz.	
		10 Sponges, X-Ray Detectable, 4" x 4"	
		1 Placenta Basin with Lid	

♦ Except as noted, all non-sterile products are for custom trays and available only through your custom tray supplier.

Copyright 2002 Kimberly-Clark Corporation. All Rights Reserved. Your visit to this site and use of the information hereon is subject to the terms of our Legal Statement. \*Registered Trademark or Trademark of Kimberly-Clark Corporation or its affiliates.