



100 SERIES SHOWER CHAIRS

MJM ADVANTAGES:

- *Healthcare Grade PVC Polymer Plastics*
- *All Products Are Extensively Tested*
- *Push/Pull Handle*
- *Deluxe Rust Proof Threaded Stem Casters*
- *Reinforced At Stress Related Areas*
- *Deluxe Elongated Open Front Seat*
- *Anti-Slip Safety Handgrips*
- *Fast Drying Removable Mesh Sling*
- *Personal Pouch Ideal to Store Hygiene Items While Showering*



**5 YEAR LIMITED
WARRANTY ON PVC FRAME**



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

Yellow	Tan	Mauve	F. Green	Red	Black	Grey	L. Blue	R. Blue	White	Navy	Mint
--------	-----	-------	----------	-----	-------	------	---------	---------	-------	------	------

ALL SOLID VINYL FABRICS ARE ANTIMICROBIAL, FLAME RESISTANT, HIGH TEAR AND TENSILE STRENGTH AND EXCELLENT CLEANABILITY

VINYL ONLY



With Healthcare Facility Support, MJM has Developed an Innovative Built in Push/Pull Handle that will Assist and Facilitate All of MJM's Shower Chairs to Maneuver with Ease.

118-3 SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, PERSONAL POCKET POUCH.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS

MJM Highly Suggest Adding One or More of Our Top Selling Shower Chair Options for SAFETY-CONVENIENCES and FUNCTION
 (1) Sliding Footrest #SF (2) Commode Pail #SQ Pail (3) Safety Belt #BB
 (4) Dual Drop Arms #DDA (5) Deluxe Elongated Soft Seat #SSDE



118-3-VS SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE EXTRA HEAVY DUTY VACUUM SEAT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Clamp on Seat: 20"
- Seat Height Clearance: 17.75"
- Seat Height Clearance with Commode Pail Rails: 16.5"
- Capacity: 300 lbs
- Shipped: UPS

Full Support Snap On Seat Provides Better Comfort

- Snap On Seat Allows for Easy Cleaning and Disinfecting
- Unique Built-In Splash Guard

MJM Use HEALTHCARE GRADE Mesh / Vinyl Fabrics That Are Flame Resistant, Mildew Resistant and Have Biocide or Antibacterial Additives.



118-3-FS SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, FLARED BASE FOR EXTREME STABILITY.

- Flared Base: 4" Wider Base than Standard 118-3
- Internal Width: 18"
- External Width: 26"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 25" W x 28" D
- Seat Height: 20"
- Capacity: 300 lbs.
- Shipped: UPS

MJM'S HEALTHCARE GRADE POLYMER PLASTICS INHIBIT/CONTROL THE GROWTH OF BACTERIA, IS PHTHALATES AND LEAD FREE.

18

118-3-FS-R0S
REVERSE OPEN SEAT



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES



Fern Dance Floral



Montego Floral



Sierra Sands Tweed

DESIGNER MESH FABRICS

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

118-3-KD

KNOCKDOWN SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, USER FRIENDLY SELF STORING FOLDING FOOTREST PREVENTS FOOT DRAGGING AND ENTRAPMENT, REINFORCED AT ALL STRESS RELATED AREAS, SQUARE PAIL.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22.5" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS



OPTIONAL:
118-KD-BAG

**SHOWER CHAIR
IN THE BAG**



118-3-SAFE

SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS. OUTRIGGERS HELP PREVENT UNIT FROM TIPPING OVER.

- Internal Width: 18"
- External Width: 22"
- External Depth, Front to Back at Bottom Section of Outriggers: 30.25"
- External Width, Left to Right at Bottom Section of Outriggers: 31.75"
- Threaded Stem Casters: Rust Proof 3" Twin
- Overall Size: 40.5" H x 31.75" W x 30.25" D
- Seat Height: 20"
- Capacity: 300 lbs.

**360° Outrigger
for Safe and
Stable Swivel
Movement
8 Casters (4 Locking)**



MJM Highly Suggest Adding One or More of Our Top Selling Shower Chair Options for SAFETY-CONVENIENCES and FUNCTION
(1) Sliding Footrest #SF (2) Commode Pail #SQ Pail (3) Safety Belt #BB
(4) Dual Drop Arms #DDA (5) Deluxe Elongated Soft Seat #SSDE

118-3-ADJ

ADJUSTABLE SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, STANDARD DELUXE ADJUSTABLE HEIGHTS.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Fixed Height: 40.5" H x 22.5" W x 25.25" D
- Adjusted Height: 45.5" H x 22.5" W x 18" D
- Height Top of Seat: 20"
- Height Top of Seat Adjustable: 25"
- Height Clearance Under Seat: 17.5"-22.5"
- Height Adjustment Increments: 6 each of 1"
- Capacity: 250 lbs.

**ADJUSTABLE HEIGHT ROLLING SHOWER
COMMODE CHAIR WILL ACCOMMODATE
STANDARD TO HANDICAP TOILETS**



www.mjminternational.com



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

Yellow	Tan	Mauve	F. Green	Red	Black	Grey	L. Blue	R. Blue	White	Navy	Mint

ALL SOLID VINYL FABRICS ARE ANTIMICROBIAL, FLAME RESISTANT, HIGH TEAR AND TENSILE STRENGTH AND EXCELLENT CLEANABILITY



118-3-OF SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT AND INCREASES CAREGIVERS ACCESS, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Height: 20"
- Capacity: 300 lbs.
- Shipped: UPS

**MJM'S HEALTHCARE GRADE
POLYMER PLASTICS INHIBIT/CONTROL THE GROWTH
OF BACTERIA, IS PHTHALATES AND LEAD FREE.**



118-3-H

MJM ADVANTAGES: HEALTHCARE GRADE, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, HIGHBACK, ASSISTS PATIENTS THAT NEED ADDITIONAL BACK SUPPORT.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 47.5" H x 22" W x 25.25" D
- Seat Height: 20"
- Top of Seat to of Back: 27.5"
- Capacity: 300 lbs.
- Shipped: UPS

MJM Use HEALTHCARE GRADE Mesh / Vinyl Fabrics That Are Flame Resistant, Mildew Resistant and Have Biocide or Antibacterial Additives.

MJM Highly Suggest Adding One or More of Our Top Selling Shower Chair Options for SAFETY-CONVENIENCES and FUNCTION
(1) Sliding Footrest #SF (2) Commode Pail #SQ Pail (3) Safety Belt #BB
(4) Dual Drop Arms #DDA (5) Deluxe Elongated Soft Seat #SSDE



118-3-TS SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, SHOWER CHAIR W/TILT SEAT W/OPTIONAL 12 QT. COMMODE PAIL (W/HANDLE).

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Front Seat Height: 20"
- Back Seat Height: 18"
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 16.625"
- Seat Height Clearance with Commode Pail Rails: 14.75"
- Capacity: 300 lbs.
- Shipped: UPS



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES



Fern Dance Floral



Montego Floral



Sierra Sands Tweed

DESIGNER MESH FABRICS

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

118-3-FD

MJM ADVANTAGES: HEALTHCARE GRADE, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, STANDARD DELUXE W/FOLDING CAPACITY.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 24" D
- Seat Height: 20.25"
- Seat Height Clearance: 17.5"
- Folded Width: 6"
- Folded Length: 39"
- Capacity: 200 lbs.
- Shipped: UPS

**NO MORE CROWDED
SHOWER ROOMS**



Available In 16"
116-3-FD

MJM Highly Suggest Adding One or More of Our Top Selling Shower Chair Options for SAFETY-CONVENIENCES and FUNCTION
(1) Sliding Footrest #SF (2) Commode Pail #SQ Pail (3) Safety Belt #BB
(4) Dual Drop Arms #DDA (5) Deluxe Elongated Soft Seat #SSDE

118-3-F-LSB-18

MJM ADVANTAGES: HEALTHCARE GRADE, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, SHOWER CHAIR W/FLAT STOCK SEAT AND LAP SECURITY BAR.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Height: 19.75"
- Capacity: 300 lbs.
- Shipped: UPS



MJM Highly Suggest Adding One or More of Our Top Selling Shower Chair Options for SAFETY-CONVENIENCES and FUNCTION
(1) Sliding Footrest #SF (2) Commode Pail #SQ Pail (3) Safety Belt #BB
(4) Dual Drop Arms #DDA (5) Deluxe Elongated Soft Seat #SSDE

118-3-SQ-PAIL

SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, DELUXE SHOWER CHAIR W/OPTIONAL SQUARE PAIL.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS



**5 YEAR LIMITED
WARRANTY ON PVC FRAME**



www.mjminternational.com



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

Yellow	Tan	Mauve	F. Green	Red	Black	Grey	L. Blue	R. Blue	White	Navy	Mint

ALL SOLID VINYL FABRICS ARE ANTIMICROBIAL, FLAME RESISTANT, HIGH TEAR AND TENSILE STRENGTH AND EXCELLENT CLEANABILITY



SEAT IN DOWN POSITION

118-3-FLS-SQ-PAIL

SHOWER CHAIR W/FLIP FRONT SEAT

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FLIP FRONT SEAT ENHANCES RESIDENT COMFORT AND INCREASES CARE-GIVERS ACCESS, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS

MJM Highly Suggest Adding One or More of Our Top Selling Shower Chair Options for SAFETY-CONVENIENCES and FUNCTION
 (1) Sliding Footrest #SF (2) Commode Pail #SQ Pail (3) Safety Belt #BB
 (4) Dual Drop Arms #DDA (5) Deluxe Elongated Soft Seat #SSDE

118-3-FS-SQ-PAIL

SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, FLARED STABILITY BASE PROVIDES FOR EXTREME STABILITY W/OPTIONAL SQUARE PAIL.

- Internal Width: 18"
- External Width: 26"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 25" W x 28" D
- Seat Height: 20"
- Capacity: 300 lbs.
- Shipped: UPS



MJM Use HEALTHCARE GRADE Mesh / Vinyl Fabrics That Are Flame Resistant, Mildew Resistant and Have Biocide or Antibacterial Additives.

MJM'S HEALTHCARE GRADE POLYMER PLASTICS INHIBIT/CONTROL THE GROWTH OF BACTERIA, IS PHTHALATES AND LEAD FREE.

118-3-SF-SQ-PAIL

SHOWER CHAIR W/SLIDING FOOTREST AND COMMODE PAIL

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, USER FRIENDLY SELF STORING SLIDING FOOTREST PREVENTS FOOT DRAGGING AND ENTRAPMENT, REINFORCED AT ALL STRESS RELATED AREAS, SQUARE PAIL.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS



118-3-FF-SQ-PAIL

Also Available with Folding Footrest



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES



Fern Dance Floral



Montego Floral



Sierra Sands Tweed

DESIGNER MESH FABRICS

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

118-3-SF

SHOWER CHAIR W/SLIDING FOOTREST

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, USER FRIENDLY SELF STORING FOLDING FOOTREST PREVENTS FOOT DRAGGING AND ENTRAPMENT, REINFORCED AT ALL STRESS RELATED AREAS.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS



118-3-FF

Also Available with Folding Footrest



118-3-SSDE-IF

SHOWER CHAIR W/INDIVIDUAL FOOTREST

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SOFT SEAT ENHANCES RESIDENT COMFORT, INDIVIDUAL FOOTREST, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 25.5" W x 31.75" D
- Seat Height: 21.25"
- Capacity: 300 lbs.
- Shipped: UPS

MJM Deluxe Elongated Soft Seat Option Adds Absolute Comfort to Any of Our Shower Chairs. MJM Soft Seat is Designed with a Combination Durable Marine and Health Care Grade Material Cover and Resilient Foam and 3/4" Treated Plywood for the Ultimate Durability and Patient Comfort.



118-3-SSC-IF



118-3-VS-SFS-DDA-10QTC

SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE EXTRA HEAVY DUTY VACUUM SEAT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, USER FRIENDLY SELF STORING SLIDING FOOTREST SUPPORT PREVENTS FOOT DRAGGING AND ENTRAPMENT, REINFORCED AT ALL STRESS RELATED AREAS, DUAL DROP ARM.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 23.5" D
- Seat Depth: 18"
- Seat Height with Clamp on Seat: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.75"
- Seat Height Clearance with Commode Pail Rails: 16.5"
- 10 QT SQ Commode Pail
- Sliding Footrest Support
- Double Drop Arm
- Capacity: 300 lbs

Full Support Snap On Seat Provides Better Comfort

- Snap On Seat Allows for Easy Cleaning and Disinfecting
- Unique Built-In Splash Guard

MJM Use HEALTHCARE GRADE Mesh / Vinyl Fabrics That Are Flame Resistant, Mildew Resistant and Have Biocide or Antibacterial Additives.



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH											VINYL ONLY
Yellow	Tan	Mauve	F. Green	Red	Black	Grey	L. Blue	R. Blue	White	Navy	Mint
ALL SOLID VINYL FABRICS ARE ANTIMICROBIAL, FLAME RESISTANT, HIGH TEAR AND TENSILE STRENGTH AND EXCELLENT CLEANABILITY											



118-3-SSDD

SHOWER CHAIR W/SOFT SEAT DELUXE DUAL USAGE W/REMOVABLE CENTER SECTION

MJM ADVANTAGES: HEALTHCARE GRADE, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, PERSONAL POCKET POUCH.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS

MJM Use HEALTHCARE GRADE Mesh / Vinyl Fabrics That Are Flame Resistant, Mildew Resistant and Have Biocide or Antibacterial Additives.



118-3-SSDE-BL

SHOWER CHAIR W/SOFT SEAT DELUXE ELONGATED

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, PERSONAL POCKET POUCH.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS

118-3-SSDE
SHOWER CHAIR W/SOFT SEAT DELUXE ELONGATED

MJM Deluxe Elongated Soft Seat Option Adds Absolute Comfort to Any of Our Shower Chairs. MJM Soft Seat is Designed with a Combination Durable Marine and Health Care Grade Material Cover and Resilient Foam and 3/4" Treated Plywood for the Ultimate Durability and Patient Comfort.

118-3-SSC

SHOWER CHAIR W/SOFT SEAT COMPLETE

MJM ADVANTAGES: HEALTHCARE GRADE, SOFT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, PERSONAL POCKET POUCH.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75"
- Capacity: 300 lbs.
- Shipped: UPS

MJM Deluxe Elongated Soft Seat Option Adds Absolute Comfort to Any of Our Shower Chairs. MJM Soft Seat is Designed with a Combination Durable Marine and Health Care Grade Material Cover and Resilient Foam and 3/4" Treated Plywood for the Ultimate Durability and Patient Comfort.



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES



Fern Dance Floral



Montego Floral



Sierra Sands Tweed

DESIGNER MESH FABRICS

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

118-3-0S-SS0S

SHOWER CHAIR W/SOFT SEAT OPEN SEAT

MJM ADVANTAGES: HEALTHCARE GRADE, SPLIT SOFT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, PERSONAL POCKET POUCH.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 25.25" D
- Seat Depth: 18"
- Seat Height with Open Front: 20"
- Seat Height with Split Soft Seat: 21.25"
- Capacity: 300 lbs.
- Shipped: UPS



TOTAL HYGIENE CHAIR



Patent Pending

The Total Hygiene Chair is a Must in Every Healthcare Facility. The Dual Split Soft Seat and the Dual Railing System Gives and Helps Aids Washing/Cleaning the Perineal Area on Patients.

RIO 118-3-TL-SFS-SQPAIL-SSDE-SADJ-BB-18

ADJUSTABLE SHOWER COMMODE CHAIR

IDEAL FOR ADULT

SEAT HEIGHT ADJUST FROM 16" – 25" FITS PATIENTS 5' TO 6'3"



MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE SOFT SEAT, USER FRIENDLY SELF STORING SLIDING FOOTREST SUPPORT PREVENTS FOOT DRAGGING AND ENTRAPMENT, REINFORCED AT ALL STRESS RELATED AREAS, DELUXE ADJUSTABLE SHOWER COMMODE CHAIR IS CAREGIVER AND PATIENT FRIENDLY. THE RIO CHAIR IS IDEAL, BECAUSE IT WILL ACCOMMODATE DIFFERENT PATIENTS HEIGHTS, ERGONOMIC PUSH PULL HANDLE, BREATHABLE MESH BACK, COMES WITH PERSONAL POUCH TO STORE TOOTH BRUSH, COMB ETC., SAFETY BELT, WITH SIDE RELEASE BUCKLES FOR QUICK RELEASE, SAFETY GRIPPERS FOR ARM – HAND SECTION, FOUR RUST PROOF TOTAL LOCK CASTER, DELUXE OPEN FRONT SOFT SEAT, IDEAL FOR BOTH MALE AND FEMALE, SQUARE COMMODE PAIL WITH BUILT IN HANDLE, COMMODE PAIL COVERS THE ENTIRE SOFT SEAT, HEIGHT ADJUSTABLE SEAT FROM 16" – 25" AT TOP OF SEAT (ADJUST WITH SPRING LOADED KNOB), ADJUSTMENT RAILS ARE BUILT INTO THE CHAIR FRAME AND ARE MADE OF HIGH IMPACT PLASTIC (EASY TO CLEAN AND DISINFECT).

- Internal Width: 18"
- External Width: 27"
- Seat Depth: 18"
- External Depth of Frame: 20"
- Sliding Footrest Support



MJM Deluxe Elongated Soft Seat Option Adds Absolute Comfort to Any of Our Shower Chairs. MJM Soft Seat is Designed with a Combination Durable Marine and Health Care Grade Material Cover and Resilient Foam and 3/4" Treated Plywood for the Ultimate Durability and Patient Comfort.



Lowest Seat Height - 16"

Highest Seat Height - 25"



www.mjminternational.com

25



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

Yellow	Tan	Mauve	F. Green	Red	Black	Grey	L. Blue	R. Blue	White	Navy	Mint
--------	-----	-------	----------	-----	-------	------	---------	---------	-------	------	------

ALL SOLID VINYL FABRICS ARE ANTIMICROBIAL, FLAME RESISTANT, HIGH TEAR AND TENSILE STRENGTH AND EXCELLENT CLEANABILITY

VINYL ONLY



118-3-AF

MJM ADVANTAGES: HEALTHCARE GRADE, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, FLAT STOCK SEAT WITH DRAIN HOLES, FLAT STOCK SEAT SUPPORTS ENTIRE BUTTOCKS, STANDARD DELUXE AMPUTEE FLATSTOCK.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 27.75" D
- Elevated Leg Extension (Padded): 14.5" L
- Seat Height: 19.75"
- Capacity: 300 lbs.
- Shipped: UPS

THE ONLY PVC SHOWER CHAIR FOR UNI-LATERAL OR BI-LATERAL BELOW KNEE AMPUTEE

FLATSTOCK SEATS WILL PROVIDE A WIDER AND MORE COMFORTABLE SEAT FOR RESIDENTS THAT DO NOT NEED OR LIKE THE TRADITIONAL PVC SHOWER CHAIRS WITH A COMMODE SEAT.



118-3-A

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, THE ONLY PVC SHOWER COMMODE CHAIR FOR UNI-LATERAL OR BI-LATERAL BELOW KNEE AMPUTEE.

- Internal Width: 18"
- External Width: 22"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 22" W x 27.75" D
- Elevated Leg Extension (Padded): 14.5" L
- Seat Height: 21"
- Capacity: 300 lbs.
- Shipped: UPS

MJM Use HEALTHCARE GRADE Mesh / Vinyl Fabrics That Are Flame Resistant, Mildew Resistant and Have Biocide or Antibacterial Additives.

MJM'S HEALTHCARE GRADE POLYMER PLASTICS INHIBIT/CONTROL THE GROWTH OF BACTERIA, IS PHTHALATES AND LEAD FREE.



MJM'S SUPERIOR SHOWER CHAIRS OFFER COMFORT AND SAFETY

100 SERIES



Fern Dance Floral



Montego Floral



Sierra Sands Tweed

DESIGNER MESH FABRICS

BREATHABLE STANDARD MESH IS FLAME RESISTANT, MILDEW RESISTANCE AND HIGH TEAR AND TENSILE STRENGTH

122-3-KD

KNOCKDOWN SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, WIDE DELUXE SHOWER CHAIR.

- Internal Width: 22"
- External Width: 26"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 26.5" W x 27" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75
- Capacity: 375 lbs.
- Shipped: UPS



OPTIONAL:
122-KD-BAG
SHOWER CHAIR
IN THE BAG



122-3

SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, WIDE DELUXE SHOWER CHAIR.

- Internal Width: 22"
- External Width: 26"
- Threaded Stem Casters: Rust Proof 3" Twin
- Over All Size: 40.5" H x 26" W x 27" D
- Seat Height: 21"
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 20"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 21.25"
- Seat Height Clearance: 17.625"
- Seat Height Clearance with Commode Pail Rails: 16.75
- Capacity: 375 lbs.
- Shipped: UPS

MJM Use HEALTHCARE GRADE Mesh / Vinyl Fabrics That Are Flame Resistant, Mildew Resistant and Have Biocide or Antibacterial Additives.



MJM Highly Suggest Adding One or More of Our Top Selling Shower Chair Options for SAFETY-CONVENIENCES and FUNCTION

- (1) Sliding Footrest #SF (2) Commode Pail #SQ Pail (3) Safety Belt #BB
(4) Dual Drop Arms #DDA (5) Deluxe Elongated Soft Seat #SSDE

126-4-NB

BARIATRIC SHOWER CHAIR

MJM ADVANTAGES: HEALTHCARE GRADE, DELUXE ELONGATED OPEN FRONT SEAT ENHANCES RESIDENT COMFORT, ANTI-SLIP HANDGRIPS, PUSH/PULL HANDLE AND FAST DRYING REMOVABLE MESH BACK/SLING, REINFORCED AT ALL STRESS RELATED AREAS, EXTRA-WIDE DELUXE SHOWER CHAIR.

- Internal Width: 26"
- External Width: 30"
- Threaded Stem Casters: Rust Proof 4" Twin
- Over All Size: 40" H x 30" W x 29.5" D
- Seat Depth: 18"
- Seat Height with Deluxe Elongated Open Front: 19.25"
- Seat Height with Deluxe Elongated Open Front Soft Seat: 20.75"
- Seat Height Clearance: 17"
- Seat Height Clearance with Commode Pail Rails: 16"
- Capacity: 425 lbs.
- Shipped: UPS



www.mjminternational.com



DELUXE USER-FRIENDLY OPTIONS FOR RESIDENT COMFORT

100 SERIES

FF: FOLDING FOOTREST

- From Floor To Top of FF: 8.5" (w/Standard 3" Casters)
- Plate Size: 14.5" x 8.5"



IF: INDIVIDUAL FOLDING FOOTREST



SF: SLIDING FOOTREST

- From Floor To Top of SF: 6.5" (w/Standard 3" Casters)
- Plate Size: 14" x 9"



DDA: DOUBLE DROP ARM

The Arms Swing Clear Of The Frame To Allow For Easy Patient Transfers.



LSB: Lap Security Bar

AVAILABLE IN THE FOLLOWING SIZES:

- LSB-15: 15" Internal Shower Chair
- LSB-18: 18" Internal Shower Chair
- LSB-22: 22" Internal Shower Chair
- LSB-26: 26" Internal Shower Chair
- LSB-30: 30" Internal Shower Chair



BB-18

BUCKLE SAFETY BELT

SAFETY BELTS INCREASES RESIDENT SAFETY, ADJUSTABLE SAFETY BELTS FOR MAXIMUM COMFORT.

SIDE RELEASE BUCKLE SIZES:

- BB-15: 15" Internal Shower Chair
- BB-18: 18" Internal Shower Chair
- BB-22: 22" Internal Shower Chair
- BB-26: 26" Internal Shower Chair
- BB-30: 30" Internal Shower Chair

VELCRO CLOSURE SIZES:

- BV-15: 15" Internal Shower Chair
- BV-18: 18" Internal Shower Chair
- BV-22: 22" Internal Shower Chair
- BV-26: 26" Internal Shower Chair
- BV-30: 30" Internal Shower Chair



MJM Use HEALTHCARE GRADE Mesh / Vinyl Fabrics That Are Flame Resistant, Mildew Resistant and Have Biocide or Antibacterial Additives.

PS: PRIVACY SOLID VINYL SKIRT

FULL LENGTH PRIVACY SKIRT

AVAILABLE IN THE FOLLOWING SIZES:

- PS-15: 15" Internal Shower Chair
- PS-18: 18" Internal Shower Chair
- PS-22: 22" Internal Shower Chair
- PS-26: 26" Internal Shower Chair
- PS-30: 30" Internal Shower Chair

ALSO AVAILABLE:

3/4 LENGTH PRIVACY SKIRT

AVAILABLE IN THE FOLLOWING SIZES:

- PS-15-S: 15" Internal Shower Chair
- PS-18-S: 18" Internal Shower Chair
- PS-22-S: 22" Internal Shower Chair
- PS-26-S: 26" Internal Shower Chair
- PS-30-S: 30" Internal Shower Chair



COV-18

THROW COVER

Conceals The Entire Frame.

AVAILABLE IN THE FOLLOWING SIZES:

- COV-15: 15" Internal Shower Chair
- COV-18: 18" Internal Shower Chair
- COV-22: 22" Internal Shower Chair
- COV-26: 26" Internal Shower Chair
- COV-30: 30" Internal Shower Chair



DELUXE USER-FRIENDLY OPTIONS FOR RESIDENT COMFORT

100 SERIES



PUSH/PULL HANDLE



DDA: DUAL
DROP ARM

LSB-18: LAP SECURITY BAR

PS-18: PRIVACY SKIRT

BB-18: BUCKLE SAFETY BELT

DEO: DELUXE ELONGATED
COMMODOE OPEN
FRONT SEAT

SQ-PAIL: COMMODOE PAIL
(Covers The Entire Underside Of The
Deluxe Elongated Commode Open Front
Seat To Help Prevent Splashing And
Spillage.)

SFS: SLIDING FOOTREST
SUPPORT

MJM'S HEALTHCARE GRADE
POLYMER PLASTICS INHIBIT/CONTROL THE GROWTH
OF BACTERIA, IS PHTHALATES AND LEAD FREE.

**HEAVY DUTY THREADED
STEM RUST PROOF CASTERS**

5 YEAR WARRANTY ON FRAME
ALL PRODUCTS EXTENSIVELY TESTED

**REINFORCED AT ALL STRESS
RELATED AREAS**



**5 YEAR LIMITED
WARRANTY ON PVC FRAME**



www.mjminternational.com

