

Latex-Free

Orthopedic Stockinettes

Tubular stockinettes offer conforming, soft fabric to protect patient skin under plaster or synthetic casts. Also used to provide consistent pressure for wound dressings. Continuous rolls are packaged in a dispensing box. Knitted from cotton or synthetic materials in 2, 3 or 4 inch widths.

- Non-sterile
- Dispensing box



Item #	Size	Case
COTTON P072502 P072503 P072504	2" x 25 yd 3" x 25 yd 4" x 25 yd	12 12 10
SYNTHETIC P072522 P072523 P072524	2" x 25 yd 3" x 25 yd 4" x 25 yd	12 12 10

Distributed By: