

Pro
ADVANTAGE[®]

by **NDC**

Powder-Free and Powdered

Exam Gloves

A comprehensive line of hand protection with various fit and function options. These exam gloves offer safety, flexibility, and consistency in the healthcare setting.

Latex

These gloves are comfortable, low protein, and made from a natural fiber. Available in powdered and powder-free.



Soft Nitrile

The microtextured, soft nitrile gloves provide enhanced comfort and tactile sensitivity. These soft and flexible latex-free gloves stretch and conform to the hands, eliminating the stiffness sometimes associated with nitrile gloves. Available in powder-free only.

Vinyl

The non-latex exam gloves provide a durable solution for professionals and their patients who risk sensitivity to natural rubber latex. These gloves offer exceptional strength and protection at a value price. Available in powdered and powder-free.



Powder-Free and Powdered Exam Gloves

Item #	Size	Color	Box
--------	------	-------	-----

LATEX EXAM GLOVES, POWDER-FREE

P359101	X-Small	White	100
P359102	Small	White	100
P359103	Medium	White	100
P359104	Large	White	100
P359105	X-Large	White	90



LATEX EXAM GLOVES, POWDERED

P359201	X-Small	White	100
P359202	Small	White	100
P359203	Medium	White	100
P359204	Large	White	100
P359205	X-Large	White	90

SOFT NITRILE EXAM GLOVES, POWDER-FREE

P359011	X-Small	Blue	200
P359012	Small	Blue	200
P359013	Medium	Blue	200
P359014	Large	Blue	200
P359015	X-Large	Blue	200



VINYL EXAM GLOVES, POWDER-FREE

P359401	X-Small	Clear	100
P359402	Small	Clear	100
P359403	Medium	Clear	100
P359404	Large	Clear	100
P359405	X-Large	Clear	100



VINYL EXAM GLOVES, POWDERED

P359501	X-Small	Clear	100
P359502	Small	Clear	100
P359503	Medium	Clear	100
P359504	Large	Clear	100
P359505	X-Large	Clear	100

Distributed By:

Pro
ADVANTAGE[®]
by **NDC**