

Digital

Thermometer Kit

This dual scale thermometer provides a digital readout in approximately 60 seconds. Beeps when measurement is complete. Can be used for oral, rectal or axillary measurement. Convenient automatic shut off and last reading memory. Moisture resistant. Meets ASTM E1112 standards. Comes with clear storage case and five sheaths. Two year warranty.

- Beeps when measurement is complete
- Auto off conserves battery life
- Recalls last reading



- F°/C° Switchable for any application
- Water-resistant to allow cleaning without risk of damage

Item #	Description
P541222	Digital Thermometer Kit

Distributed By: