

ri-former® Wall mounted diagnostic station “PLUG AND PLAY”



A truly modular “**PLUG AND PLAY**” proposition with virtually limitless potential for diagnostic screening.
Maintenance-free diagnosis at your finger tips!

- Maintenance-free high performance LED close to daylight (5,500k) for more efficient diagnosis.
- Fiber-Optic, Xenon or LED Illumination Options.
- Optimal Beaming and Transmission of Light.
- Bi-Directional Swiveling.
- Sealing Optical Lens with 3-Fold and 4-Fold Magnification.
- Easy to Operate **SpecEjec®** Touch-less Specula Ejection Device.
- Anti-Theft Proof.
- Suitable for Pneumatics Tests.

- Maintenance-free high performance LED close to daylight (5,500K) for more efficient diagnosis.
- Fiber-Optic, Xenon or LED Illumination Options.
- Diopter Wheel with 29 or 89 Corrective Lenses,
- 6 Aperture Hand-Wheel [Semi-Circle, Small/Medium/Large Circle, Fixation Star, and Slit].
- 4 Filter Wheel For All Apertures [Symbol Display, Red-Free Filter, Blue Filter and Polarization Filter].
- Anti-Theft Proof, Dust-Proof.
- Parallel Beam Path [Coaxial].

- Fast precise measurement of the body temperature in the ear (1-2 seconds).
- Digital display for optimum readability of the measurement values in °C or in °F.
- Memory functions for 12 measurements.



- Clinically validate for accurate results (±3mmHg).
- Latex-Free Velcro Cuff with 10' of Latex-Free Coiled Tubing.
- Total Swivel Radius 130°.
- Shock resistant.
- Large and Sturdy for Tough Everyday Situations.

- Extendable up to 5 handles.
- 3 years warranty including LED.
- Made in Germany.

- Specula dispenser with 5 tubes for up to 5 specula sizes.
- Hygienic and practical solution dispenses a new speculum for each examination.
- Simple to clean and disinfect.
- Disposable and reusable specula available Pneumatic Bulb.

**Human Operation
Otoscope Head**
3.5V Xenon Illumination
Anti-Theft Proof



**Tongue Blade
Holder Head**
Fiber-Optic 3.5V
LED Illumination
Anti-Theft Proof



Nasal Speculum Head
Fiber-Optic 3.5V
LED Illumination
Anti-Theft Proof



Dermatoscope Head
Fiber-Optic 3.5V
LED Illumination
Anti-Theft Proof



Light Conductor Head
Fiber-Optic 3.5V
LED Illumination
Anti-Theft Proof



Pneumatic Bulb



Riester
Made in Germany

Riester
Authorized Distributor
USA