

CAPJECT®

Need samples?

We've now made it faster & easier
for you to obtain samples and literature of our products
through your local TMP rep.
Just log-on to www.TerumoTMP.com and select
"Find Your Rep" from our homepage.

CAPJECT®

Capillary Micro Collection Tubes
& Safety Lancet • Blade



The Perfect Pairing of Ease and Precision.



CAPJECT® Capillary Micro Collection Tubes & Safety Lancet • Blade

©2008 All rights reserved. Capject is a registered trademark of Terumo Medical Corporation, Somerset, NJ. Circle T and Terumo are registered trademarks of Terumo Corporation, Japan. BD and Genie are trademarks of Becton Dickinson Co. and are not affiliated with Terumo or Capject.

Capiject Safety Lancets.

The same precision, convenience, and reliability as BD Genie®. Because they are the very same safety lancets.

A legacy of proven design and trusted performance.

Unmatched design

- Ergonomically shaped for comfortable use
- Promotes confident operation and handling

Saves time

- Single-hand, one-step operation
- Reduced training time
- No special staging, ready to use right out of the box

Greater accuracy

- Correct product selection is easier with color codes
- Top mounted activation button maximizes visibility for precise placement

Promotes safety

- Specifically designed for single use only, preventing reuse
- Automatic retraction of lancet minimizes exposure of sharp and possible contamination

Catalog Number	Lancet Color	Width/Gauge	Depth	Blood Volume
200101	Pink	1.5mm	1.0mm	Low Blood flow
200102	Green	1.5mm	1.5mm	Medium Blood flow
200103	Blue	1.5mm	2.0mm	High Blood flow
200104	Orange	23G	2.25mm	Single Blood Drop
200105	Purple	28G	1.25mm	Single Blood Drop

Packaging: • 200 Lancets per box
• 10 boxes per case—2,000 Lancets



Capiject Capillary Micro Collection Tubes.

Unmatched design and unparalleled performance.

With a full compliment of collection tubes for serum, plasma, and whole blood determinations, a small but adequate amount of blood can be easily obtained by skin puncture with a Capiject Safety Lancet.

Capiject® makes it easy.

- Extra-wide, built in collection tip for enhanced specimen flow
- Single piece, tapered design for easy collection and handling
- Color-coded fill lines and caps for volume accuracy and easy identification
- Lot number and expiration date printed on each tube for easy identification
- Silicone coated anticoagulant tubes designed to reduce clotting
- Olefin-based gel for clean separations
- Tube nests in cap allowing tube to stand alone
- Certified for lead testing
- Latex-free

Catalog Number	Tube Color	Tube Additive	Application	Lead Testing Certified	Fill Volume µL
T-MG	Cherry Red	Gel Barrier with Clot Activator (silica particles and inert gel)	Serum determinations in clinical chemistry	X	500
T-MQK	Lavender	K2-EDTA (0.78 mg)	Whole blood hematology procedures	X	500
T-M	Brick Red	No Additive/Clot Tube	Serum determinations in clinical chemistry, serology, and blood banking	X	500
T-MLH	Green	Lithium Heparin (12.5 IU)	Plasma determinations in clinical chemistry	X	500
T-MLHG	Mint Green	Lithium Heparin with Gel (12.5 IU)	Plasma determinations in clinical chemistry	X	500
T-MLHGA	Mint Green	Amber tube with Lithium Heparin and Gel Barrier (12.5 IU)	Light-sensitive plasma determinations in clinical chemistry	X	500
T-MPS	Gray	Sodium Fluoride (0.72 mg) Potassium Oxalate (0.58 mg)	Glucose determinations	X	250
T-MGA	Gold	Amber Tube with Gel Barrier and Clot Activator (silica particles and inert gel)	Light-sensitive serum analyte determinations (Bilirubin determinations)	X	500

For in vitro diagnostic use

Packaging: • 50 Capiject Capillary Blood Collection Tubes per bag
• 5 bags per box—250 Capiject Tubes
• 4 boxes per case—1,000 Capiject Tubes