

Modular Wall Cabinet

Model 6074

UMF Medical's 6074 Modular Wall Cabinet has 2-doors, and an adjustable stainless steel shelf. The unit includes wall mounting hardware. It is the perfect complement to any exam room or working area in need of modular cabinet options.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel construction
- › Commercial grade, adjustable, spring assisted door hinges
- › Adjustable shelf
- › Wall mounting bracket
- › 3-year limited warranty

Specifications

Width: 24" (61.0cm)
Height: 20" (50.8cm)
Depth: 11" (27.9cm)
Net Weight: 40lb (18.1kg)

Shipping Weight: 70lb (31.8kg)
Shipping Dimensions:
L36xW28xH31
(91.4cm x 71.1cm x 78.7cm)

Accessories

502 - Double door lock (*Factory Installed Only*)
507 - Additional keys (x2)