

Convenient Storage Cabinet



Shown with optional basket and pad

Large Storage Cabinet w/ Two Doors
Basket Can Be Placed Horizontally or
Elevated On One End
3" Ball Bearing Swivel Casters For Safe
Transport
All-Welded Stainless Steel

The UMF Medical SS8528 Stainless Steel Bassinet with two doors is built on a heavy duty 1" O.D. (2.54cm) tubular stainless steel framework. The basket supports with insert bends allow basket to be moved to horizontal or elevated positions. It is mounted on 3" (7.6cm) ball bearing swivel casters.

Basket and mattress sold separately.

3 year limited warranty included.

Standard Features

- Heavy duty tubular stainless steel framework
- Full width storage area with two doors
- 3" (7.6cm) ball bearing swivel casters
- 3 year limited warranty

Specifications

- Height: 35.375" (89.9cm)
- Width: 31.25" (79.4cm)
- Depth: 17" (43.2cm)

Options

- 132 set - 3" swivel casters with brake
- 8545 - Basket
- 8546 - Pad