

## Green Series™ Medical Exam Lights

Patented LED technology that is whiter, brighter and greener than halogen.

### Green Series Lights

- Touchless controls
- Lower cost of ownership
- Improve patient exams



Learn more  
on page 181



## Exam Lights

Green Series™ Exam Light IV.....	178
Green Series 300 General Exam Light.....	179

## Minor Procedure Lighting

Green Series™ 600 Minor Procedure Light.....	180
Green Series 900 Procedure Light.....	181

## Diagnostic Headlights

Green Series™ Portable Headlights.....	182
--	-----

## Procedure Headlights

Green Series™ Procedure Headlights.....	183
---	-----

## Skin Surface Microscopes

3.5 V EpiScope® Skin Surface Microscope.....	184
--	-----

## PenLite

Professional PenLite.....	185
---------------------------	-----



**48810**  
GS General Exam  
Light IV with Table/  
Wall Mount

**48950**  
Mobile Stand for GS  
Exam Light IV

## Green Series™ Exam Light IV

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas, the GS Exam Light IV provides you with 3X the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes. And with a compact, unobtrusive head, it ensures pinpoint accuracy and provides true tissue color rendition during an exam. Utilizing white LEDs, the GS Exam Light IV provides a whiter, brighter light than halogen and consumes less energy.

- One LED—3X the typical light output; pinpoint accuracy; true tissue color rendition
- Highly flexible pipe for easy maneuverability
- Lightweight pipe is twice as long as conventional pipes
- Touchless On/Off
- Multiple mounting options for floor space maximization
- Disposable sheath available to help reduce the risk of cross-contamination
- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble

**48810** GS General Exam Light IV with Table/Wall Mount

### Green Series Exam Light IV: Accessories

**48950** Mobile Stand for GS Exam Light IV

**48955** Table/Wall Mount for GS Exam Light IV

**52640** Disposable Sheath for GS Exam Light IV (25/box, 5 boxes/case)



**48810**  
GS General Exam  
Light IV with Table/  
Wall Mount



**48950**  
Mobile Stand for  
GS Exam Light IV



**48955**  
Table/Wall Mount for  
GS Exam Light IV

## Green Series™ 300 General Exam Light

Utilizing a single white LED, the GS 300 Exam Light provides a whiter, brighter light than halogen and consumes less energy.

- Perfect combination of intense light output and cool operation in a compact design.
- Affordable price and rugged performance
- Touchless On/Off
- Multiple mounting options for floor space maximization
- Low cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy assembly

**44400** GS 300 General Exam Light with Mobile Stand

**44410** GS 300 General Exam Light with Table/Wall Mount

### Green Series 300 General Exam Light: Accessories

**48955** Table/Wall Mount for GS 300 Exam Light



**44400**  
GS 300 with  
Mobile Stand



**48955**  
Table/Wall Mount  
for GS 300



**44410**  
GS 300 with Table/  
Wall Mount

# MINOR PROCEDURE LIGHTING

## LIGHTING

### Green Series™ 600 Minor Procedure Light

Featuring three LEDs and a wider head, the GS 600 is designed specifically for demanding viewing conditions.

- Three White LEDs for greater intensity
- Touchless power controls and adjustable light output
- Disposable sheath available to help reduce the risk of cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Low cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble and clean



**44600**  
GS 600 with  
Mobile Stand

**44600** GS 600 Minor Procedure Light with Mobile Stand

**44610** GS 600 Minor Procedure Light with Table/Wall Mount

### Green Series 600 Minor Procedure Light: Accessories

**48955** Table/Wall Mount for GS 600

**52630** Disposable Sheath for GS 600 (50/box, 5 boxes/case)



**44610**  
GS 600 with Table/Wall Mount



**48955**  
Table/Wall Mount for GS 600

### Green Series™ 900 Procedure Light

This procedure light, featuring six LEDs, combines outstanding durability, extraordinary maneuverability, unparalleled spot quality and multiple mounting options.

- Six white LEDs for greater intensity
- Touchless On/Off with dimming option
- Multiple mounting options for floor space maximization
- Disposable sheath available to help reduce the risk of cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Lower cost of ownership due to elimination of bulb replacement and lower energy consumption

**44900** GS 900 Procedure Light with Mobile Stand

**44900-C** GS 900 Procedure Light with Ceiling Mount

**44900-W** GS 900 Procedure Light with Wall Mount



**44900**  
GS 900 Procedure Light  
with Mobile Stand

### Green Series 900 Procedure Light: Accessories

**52630** Disposable Sheath for GS 900 (50/box, 5 boxes/case)



**44900-C**  
GS 900 Procedure Light  
with Ceiling Mount



**44900-W**  
GS 900 Procedure Light  
with Wall Mount

# DIAGNOSTIC HEADLIGHTS

## LIGHTING



## Green Series™ Portable Headlights

Provide illumination precisely where it's needed with cool, bright light that ensures comfort for both doctor and patient.

- Weighing approximately four ounces, it's one of the most comfortable and lightweight headlights available
- Fits inside a shirt pocket and allows complete freedom of movement
- Supplies bright, white, shadow-free light for true tissue color
- Green Series lamp has 10,000-hour life expectancy with minimal degradation of light output

**46070** Green Series Portable Headlight with Direct Power Supply/Charger

**46070V** Veterinary Green Series Portable Headlight with Direct Power Supply/Charger

### Rigid Headband

**46070R (Domestic)** Green Series Portable Headlight with Direct Power Supply/Charger

### Solid State Portable Headlights: Accessories

**05232-U** Soft Carrying Case

**49642** Soft Terrycloth Headband (case of 10)

**72270** Rechargeable Battery

**74170** Direct Power Source/Charger

### Green Series™ Procedure Headlights

Provide cool, bright white light that is ideal for a variety of office-based procedures. When you need targeted, well-defined illumination without the heat of traditional headlights, look to our Green Series headlights for superior quality and reliability.

- Weighing approximately six ounces, it is extremely comfortable to wear for simple and complex procedures
- Supplies bright, white, shadow-free light for true tissue color
- Rechargeable battery pack makes this a completely portable, compact unit that allows freedom of movement
- Enjoy the flexibility of operating the headlight using the direct power source or the rechargeable “belt-clip” portable power source
- Solid state lamp has 10,000-hour life expectancy with minimal degradation of light output

**49020** Green Series Procedure Headlight with Direct Power Supply

### Green Series Procedure Headlights: Accessories

**49099** Carrying Case

**75260** Portable Power Source with Charger

**49000** Procedure Headlight (no power supply)

**49097** Set, Front and Rear Pad

**72250** Replacement Cells for Battery Pack

**74180** Charger only for Portable Power Source

**75270** Direct Power Supply

**75200** Portable Power Source only

Operate the Headlight from a Belt-Clip or Direct Power Source

Weighs Approximately 6 oz



Solid-State Lamp

Bright, White, Shadow-Free Light

**49020**  
Green Series Procedure Headlight with Direct Power Supply



**49099**  
Carrying Case



**75200**  
Portable Power Source only



**74180**  
Charger only for Portable Power Source





**47300**  
EpiScope® only

### 3.5 V EpiScope® Skin Surface Microscope

Helps distinguish pigmented skin lesions such as basal cell carcinoma vs. malignant melanoma by providing the ability to magnify, illuminate and measure.

- Provides 10X times magnification and illumination for crisp, clear viewing
- Reticle contact plate facilitates measurement and tracking of lesion growth
- Fits any 3.5 V Welch Allyn power source
- Halogen illumination for true tissue color
- Free reference guide available for patient education

**47300** EpiScope only

**47351** EpiScope Set with Rechargeable Handle and Soft Carrying Case

**47351-C** EpiScope Set with Convertible Handle



**05232-U**  
Soft Carrying Case

### 3.5 V EpiScope Skin Surface Microscope: Accessories

**05232-U** Soft Carrying Case

**47310** Reticle Contact Plate for EpiScope

**47320** Plano Contact Plate for EpiScope

**03000-U** 3.5 V Halogen Lamp

## Professional PenLite

For general non-contact illumination.

- Bright, white halogen light
- Durable solid brass construction

**76600** Halogen Professional PenLite

### Professional PenLite: Accessories

**07600-U** 2.5 V Halogen Lamp



Halogen Light

**76600**  
Halogen Professional PenLite